# CanAg CEA EIA

REF 401-10

IVD

Instructions for use. 2022-06

Read highlighted changes

EN	EXPLANATION OF SYMBOLS
BG	ОБЯСНЕНИЕ НА СИМВОЛИТЕ
CS	VÝZNAM SYMBOLŮ
DA	SYMBOLFORKLARING
DE	ERKLÄRUNG DER SYMBOLE
EL	ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES	SIGNIFICADO DE LOS SÍMBOLOS
ET	SÜMBOLITE SELGITUS
FR	EXPLICATION DES SYMBOLES
HR	OBJAŠNENJE SIMBOLA
HU	JELMAGYARÁZAT
IT	SPIEGAZIONE DEI SIMBOLI
LT	SIMBOLIŲ PAAIŠKINIMAI
LV	SIMBOLU SKAIDROJUMS
NL	VERKLARING DER SYMBOLEN
NO	SYMBOLFORKLARING
PL	OBJAŚNIENIE SYMBOLI
PT	EXPLICAÇÃO DOS SÍMBOLOS
RO	SEMNIȚAȚIA SIMBOLURILOR
RU	ОБОЗНАЧЕНИЯ
SV	SYMBOLFÖRKLARING
SK	VÝZNAM SYMBOLOV
SL	RAZLAGA SIMBOLOV
SR	OBJAŠNENJE SIMBOLA
TR	SEMBOLLERİN AÇIKLAMALARI

Use By/Годно до/Použitelné do/  
Holdbar til/Verwendbar bis/  
Ημερομηνία λήξης/Fecha  
de caducidad/Kölblik kuni/  
Utiliser jusque/Rok valjanosti/  
Felhasználható/Utilizzare entro/  
Sunautoti iki/Izlietot līdz/Houdbaar  
tot/Brukes innen/Uzyć przed/  
Prazo de validade/Expírã la/  
Использовать до/Använd före/  
Použite né do/ Uporabno do/  
Upotrebljivo do/Son Kullanna Tarihi

LOT

Batch code/Номер на партида/  
Číslo šarže/Lotnummer/  
Chargenbezeichnung/Αριθμός  
Παρτίδας/Código de lote/Partii  
kood/Code du lot/Kod serije/  
Sarzsám/Codice del lotto/  
Partijas kods/Partijas kods/Lot  
nummer/Partikode/Kod partii/  
Código do lote/Număr de lot/  
Номер лота/Lotnummer/Číslo  
šarže/Številka serije/Kod partije/  
Parti Kodu

Date of manufacture/Дата на производство/Datum výroby/  
Produktionsdato/Herstellungsdatum/  
Ημερομηνία παραγωγής/Fecha de fabricación/Valmistamise kuupäev/  
Date de fabrication/Datum proizvodnje/  
Gyártási idő/Data di produzione/  
Pagaminimo data/Ražošanas datums/  
Productiedatum/Fremstillingsdato/  
Data produkcji/Data de fabrico/Data fabricației/Дата производства/  
Tillverkningsdatum/Dátum výroby/Datum izdelave/Datum proizvodnje/Úretim tarihi

Temperature limitation/  
Температурни граници/  
Теплотни омеzeи/  
Temperaturbegrænsning/  
Temperaturbegrenzung/  
Περιορισμοί θερμοκρασίας/  
Limites de temperatura/  
Temperatuuri piirang/  
Limite de température/  
Temperaturno ograničenje/  
Hőmérsékletre vonatkozó korlátozás/  
Limiti di temperatura/  
Temperatūriniai apribojimai/  
Temperatūras ierobežojums/  
Temperaturbepèrking/  
Temperaturbegrensninger/  
Temperaturey granicne/  
Limite de temperatura/  
Limite de temperatură/  
Температурный режим/  
Temperaturbegrænsning/  
Теплотне обмеzenie  
Omejitve temperature/  
Temperaturno ograničenje/  
Sıcaklık sınırlaması/

## IVD

In Vitro Diagnostic Medical Device/  
Медицински уред за диагностика  
ин vitro/Diagnostický zdravotnícký  
prostředek in vitro/Medicinsk udstyr til  
in vitro-diagnostik/In-vitro-Diagnostikum/  
Ιατροτεχνολογικό προϊόν για διάγνωση  
In Vitro/Dispositivo médico para  
diagnóstico in vitro/In vitro diagnostiline  
meditsiiniseade/Dispositif médical de  
diagnostic in vitro/Diagnostički medicinski  
uređaj In Vitro/In vitro orvosdiagnostikai  
eszköz/Dispositivo medico per test  
diagnostici in vitro/In Vitro Diagnostinė  
Medicinos Priemonė/Medicínska ierice  
in vitro diagnostikai/In vitro-diagnostisch  
medisch instrument/In vitro diagnostisk  
medisinsk utstyr/Wyrób medyczny do  
diagnostyki in vitro/Dispositivo Médico  
de Diagnóstico In Vitro/Dispozitiv medical  
pentru diagnostic in vitro/Только для  
диагностики In Vitro/Endast för in  
vitro-diagnostik/ Zdravotnícka pomôcka na  
diagnostiku in vitro/In vitro diagnostični  
pripomoček/Diagnostički medicinski  
uređaj In Vitro/<96> testleri için yeterlilik  
içerir

Contains sufficient for <96> tests/Съдържа  
достатъчно количество за тестове  
<96>/Lze použít pro <96> testů/Ineholder  
tilstrækkeligt/Inhalt ausreichend für <96>  
Prüfungen/Περισχόμενο επαρκές για  
«96» εξετάσεις/Contenido suficiente para  
<96> ensayos/Kogusest piisab <96> testi  
lääbiviimiseks/Contenu suffisant pour «96»  
tests/Sadržaj dovoljno za <96> testova/A  
doboz tartalma <96> vizsgálat elvégzéséhez  
elegendő/Contenuto sufficiente per «96»  
saggi/Turiny's skirtas atlikti <96> tyrimus/  
Satur's pietiekams <96> testiem/Inhoud  
voldoende voor «96» testen/til «96» test/  
Tilstrækkelig innhold for <96> prøver/  
Wystarczy na wykonanie <96> testów/  
Conteúdo suficiente para «96» ensaios/  
Conținut suficient pentru 96 de teste/  
Содержит достаточные количества для  
«96» определений/Innehåller tillräckligt  
till «96» antal tester/Obsah postačuje na  
tento počet testov: <96>/Vsebinsa zadostuje  
za <96> testov/Sadržina dovoljna za <96>  
testova/<96> testleri için yeterlilik içerir

## REF

Catalogue number/Каталожен номер/  
Katalogové číslo/Katalognummer/  
Bestellnummer/Αριθμός καταλόγου/  
Número de catálogo/Katalogi number/  
Numéro de catalogue/Kataloški broj/  
Katalógusszám/Numero di catalogo/  
Katalogo numeris/Numurs katalogā/  
Catalogusnummer/Katalognummer/  
Numer katalogowy/Número do catálogo/  
Număr de catalog/Номер по каталогу/  
Produktnummer/Katalógové číslo/  
Kataloška številka/Kataloški broj/  
Katalog numarası

Consult Instructions for Use/  
Прочетете инструкцията за  
употреба/Konzultujte s návodem  
k použití/Se brugsanvisning/Siehe  
Gebrauchsanweisung/Συμβουλευτείτε  
της Οδηγίας σχετικά με τη χρήση/  
Consulte las instrucciones de uso/  
Vt kasutusjuhendit/Consulter le mode  
d'emploi/Pročítajte upute za uporabu/  
Olvassa el a használati utasítást/  
Consultare le istruzioni per l'uso/Dél  
naudojimo žiūrėkite instrukcijas/Izlasiet  
lietošanas instrukciju/Raadpleeg de  
instructies voor gebruik/Les instruksene  
for bruk/Sprawdzić w instrukcji użycia/  
Consulte as Instruções de Utilização/  
Consultați instrucțiunile de utilizare/  
Обратитесь к инструкции по  
применению/Se bruksanvisning/  
Prečítajte si návod na používanie/  
Pročítajte uputstvo za upotrebu/  
Kullanım Talimatlarını Bakınız

## CONT

Contents of kit/Съдържание на набора/  
Obsah soupravy/Kittets indhold/Inhalt  
des Kits/Περιεχόμενα του κιτ/Contenido  
del kit/Komplekt sisaldab/Contenu du  
kit/Sadržaj opreme/A készlet tartalma/  
Contenuto del kit/Rinkinio turinys/  
Komplekta saturs/Inhoud van de set/  
Settets innhold/Zawartość zestawu/  
Conteúdo do kit/Conținutul setului/  
Компоненты набора/Kit innehåll/  
Obsah súpravy/Vsebina kompleta/Sadržaj  
opreme/Kitin içindekiler

Biological risks/Биологическа  
опасност/Biológická rizika/Biologisk  
fare/Biologische Gefahren/Βιολογικοί  
κίνδυνοι/Riesgos biológicos/  
Biolooigilised ohud/Risques biologiques/  
Biolóškli rizici/Biológiai kockázatok/Rischi  
biologici/Biologinis pavojus/Biolóškisks  
risks/Biologische risico's/Biologische  
risikoer/Zagroženie biologiczne/Riscos  
biológicos/ Biologisk risk/Pericole  
biologice/Биологическая опасность/  
Biologicky rizikové/Biológické riziká/  
Biolóškli rizici/Biyolojik riskler

## ORIG HUM

Human/C човешки производ/Ľidské/  
Human/Human/δείγματα αναφοράς/  
Humano/Inimāritolu/Humaine/Ljudskog  
porjekla/Humán/Origine Umana/  
Žmogaus kilmės/Cilvēku izcelsmes/  
Human/Menneske/Ludzka/Humano/  
Origine umână/Человеческого  
происхождения/Human/Ludské/  
Humanega izvora/Ljudskog porekla/Ľnsan

## ORIG MOU

From mouse/C миши производ/Myši/  
Fra mus/Maus/από πογτίκι/de ratón/  
Hiirtelt/De souris/Mišijeg porjekla/  
Egérböli/Murino/Pelés kilmés/No peles/  
Van muizen/Fra mus/Mysia/Do rato/De  
la șoareci/Мышиного происхождения/  
Från mus/Myšije/Mišjega izvora/Mišijeg  
porekla/Fareden

## ORIG BOV

Bovine/C говежди производ/  
Hovēži/Bovin/Rind/από βοοειδή/  
Bovino/Veistelt/Bovine/Rogate stoke/  
Szarvasmarha/Bovino/Jaučio/No  
liellopa/Bovien/Bovin/Wolowy/Bovino/  
Origine bovină/крупного рогатого  
скота/Från ko/Hovädzie/Rogvega  
izvora/Rogate krupne stoke/Bovin

Reconstitute with/Пазтваряне с/  
Rozfeďte pomoci/Rekonstitueres med/  
Rekonstituieren mit/Ανασύσταση με/  
Reconstituir con/Lahjendamine/  
Reconstituer avec/Rekonstituiraite s/  
Feloldáshoz/Ricostituire con/Atkurti,  
ištirpdant su/Atšķaidīt ar/Reconstituire  
met/Rekonstitueres med/Odtworzyć  
za pomocą/Reconstituir com/A  
se reconstitui cu/Пастворить в/  
Rekonstituera med/Rozriedte pomocou/  
Rekonstituiraite z/s/Ponovno formiranje  
sa/Yeniden oluşturalur

Manufacturer/Производител/Výrobce/  
Producent/Hersteller/Κατασκευαστής/  
Fabricante/Tootja/Fabricant/Proizvođač/  
Gyártó/Fabbricante/Gamintojas/  
Ražotājs/Fabrikant/Produsent/  
Producent/Fabricante/Producător/  
Производитель/Tilverkare/ Výrobca/  
Izdelovalec/Proizvođač/Üretici

# INSTRUCTIONS FOR USE

EN

## INSTRUCTIONS FOR USE

Please visit our website [www.fdi.com/ifu](http://www.fdi.com/ifu) to obtain the Instructions For Use (IFU) in additional languages.

To ensure that you download the correct IFU for your kit lot, please select the revision corresponding to the issue date printed on the front page of the IFU provided with this kit.

Please follow the IFU carefully. Instructions for safe handling are found in the WARNINGS AND PRECAUTIONS section. Material Safety Data Sheets (MSDS) are available on our website [www.fdi.com](http://www.fdi.com). If you do not have access to the internet, please contact your local distributor, or Fujirebio Diagnostics AB for assistance.

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# CanAg CEA EIA

Instructions for use

Enzyme immunometric assay kit  
For 96 determinations

## INTENDED USE

The CanAg CEA EIA kit is intended for the quantitative determination of the cancer associated antigen CEA in serum.

## SUMMARY AND EXPLANATION OF THE ASSAY

Carcinoembryonic antigen (CEA) is a glycoprotein, which was first identified in patients with colonic carcinoma and in epithelial tumours of endodermal origin (gastrointestinal tract) by Gold and Freedman (1). The CEA molecule is quite heterogeneous due to the carbohydrate contents (50-60%) and depending on the purification procedure employed. It is soluble in perchloric acid and has a molecular weight of about 175.000–200.000 Daltons (2). Immunological and genetic characterization of CEA has identified a family of CEA-like molecules sharing common antigenic determinants. The most relevant CEA-like molecule is NCA (non-specific cross-reacting antigen) synthesized both by normal and pathological tissues. The problem of cross-reacting CEA-like molecules when assaying CEA is possible to overcome by the use of monoclonal antibodies. The CanAg CEA EIA is based on two mouse monoclonal antibodies against the Gold epitopes IV and V (3, 4).

CEA is secreted from tumour cells and is a widely used serological marker of gastrointestinal carcinomas, lung cancer and breast cancer. In colorectal cancer, the clinical use of CEA testing for monitoring response to therapy and for documenting progressive disease is well established (5, 6). CEA may also be present in benign gastrointestinal inflammatory diseases or in hepatobiliary diseases. These observations make it necessary to emphasize that the CEA assay should not be used as a cancer-screening test.

## PRINCIPLE OF THE TEST

The CanAg CEA EIA is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. Calibrators, controls and patient samples are incubated together with biotinylated Anti-CEA monoclonal antibody and horseradish peroxidase (HRP) labelled Anti-CEA monoclonal antibody in Streptavidin coated microstrips. After washing, buffered Substrate/ Chromogen reagent (hydrogen peroxide and 3, 3', 5, 5' tetra-methylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction a blue colour will develop if antigen is present.

The intensity of the colour is proportional to the amount of CEA present in the samples.

The colour intensity is determined in a microplate spectrophotometer at 620 nm (or optionally at 405 nm after addition of Stop Solution). Calibration curves are constructed for each assay by plotting absorbance value versus the concentration for each calibrator. The CEA concentrations of patient samples are then read from the calibration curve.

## REAGENTS

- Each CanAg CEA EIA kit contains reagents for 96 tests.
- The expiry date of the kit is stated on the label on the outside of the kit box.
- Do not use the kit beyond the expiry date.
- Do not mix reagents from different kit lots.
- Store the kit at 2–8° C. Do not freeze.
- Opened reagents are stable according to the table below provided they are not contaminated, stored in resealed original containers and handled as prescribed. Return to 2–8° C immediately after use.

Component	Quantity	Storage and stability after first opening
-----------	----------	---

### MICROPLA

<b>Microplate</b>	1 Plate	2–8° C until expiry date stated on the plate
-------------------	---------	--

12 x 8 wells coated with Streptavidin. After opening, immediately return unused strips to the aluminium pouch, containing desiccant. Reseal carefully to keep dry.

<b>CEA Calibrators</b>	6 vials	2–8° C until expiry date stated on the vials
------------------------	---------	--

CAL	CEA	0
-----	-----	---

0 µg/L 1 x 8 mL

CAL	CEA	2
-----	-----	---

2 µg/L 1 x 0.75 mL

CAL	CEA	5
-----	-----	---

5 µg/L 1 x 0.75 mL

CAL	CEA	15
-----	-----	----

15 µg/L 1 x 0.75 mL

CAL	CEA	50
-----	-----	----

50 µg/L 1 x 0.75 mL

CAL	CEA	75
-----	-----	----

75 µg/L 1 x 0.75 mL

Human CEA in a Tris-HCl buffered salt solution containing bovine serum albumin, an inert yellow dye and 0.01 % methyl-isothiazolone (MIT) as preservative. Ready for use. 

CAL	CEA	0
-----	-----	---

 should also be used for dilution of samples.

Component	Quantity	Storage and stability after first opening			
<b>CEA Controls</b>	2 vials	2–8° C until expiry date stated on the vials			
<table border="1" style="display: inline-table;"><tr><td>CONTROL</td><td>CEA</td><td>1</td></tr></table>	CONTROL	CEA	1	1 x 0.75 mL	
CONTROL	CEA	1			
<table border="1" style="display: inline-table;"><tr><td>CONTROL</td><td>CEA</td><td>2</td></tr></table>	CONTROL	CEA	2	1 x 0.75 mL	
CONTROL	CEA	2			

Human CEA in a Tris-HCl buffered salt solution containing bovine serum albumin, and 0.01 % methyl-isothiazolone (MIT) as preservative. Ready for use.

<b>BIOTIN</b>	<b>Anti-CEA</b>
---------------	-----------------

<b>Biotin Anti-CEA</b>	1 x 15 mL	2–8° C until expiry date stated on the vial
------------------------	-----------	---

Biotin Anti-CEA monoclonal antibody from mouse, approximately 3 µg/mL. Contains phosphate buffered saline (pH 7.2), bovine serum albumin, bovine immunoglobulin, blocking agents, Tween 20, an inert blue dye and 0.01 % methyl-isothiazolone (MIT) as preservative. To be mixed with Tracer, HRP Anti-CEA before use.

<b>CONJ</b>	<b>Anti-CEA</b>
-------------	-----------------

<b>Tracer, HRP Anti-CEA</b>	1 x 0.75 mL	2–8° C until expiry date stated on the vial
-----------------------------	-------------	---

Stock solution of HRP Anti-CEA monoclonal antibody from mouse, approximately 60 µg/mL. Contains preservatives. To be mixed with Biotin Anti-CEA before use.

<b>SUBS</b>	<b>TMB</b>
-------------	------------

<b>TMB HRP-Substrate</b>	1 x 12 mL	2–8° C until expiry date stated on the vial
--------------------------	-----------	---

Contains buffered hydrogen peroxide and 3, 3', 5, 5' tetramethyl-benzidine (TMB). Ready for use.

Component	Quantity	Storage and stability after first opening
-----------	----------	---

**STOP**

<b>STOP Solution</b>	1 x 15 mL	2–8° C until expiry date stated on the vial
----------------------	-----------	---

Contains 0.12 M hydrochloric acid. Ready for use.

**WASHBUF 25X**

<b>Wash Concentrate</b>	1 x 50 mL	2–8° C until expiry date stated on the bottle
-------------------------	-----------	---

A Tris-HCl buffered salt solution with Tween 20. Contains Germall II as preservative. To be diluted with water 25 times before use.

### Indications of instability

The TMB HRP-Substrate should be colourless or slightly bluish. A blue colour indicates that the reagent has been contaminated and should be discarded.

### WARNINGS AND PRECAUTIONS

#### For in vitro diagnostic use.

- For professional use only.
- Please refer to the US Department of Health and Human Services (Bethesda, Md., US) publication No. (CDC) 88-8395 on laboratory safety or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Follow local guidelines for disposal of all waste material.

#### Caution

Material used in the preparation of human source reagent has been tested and found to be Non Reactive for HIV-1/2 Antibody, HCV Antibody and Hepatitis B Surface Antigen (HBsAg). Since no method can completely rule out the presence of blood borne diseases, the handling and disposal of human source reagents from this product should be made as if they were potentially infectious.

#### CLP (1272/2008) HAZARD CLASSIFICATION

Information about CLP (1272/2008) HAZARD CLASSIFICATION can be found at the end of this document.

# Protocol Sheet

## CanAg CEA EIA REF 401-10

Mix the components directly before use. Use shaking conditions according to the Instructions.

Step	Vial/Plate	Procedure																																	
1. Prepare Wash Solution	WASHBUF 25X	Dilute 50 mL of Wash Concentrate with 1200 mL of distilled water or deionized water.																																	
Prepare Antibody Solution	CONJ Anti-CEA	Mix 50 $\mu$ L of Tracer, HRP Anti-CEA, with 1 mL of Biotin Anti-CEA per strip:																																	
	BIOTIN Anti-CEA																																		
		<table border="1"> <thead> <tr> <th>No. of Strips</th> <th>Tracer, HRP Anti-CEA (<math>\mu</math>L)</th> <th>Biotin Anti-CEA (mL)</th> </tr> </thead> <tbody> <tr><td>1</td><td>50</td><td>1</td></tr> <tr><td>2</td><td>100</td><td>2</td></tr> <tr><td>3</td><td>150</td><td>3</td></tr> <tr><td>4</td><td>200</td><td>4</td></tr> <tr><td>5</td><td>250</td><td>5</td></tr> <tr><td>6</td><td>300</td><td>6</td></tr> <tr><td>7</td><td>350</td><td>7</td></tr> <tr><td>8</td><td>400</td><td>8</td></tr> <tr><td>9</td><td>450</td><td>9</td></tr> <tr><td>10</td><td>500</td><td>10</td></tr> </tbody> </table>	No. of Strips	Tracer, HRP Anti-CEA ( $\mu$ L)	Biotin Anti-CEA (mL)	1	50	1	2	100	2	3	150	3	4	200	4	5	250	5	6	300	6	7	350	7	8	400	8	9	450	9	10	500	10
No. of Strips	Tracer, HRP Anti-CEA ( $\mu$ L)	Biotin Anti-CEA (mL)																																	
1	50	1																																	
2	100	2																																	
3	150	3																																	
4	200	4																																	
5	250	5																																	
6	300	6																																	
7	350	7																																	
8	400	8																																	
9	450	9																																	
10	500	10																																	

## SPECIMEN COLLECTION AND HANDLING

The CanAg CEA EIA is intended for use with serum. Collect blood by venipuncture and separate the serum according to common procedures. Samples can be stored at 2–8° C for 2 days. For longer periods it is recommended to store the samples at –20° C or below. Avoid repeated freezing and thawing of the samples. Allow frozen samples to thaw slowly, preferably at 2–8° C over night and then bring the samples to room temperature before analysis.

## PROCEDURE

### Materials required but not supplied with the kit

**1. Microplate shaker**

Shaking should be medium to vigorous, approximately 700-1100 oscillations/min.

**2. Microplate washer**

Automatic plate washer capable of performing 1 and 6 washing cycles with a minimal fill volume of 350 µL/well/washcycle.

An 8-channel pipette with disposable plastic tips for delivery of 350 µL is recommended if an automatic microplate washer is not used.

**3. Microplate spectrophotometer**

With a wavelength of 620 nm and/or 405 nm and an absorbance range of 0 to 3.0.

**4. Precision pipettes**

With disposable plastic tips to deliver microlitre and millilitre volumes.

An 8-channel pipette or respenser pipette with disposable plastic tips for delivery of 100 µL is useful but not essential.

**5. Distilled or deionized water**

For preparation of Wash Solution.

### Procedural notes

1. A thorough understanding of this package insert is necessary to ensure proper use of the CanAg CEA EIA kit. The reagents supplied with the kit are intended for use as an integral unit. Do not mix identical reagents from kits having different lot numbers. Do not use the kit reagents after the expiry date printed on the outside of the kit box.
2. Reagents should be allowed to reach room temperature (20–25° C) prior to use. The assay should only be performed at temperatures between 20–25° C to obtain accurate results. Frozen specimens should be brought to room temperature slowly and must be gently but thoroughly mixed after thawing.

3. Before starting to pipette calibrators, controls and patient specimens it is advisable to mark the strips to be able to clearly identify the samples during and after the assay.
4. The requirement for efficient and thorough washing for separation of bound and unbound antigen and reagents from the solid-phase bound antibody-antigen complexes is one of the most important steps in an EIA. In order to ensure efficient washing make sure that all wells are completely filled to the top edge with wash solution during each wash cycle, that wash solution is dispensed at a good flow rate, that the aspiration of the wells between and after the wash cycles is complete and that the wells are empty. If there is liquid left, invert the plate and tap it carefully against absorbent paper.
  - Automatic strip washer: Follow the manufacturer's instructions for cleaning and maintenance diligently and wash the required number of wash cycles prior to and after each incubation step. It's highly recommended to use *strip* process mode and *overflow* wash mode with a dispensing volume of 800  $\mu$ L. The aspiration/wash device should not be left standing with the Wash Solution for long periods, as the needles may get clogged resulting in poor liquid delivery and aspiration.
5. The TMB HRP-Substrate is very sensitive for contamination. For optimal stability of the TMB HRP-Substrate, pour the required amount from the vial to a carefully cleaned reservoir or preferably a disposable plastic tray to avoid contamination of the reagent. Be sure to use clean disposable plastic pipette tips (or respenser pipette tip).
6. Be sure to use clean disposable plastic pipette tips and a proper pipetting technique when handling samples and reagents. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid. A proper pipetting technique is of particular importance when handling the TMB HRP-Substrate Solution.

Preparation of reagents	Stability of prepared reagent
<b>Wash Solution</b>	2 weeks at 2–25° C in a sealed container

Pour the 50 mL Wash Concentrate into a clean container and dilute 25- fold by adding 1200 mL of distilled or deionized water to give a buffered Wash Solution.

<b>Antibody Solution</b>	3 weeks at 2–8° C
Prepare the required quantity of Antibody Solution by mixing 50 µL of Tracer, HRP Anti-CEA with 1 mL of Biotin Anti-CEA per strip (see table below and the Protocol Sheet).	

No. of Strips	Tracer, HRP Anti-CEA (µL)	Biotin Anti-CEA (mL)
1	50	1
2	100	2
3	150	3
4	200	4
5	250	5
6	300	6
7	350	7
8	400	8
9	450	9
10	500	10
11	550	11
12	600	12

Be sure to use a clean plastic or glass bottle for preparation of the Antibody Solution.

**Alternative:** Pour the content of the Tracer, HRP Anti-CEA into the vial of Biotin Anti-CEA and mix gently. Make sure that all of the Tracer, HRP Anti-CEA is transferred to the vial of Biotin Anti-CEA.

**NOTE:** The Antibody Solution is stable for 3 weeks at 2–8° C. Do not prepare more Antibody Solution than will be used within this period and make sure that it is stored properly.

### Assay procedure

Perform each determination in duplicate for calibrators, controls and patient samples. A calibration curve should be run with each assay. All reagents and samples must be brought to room temperature (20–25° C) before use.

1. Start to prepare Wash Solution and Antibody Solution. It is important to use clean containers. Follow the instructions carefully.

- Transfer the required number of microplate strips to a strip frame. (Immediately return the remaining strips to the aluminium pouch containing a desiccant and reseal carefully). Wash each strip once with the Wash Solution. Do not wash more strips than can be handled within 30 min.
- Pipette 25  $\mu$ L of the CEA Calibrators (CAL 0, 2, 5, 15, 50, 75), controls (c) and patient samples (unknowns-Unk) into the strip wells according to the following scheme:

	1	2	3	4	5	6	7 etc
A	Cal 0	Cal 50	Unk 1				
B	Cal 0	Cal 50	Unk 1				
C	Cal 2	Cal 75	Unk 2				
D	Cal 2	Cal 75	Unk 2				
E	Cal 5	C1	etc.				
F	Cal 5	C1					
G	Cal 15	C2					
H	Cal 15	C2					

- Add 100  $\mu$ L of Antibody Solution to each well using a 100  $\mu$ L precision pipette (or an 8-channel 100  $\mu$ L precision pipette). Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or the surface of the liquid.
- Incubate the frame containing the strips for 1 hour ( $\pm$  5 min) at room temperature (20–25°C) with constant shaking of the plate using a microplate shaker.
- Wash each strip 6 times, using the wash procedure described in Procedural notes item 4.
- Add 100  $\mu$ L of TMB HRP-Substrate to each well using the same pipetting procedure as in item 4. The TMB HRP-Substrate should be added to the wells as quickly as possible and the time between the addition to the first and last well should not exceed 5 min.
- Incubate for 30 min ( $\pm$  5 min) at room temperature with constant shaking. Avoid direct sunlight.
- Immediately read the absorbance at 620 nm in a microplate spectrophotometer.

## Option

If the laboratory does not have access to a microplate spectrophotometer capable of reading at 620 nm, the absorbance can be determined as follows:

Alt. 9. Add 100  $\mu\text{L}$  of Stop Solution. Mix and read the absorbance at 405 nm in a microplate spectrophotometer within 5 minutes after addition of Stop Solution.

## Measurement range

The CanAg CEA EIA measures concentrations between 0.25 and 75  $\mu\text{g/L}$ . If CEA concentrations above the measuring range are to be expected, it is recommended to dilute samples with CEA Calibrator 0 prior to analysis.

## Quality control

CEA Control 1 and 2 may be used for validation of the assay series. Ranges of expected results are indicated on the vial labels. If values outside of the specified range are obtained, a complete check of reagents and reader performance should be made and the analysis repeated. It is recommended that each laboratory in addition prepare its own serum pools at different levels, which can be used as internal controls in order to assure the accuracy of the assay.

## Reference material

The 1<sup>st</sup> International Reference Preparation IRP 73/601 may be used as a reference standard. Values for CEA Calibrators and Controls were assigned against a set of in-house reference standards whose values are traceable to IRP 73/601 using the conversion factor 13.5, i.e. 1  $\mu\text{g/L}$  corresponds to 13.5 IU/L.

## CALCULATION OF RESULTS

If a microplate spectrophotometer reader with built-in data calculation program is used, refer to the manual for the plate reader and create a program using the concentration stated on the labels of each of the CEA Calibrators.

For automatic calculation of CEA results it is recommended to use either of the following methods:

- Cubic spline curve fit method. Calibrator 0 should be included in the curve with the value 0  $\mu\text{g/L}$ .
- Spline smoothed curve fit method. Calibrator 0 should be used as plate blank.
- Interpolation with point-to-point evaluation. Calibrator 0 should be included in the curve with the value 0  $\mu\text{g/L}$ .
- Quadratic curve fit method. Calibrator 0 should be included in the curve with the value 0  $\mu\text{g/L}$ .

**NOTE:** 4-parametric or linear regression should not be used.

For manual evaluation, a calibration curve is constructed by plotting the absorbance (A) values obtained for each CEA calibrator against the corresponding CEA concentration (in  $\mu\text{g/L}$ ), see figure below. The unknown CEA concentrations can then be read from the calibration curve using the mean absorbance value of each patient specimen.

If samples in an initial analysis give CEA levels higher than  $75 \mu\text{g/L}$  the samples should be diluted 1/10 and 1/100 with CEA calibrator 0 and reanalyzed to obtain the accurate CEA concentration.

1 : 10 dilution =  $50 \mu\text{L}$  of specimen +  $450 \mu\text{L}$  of CEA 0  $\mu\text{g/L}$

1 : 100 dilution =  $50 \mu\text{L}$  of 1:10 dilution +  $450 \mu\text{L}$  of CEA 0  $\mu\text{g/L}$


The CEA concentration of the undiluted sample is then obtained as follows:

Dilution 1/10:  $10 \times$  Measured value

Dilution 1/100:  $100 \times$  Measured value

### Example of results

Specimen	Calibrator values	Mean abs value (A)	CEA ( $\mu\text{g/L}$ )
CAL CEA 0	0 $\mu\text{g/L}$	0.050	
CAL CEA 2	2 $\mu\text{g/L}$	0.131	
CAL CEA 5	5 $\mu\text{g/L}$	0.259	
CAL CEA 15	15 $\mu\text{g/L}$	0.657	
CAL CEA 50	50 $\mu\text{g/L}$	1.857	
CAL CEA 75	75 $\mu\text{g/L}$	2.519	
Specimen A		0.220	4.1
Specimen B		1.290	32.3

**Example (do not use this curve or table above to determine actual assay results).**

## LIMITATIONS OF THE PROCEDURE

The level of CEA cannot be used as absolute evidence for the presence or absence of malignant disease and the CEA test should not be used in cancer screening. The results of the test should be interpreted only in conjunction with other investigations and procedures in the diagnosis of disease and the management of patients and the CEA test should not replace any established clinical examination.

Anti-reagent antibodies (human anti-mouse antibody (HAMA) or heterophilic antibodies) in the patient sample may occasionally interfere with the assay, even though specific blocking agents are included in the buffer.

Biotin may interfere with the assay giving false low results. This should be taken into consideration for patients taking dietary supplements or receiving therapy containing high (> 5mg/day) or extremely high (300 mg/day) biotin doses. Peak serum levels have been reported to occur 1-3h post ingestion and a physiological half-life of 8-16h depending on renal function (9,10). In a study by Grimsey et al. (9) a specimen concentration of 30 ng/mL was reached 8h following intake of 10 mg of biotin. For patients on very high doses of biotin it is recommended to stop taking biotin for at least 2 days before blood draw (11).

## EXPECTED VALUES

CanAg CEA was measured in 95 healthy blood donors and in 117 healthy individuals between 60 and 64 years. The lower and upper extremes of the normal range were examined using IFCC recommended non-parametric statistical treatment. The reference interval contains the central 95% fraction of the reference distribution.

Reference limits may accordingly be estimated as the 2.5 % (lower) and 97.5 % (upper) fractiles. These limits cut off a fraction of 2.5 % of the values in each tail of the reference distribution. Non-parametric estimates:

	Mean (µg/L)	SD (µg/L)	Median (µg/L)	Range (µg/L)	Upper reference limit (Central 95 % fraction)
Healthy blood donors n=95	1.3	1.0	1.0	0.5–9.1	3.2 µg/L
Healthy individuals age 60-64, n=117	2.4	1.7	1.9	0.5–8.8	7.4 µg/L

96 % of the healthy subjects had assay values below 5 µg/L.

It is recommended that each laboratory establish their own normal range to account for such local environmental factors as diet, climate, living conditions, patient selection, etc. It should also be borne in mind that the individual patient's own baseline results provides the most important reference point for interpretation of marker results. Smoking may increase CEA levels in healthy individuals.

## PERFORMANCE CHARACTERISTICS

### Precision

Intermediate precision was calculated according to NCCLS guideline EP5-A (7) using four levels of frozen pooled human serum containing added CEA and two different CanAg CEA EIA reagent combinations. Each sample was randomly pipetted (n=2/analysis) and analysed twice each day over 20 days.

Sample	Replicates	Mean (µg/L)	Within-run SD (µg/L)	Within-run CV %	Between-day SD (µg/L)	Between-day CV %
CEA 1	80	2.78	0.07	2.5	0.08	2.7
CEA 2	80	5.97	0.15	2.6	0.11	1.8
CEA 3	80	20.8	0.44	2.1	0.36	1.7
CEA 4	80	57.3	1.57	2.7	0.87	1.5

### Detection limit

The detection limit of the CanAg CEA EIA is  $\leq 0.25$   $\mu\text{g/L}$  defined as the concentration corresponding to the mean of the absorbance values of the CEA calibrator 0 plus 2 standard deviations according to formula:

$$\frac{2 \times \text{SD CAL 0}}{\text{OD CAL 2} - \text{OD CAL 0}} \times 2 \mu\text{g/L}$$

### Recovery

Spiked serum samples were prepared by adding human CEA antigen to normal serum samples. The recovery of the added antigen was in the range 90–115 %.

### Hook effect

When reading absorbance at 405 nm, i.e. using the Optional assay procedure with addition of STOP solution, no hook effect has been noticed for samples containing up to 250 000  $\mu\text{g/L}$ . When absorbance is read at 620 nm, extremely high samples may change the colour of the substrate from blue to greenish. This may lead to a falsely low absorbance that may fall within the calibration curve range and noticed as a hook. Such a hook effect at 620 nm has been noticed for samples containing more than 2000  $\mu\text{g/L}$ .

In order to avoid reporting misleadingly low results due to apparent hook effect when absorbance is read at 620 nm it is recommended to use the Optional assay procedure and determine absorbance at 405 nm in patients analysed for the first time or in patients where very high CEA values may be expected.

### Linearity

Patient samples were serially diluted with CEA Calibrator 0 and analysed. The obtained values were between 90–120 % of the expected values.

### Specificity

The CanAg CEA EIA is based on two mouse monoclonal antibodies, the catching MAb 12-140-10 against Gold epitope IV and the detecting MAb 12-140-1 against Gold epitope V (4, 5). The NCCLS guideline EP7-P (8) was followed to determine possible sources of interference. The following substances and concentrations were tested and found not to interfere with the test.

	Concentration with no significant (± 10 %) interference
Lipemia (Intralipid®)	10 mg/mL
Bilirubin, unconjugated	0.6 mg/mL
Hemoglobin	5 mg/mL

### Biotin interference

A study was conducted to evaluate biotin interference. Low and high serum control samples were spiked to final biotin concentrations of 15, 30, 60, and 200 ng/mL. The mean CEA concentration was determined for each sample and the percent recovery for each biotin concentration was calculated using the formula: Recovery (%) = 100x (Mean CEA concentration w. biotin added/Mean CEA concentration w. diluent only added).

CEA analyte level	Biotin test conc. (ng/mL)	Expected CEA conc. (µg/L)	Observed CEA conc. (µg/L)	Recovery (%)
Low	15	3,72	3,53	95
High	15	39,5	37,8	96
Low	30	3,87	3,62	93
High	30	40,5	37,5	92
Low	60	3,90	3,13	80
High	60	41,3	32,5	79
Low	200	3,76	1,21	32
High	200	40,9	14,9	36

Based on linear regression analysis, the lowest concentration of biotin found to influence test results (≥10%) was 33 ng/mL.

### Method comparison

The CanAg CEA EIA was compared to the Wallac Delfia CEA kit. Seventy-seven human serum samples ranging in values from 0-790 µg/L were measured and linear regression analyses of the results yielded:

$$\text{CanAg CEA} = 0.90 \times \text{Delfia CEA} + 0.53 \quad r = 1.00$$

## CLP (1272/2008) HAZARD CLASSIFICATION

The following warnings and precautions apply to

SUBS	TMB
------	-----

### Hazard pictograms:

**Signal word:** Danger

**Hazard Statement:** Repr. 1B: H360D May damage the unborn child.

**Prevention statement:** P202 Do not handle until all safety precautions have been read and understood.

**Prevention:** P280 Wear protective gloves / protective clothing / eye protection / face protection.

**Precautionary statement response:** P308+P313 IF exposed or concerned get medical advice/attention.

**Precautionary statement disposal:** P501 Dispose of contents / container to an approved hazardous / special waste disposal facility in accordance with local and national regulations.

**Restricted to professional users.**

**Hazardous substances:** 2- Pyrrolidone

### Other hazards

None of the mixtures in the kit contains any substances considered to meet the criteria classifying them as PBT and/or vPvB.

## WARRANTY

The performance data presented here were obtained using the assay procedure indicated. Any change or modification of the procedure not recommended by Fujirebio Diagnostics AB may affect the results, in which event Fujirebio Diagnostics AB disclaims all warranties expressed, implied or statutory including the implied warranty of merchantability and fitness for use.

## REFERENCES

1. Gold P. and Freedman S. O. (1965) Specific carcinoembryonic antigens of the human digestive system. *J. Exp Med* 122: 467–481.
2. Thompson J.A. and Zimmerman W. (1988) The carcinoembryonic antigen gene family: Structure, expression and evolution. *Tumor Biol.*; 9: 63–83.
3. Börner, O.P. (1992) Thesis "Immunoassays for Carcinoembryonic antigen, Specificity and Interferences", ISBN 82–7633–014–2.
4. Hammarström S. et al. (1989) Antigenic sites in carcinoembryonic antigen. *Cancer Res* 49: 4852–4858.
5. Tumor Marker Expert Panel, ASCO (1996) Clinical practice guidelines for the use of tumor markers in breast and colorectal cancer. *J. Clin Oncol.*; 14: 2843–2877.
6. Fleisher M. et al. (2002) Practice guidelines and recommendations for use of tumor markers in the clinic. *National Academy of Clinical Biochemistry* 15: 26-29.
7. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
8. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).
9. Grimsey P et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. *J Pharmacokinet.* 2017; 2(4); 247-256.
10. Jenkins Colon P, Greene D.N. Biotin Interference in Clinical Immunoassays. *JALM* 2018; 2(6); 941-951.
11. Chun KY. Biotin Interference in Diagnostic Tests. *Clinical Chemistry* 2017; 63; 619-620.

---

**CanAg<sup>®</sup>** is a registered trademark of Fujirebio Diagnostics AB

**Fujirebio Diagnostics AB**

**Elof Lindälvs gata 13**

**SE-414 58 Göteborg**

**Sweden**

**Phone + 46 31 85 70 30**

**Fax + 46 31 85 70 40**

**info@fdab.com**

**www.fdab.com**