# CanAg SCC EIA

REF

800-10

IVD

Instructions for use. 2022-06

Read highlighted changes

EN	EXPLANATION OF SYMBOLS
BG	ОБЯСНЕНИЕ НА СИМВОЛИТЕ
CS	VÝZNAM SYMBOLŮ
DA	SYMBOLFORKLARING
DE	ERKLÄRUNG DER SYMBOLE
EL	ΕΠΕΞΗΓΗΣΗ ΤΩΝ ΣΥΜΒΟΛΩΝ
ES	SIGNIFICADO DE LOS SÍMBOLOS
ET	SÜMBOLITE SELGITUS
FR	EXPLICATION DES SYMBOLES
HR	OBJAŠNJENJE SIMBOLA
HU	JELMAGYARÁZAT
IT	SPIEGAZIONE DEI SIMBOLI
LT	SIMBOLIŲ PAAIŠKINIMAI
LV	SIMBOLU SKAIDROJUMS
NL	VERKLARING DER SYMBOLEN
NO	SYMBOLFORKLARING
PL	OBJAŚNIENIE SYMBOLI
PT	EXPLICAÇÃO DOS SÍMBOLOS
RO	SEMNIȚAȚIA SIMBOLURILOR
RU	ОБОЗНАЧЕНИЯ
SV	SYMBOLFÖRKLARING
SK	VÝZNAM SYMBOLOV
SL	RAZLAGA SIMBOLOV
SR	OBJAŠNJENJE SIMBOLA
TR	SEMBOLLERİN AÇIKLAMALARI

Use By/Годно до/Použitelné do/  
Holdbar til/Verwendbar bis/  
Ημερομηνία λήξης/Fecha  
de caducidad/Kölblik kuni/  
Utiliser jusque/Rok valjanosti/  
Felhasználható/Utilizzare entro/  
Sunautoti iki/Izlietot līdz/Houdbaar  
tot/Brukes innen/Użyç przed/  
Prazo de validade/Expirã la/  
Использовать до/Använd före/  
Použite né do/ Uporabno do/  
Upotrebljivo do/Son Kullanna Tarihi

LOT

Batch code/Номер на партида/  
Číslo šarže/Lotnummer/  
Chargenbezeichnung/Αριθμός  
Παρτίδας/Código de lote/Partii  
kood/Code du lot/Kod serije/  
Sarzsám/Codice del lotto/  
Partijas kods/Partijas kods/Lot  
nummer/Partikode/Kod partii/  
Código do lote/Număr de lot/  
Номер лота/Lotnummer/Číslo  
šarže/Številka serije/Kod partije/  
Parti Kodu

Date of manufacture/Дата на производство/Datum výroby/  
Produktionsdato/Herstellungsdatum/  
Ημερομηνία παραγωγής/Fecha de fabricación/Valmistamise kuupäev/  
Date de fabrication/Datum proizvodnje/  
Gyártási idő/Data di produzione/  
Pagaminimo data/Ražošanas datums/  
Productiedatum/Fremstillingsdato/  
Data produkcji/Data de fabrico/Data fabricației/Дата производства/  
Tillverkningsdatum/Dátum výroby/Datum izdelave/Datum proizvodnje/Úretim tarihi

Temperature limitation/  
Температурни граници/  
Теплотни омеzeи/  
Temperaturbegrænsning/  
Temperaturbegrenzung/  
Περιορισμοί θερμοκρασίας/  
Limites de temperatura/  
Temperatuuri piirang/  
Limite de température/  
Temperaturno ograničenje/  
Hőmérsékletre vonatkozó korlátozás/  
Limiti di temperatura/  
Temperatūriniai apribojimai/  
Temperatūras ierobežojums/  
Temperaturbepërking/  
Temperaturbegrensinger/  
Temperaturey graniczne/  
Limite de temperatura/  
Limite de temperatură/  
Температурный режим/  
Temperaturbegrænsning/  
Теплотне обмеzenie  
Omejitve temperature/  
Temperaturno ograničenje/  
Sıcaklık sınırlaması/

## IVD

In Vitro Diagnostic Medical Device/  
Медицински уред за диагностика  
ин vitro/Diagnostický zdravotnícký  
prostředek in vitro/Medicinsk udstyr til  
in vitro-diagnostik/In-vitro-Diagnostikum/  
Ιατροτεχνολογικό προϊόν για διάγνωση  
In Vitro/Dispositivo médico para  
diagnóstico in vitro/In vitro diagnostiline  
meditsiiniseade/Dispositif médical de  
diagnostic in vitro/Diagnostički medicinski  
uređaj In Vitro/In vitro orvosdiagnostikai  
eszköz/Dispositivo medico per test  
diagnostici in vitro/In Vitro Diagnostinė  
Medicinos Priemonė/Medicínska ierice  
in vitro diagnostikai/In vitro-diagnostisch  
medisch instrument/In vitro diagnostisk  
medisinsk utstyr/Wyrób medyczny do  
diagnostyki in vitro/Dispositivo Médico  
de Diagnóstico In Vitro/Dispozitiv medical  
pentru diagnostic in vitro/Только для  
диагностики In Vitro/Endast för in  
vitro-diagnostik/ Zdravotnicka pomôcka na  
diagnostiku in vitro/In vitro diagnostični  
pripomoček/Diagnostički medicinski  
uređaj In Vitro/<96> testleri için yeterlilik  
içerir

Contains sufficient for <96> tests/Съдържа  
достатъчно количество за тестове  
<96>/Lze použít pro <96> testů/Ineholder  
tilstrækkeligt/Inhalt ausreichend für <96>  
Prüfungen/Περισχόμενο επαρκές για  
«96» εξετάσεις/Contenido suficiente para  
<96> ensayos/Kogusest piisab <96> testi  
lääbiviimiseks/Contenu suffisant pour «96»  
tests/Sadržaj dovoljno za <96> testova/A  
doboz tartalma <96> vizsgálat elvégzéséhez  
elegendő/Contenuto sufficiente per «96»  
saggi/Turiny's skirtas atlikti <96> tyrimus/  
Satur's pietiekams <96> testiem/Inhoud  
voldoende voor «96» testen/til «96» test/  
Tilstrækkelig innhold for <96> prøver/  
Wystarczy na wykonanie <96> testów/  
Conteúdo suficiente para «96» ensaios/  
Conținut suficient pentru 96 de teste/  
Содержит достаточные количества для  
«96» определений/Innehåller tillräckligt  
till «96» antal tester/Obsah postačuje na  
tento počet testov: <96>/Vsebinsa zadostuje  
za <96> testov/Sadržina dovoljna za <96>  
testova/<96> testleri için yeterlilik içerir

## REF

Catalogue number/Каталожен номер/  
Katalogové číslo/Katalognummer/  
Bestellnummer/Αριθμός καταλόγου/  
Número de catálogo/Katalogoi number/  
Numéro de catalogue/Kataloški broj/  
Katalógusszám/Numero di catalogo/  
Katalogo numeris/Numurs katalogā/  
Catalogusnummer/Katalognummer/  
Numer katalogowy/Número do catálogo/  
Număr de catalog/Номер по каталогу/  
Produktnummer/Katalógové číslo/  
Kataloška številka/Kataloški broj/  
Katalog numarası

Consult Instructions for Use/  
Прочетете инструкцията за  
употреба/Konzultujte s návodem  
k použití/Se brugsanvisning/Siehe  
Gebrauchsanweisung/Συμβουλευτείτε  
της Οδηγίας σχετικά με τη χρήση/  
Consulte las instrucciones de uso/  
Vt kasutusjuhendit/Consulter le mode  
d'emploi/Pročítajte upute za uporabu/  
Olvassa el a használati utasítást/  
Consultare le istruzioni per l'uso/Dél  
naudojimo žiūrėkite instrukcijas/Izlasiet  
lietošanas instrukciju/Raadpleeg de  
instructies voor gebruik/Les instruksene  
for bruk/Sprawdzić w instrukcji użycia/  
Consulte as Instruções de Utilização/  
Consultați instrucțiunile de utilizare/  
Обратитесь к инструкции по  
применению/Se bruksanvisning/  
Prečítajte si návod na používanie/  
Pročítajte uputstvo za upotrebu/  
Kullanım Talimatlarını Bakınız

## CONT

Contents of kit/Съдържание на набора/  
Obsah soupravy/Kittets indhold/Inhalt  
des Kits/Περιεχόμενα του κιτ/Contenido  
del kit/Komplekt sisaldab/Contenu du  
kit/Sadržaj opreme/A készlet tartalma/  
Contenuto del kit/Rinkinio turinys/  
Komplekta saturs/Inhoud van de set/  
Settets innhold/Zawartość zestawu/  
Conteúdo do kit/Conținutul setului/  
Компоненты набора/Kit innehåll/  
Obsah súpravy/Vsebina kompleta/Sadržaj  
opreme/Kitin içindekiler

Biological risks/Биологическа  
опасност/Biológická rizika/Biologisk  
fare/Biologische Gefahren/Βιολογικοί  
κίνδυνοι/Riesgos biológicos/  
Biolooigilised ohud/Risques biologiques/  
Biolóskli rizici/Biológiai kockázatok/Rischi  
biologici/Biologinis pavojus/Biológiskais  
risks/Biologische risico's/Biologische  
risikoer/Zagroženie biologiczne/Riscos  
biológicos/ Biologisk risk/Pericole  
biologice/Биологическая опасность/  
Biologicky rizikové/Biológické riziká/  
Biolóskli rizici/Biyolojik riskler

## ORIG HUM

Human/C човешки производ/Ľidské/  
Human/Human/δείγματα αναφοράς/  
Humano/Inimāritolu/Humaine/Ljudskog  
porjekla/Humán/Origine Umana/  
Žmogaus kilmės/Cilvēku izcelsmes/  
Human/Menneske/Ludzka/Humano/  
Origine umână/Человеческого  
происхождения/Human/Ludské/  
Humanega izvora/Ljudskog porekla/Ľnsan

## ORIG MOU

From mouse/C миши производ/Myši/  
Fra mus/Maus/από πογτίκι/de ratón/  
Hiirtelt/De souris/Mišijeg porjekla/  
Egérböli/Murino/Pelés kilmés/No peles/  
Van muizen/Fra mus/Mysia/Do rato/De  
la șoareci/Мышиного происхождения/  
Från mus/Myšije/Mišjega izvora/Mišijeg  
porekla/Fareden

## ORIG BOV

Bovine/C говежди производ/  
Hovēži/Bovin/Rind/από βοοειδή/  
Bovino/Veistelt/Bovine/Rogate stoke/  
Szarvasmarha/Bovino/Jaučio/No  
liellopa/Bovien/Bovin/Wolowy/Bovino/  
Origine bovină/крупного рогатого  
скота/Från ko/Hovädzie/Rogveja  
izvora/Rogate krupne stoke/Bovin

Reconstitute with/Пазтваряне с/  
Rozfeďte pomocí/Rekonstitueres med/  
Rekonstituieren mit/Ανασύσταση με/  
Reconstituir con/Lahjendamine/  
Reconstituer avec/Rekonstituiraite s/  
Feloldáshoz/Ricostituire con/Atkurti,  
ištirpdant su/Atšķaidīt ar/Reconstituite  
met/Rekonstitueres med/Odtworzyć  
za pomocą/Reconstituir com/A  
se reconstitui cu/Пастворить в/  
Rekonstituera med/Rozriedte pomocou/  
Rekonstituiraite z/s/Ponovno formiranje  
sa/Yeniden oluşturalur

Manufacturer/Производитель/Výrobce/  
Producent/Hersteller/Κατασκευαστής/  
Fabricante/Tootja/Fabricant/Proizvođač/  
Gyártó/Fabbricante/Gamintojas/  
Ražotājs/Fabrikant/Produsent/  
Producent/Fabricante/Producător/  
Производитель/Tilverkare/ Výrobca/  
Izdelovalec/Proizvođač/Üretici

# INSTRUCTIONS FOR USE

EN

## INSTRUCTIONS FOR USE

Please visit our website [www.fdi.com/ifu](http://www.fdi.com/ifu) to obtain the Instructions For Use (IFU) in additional languages.

To ensure that you download the correct IFU for your kit lot, please select the revision corresponding to the issue date printed on the front page of the IFU provided with this kit.

Please follow the IFU carefully. Instructions for safe handling are found in the WARNINGS AND PRECAUTIONS section. Material Safety Data Sheets (MSDS) are available on our website [www.fdi.com](http://www.fdi.com). If you do not have access to the internet, please contact your local distributor, or Fujirebio Diagnostics AB for assistance.

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# CanAg SCC EIA

Instructions for use

Enzyme immunometric assay kit  
For 96 determinations

## IMPORTANT USER INFORMATION

SCC antigen is present in skin, sweat and saliva, and is easily distributed in aerosols (e.g. as a result of sneezing). In order to avoid false elevated values due to contamination, gloves should be used after opening the kit box and throughout the test procedure when handling reagent vials, microplate, pipette tips etc. In addition, all elevated values should be confirmed by repeat testing.

## INTENDED USE

The CanAg SCC EIA kit is intended for the quantitative determination of squamous cell carcinoma (SCC) antigen in serum as an aid in the management of patients with squamous cell carcinoma.

## SUMMARY AND EXPLANATION OF THE ASSAY

Squamous cell carcinoma antigen (SCC ag) is a group of glycoproteins with molecular weight ~45 kDa, belonging to the family of serine/cysteine -protease inhibitors (1). The protein was originally isolated by Kato and co-workers from human squamous cell carcinoma tissue and shown to consist of at least 10 subfractions differing in isoelectric point (2). More recent studies have shown that SCC antigen is composed of two distinct but highly homologous gene products, SCCA1 and SCCA2 with different inhibitor specificities (3).

SCC antigen is a serological marker of squamous cell carcinomas of the uterine cervix, vulva, lung, head & neck, and oesophagus (4-6). In squamous cell carcinoma of the uterine cervix, pre-treatment serum SCC ag may be used as an early stage prognostic factor (7) and the use of pre-treatment SCC ag have been suggested in order to select high-risk patients for adjuvant therapy (4). Further, for patients with elevated levels of SCC ag before start of treatment, the profile of SCC ag correlates with the response to radio- and chemo-therapy and measurement of SCC ag may thus be used to monitor the effect of therapy and for early detection of recurrent disease (4).

## PRINCIPLE OF THE TEST

The CanAg SCCEIA is a solid-phase, non-competitive immunoassay based upon the direct sandwich technique. Calibrators and patient samples are incubated together with biotinylated Anti-SCC monoclonal antibody and horseradish peroxidase (HRP) labelled Anti-SCC monoclonal antibody in Streptavidin coated microstrips. After washing, buffered Substrate/ Chromogen reagent (hydrogen peroxide and 3, 3', 5, 5' tetra-methylbenzidine) is added to each well and the enzyme reaction is allowed to proceed. During the enzyme reaction a blue colour will develop if antigen is present. The intensity of the colour is proportional to the amount of SCC present in the samples.

The colour intensity is determined in a microplate spectrophotometer at 620 nm (or optionally at 405 nm after addition of Stop Solution). Calibration curves are constructed for each assay by plotting absorbance value versus the concentration for each calibrator. The SCC concentrations of patient samples are then read from the calibration curve.

## REAGENTS

- Each CanAg SCC EIA kit contains reagents for 96 tests.
- The expiry date of the kit is stated on the label on the outside of the kit box.
- Do not use the kit beyond the expiry date.
- Do not mix reagents from different kit lots.
- Store the kit at 2–8 °C. Do not freeze.
- Opened reagents are stable according to the table below provided they are not contaminated, stored in resealed original containers and handled as prescribed. Return to 2–8 °C immediately after use.

Component	Quantity	Storage and stability after first opening
<b>MICROPLA</b>		
<b>Microplate</b>	1 Plate	2–8 °C until expiry date stated on the plate

12 x 8 wells coated with Streptavidin. After opening, immediately return unused strips to the aluminium pouch, containing desiccant. Reseal carefully to keep dry.

Component	Quantity	Storage and stability after first opening			
<b>SCC Calibrators</b>	5 vials, lyophilized	4 weeks at 2–8 °C 3 months at -20 °C			
<table border="1"><tr><td>CAL</td><td>SCC</td><td>A</td></tr></table>	CAL	SCC	A	1 x 0.75 mL	
CAL	SCC	A			
<table border="1"><tr><td>CAL</td><td>SCC</td><td>B</td></tr></table>	CAL	SCC	B	1 x 0.75 mL	
CAL	SCC	B			
<table border="1"><tr><td>CAL</td><td>SCC</td><td>C</td></tr></table>	CAL	SCC	C	1 x 0.75 mL	
CAL	SCC	C			
<table border="1"><tr><td>CAL</td><td>SCC</td><td>D</td></tr></table>	CAL	SCC	D	1 x 0.75 mL	
CAL	SCC	D			
<table border="1"><tr><td>CAL</td><td>SCC</td><td>E</td></tr></table>	CAL	SCC	E	1 x 0.75 mL	
CAL	SCC	E			

The lyophilised calibrators contain human SCC in a Tris-HCl buffered salt solution containing bovine serum albumin, excipient, an inert yellow dye and 0.01% methyl-isothiazolone (MIT) as preservative. To be reconstituted with water before use. **NOTE:** The exact SCC concentration is lot specific and is indicated on the label of each vial.

<b>BIOTIN</b>	<b>Anti-SCC</b>
---------------	-----------------

<b>Biotin Anti-SCC</b>	1 x 15 mL	2–8 °C until expiry date stated on the vial
------------------------	-----------	---

Biotin Anti-SCC monoclonal antibody from mouse, approximately 1 µg/mL. Contains phosphate buffered saline (pH 7.2), bovine serum albumin, bovine immunoglobulin, blocking agents, detergent, an inert blue dye and 0.01% methyl-isothiazolone (MIT) as preservative. To be mixed with Tracer, HRP Anti-SCC before use.

<b>CONJ</b>	<b>Anti-SCC</b>
-------------	-----------------

<b>Tracer, HRP Anti-SCC</b>	1 x 0.75 mL	2–8 °C until expiry date stated on the vial
-----------------------------	-------------	---

Stock solution of HRP Anti-SCC monoclonal antibody from mouse, approximately 40 µg/mL. Contains preservatives. To be mixed with Biotin Anti-SCC before use.

Component	Quantity	Storage and stability after first opening
-----------	----------	---

<b>SUBS</b>	<b>TMB</b>
-------------	------------

<b>TMB HRP-Substrate</b>	1 x 12 mL	2–8 °C until expiry date stated on the vial
--------------------------	-----------	---

Contains buffered hydrogen peroxide and 3, 3', 5, 5' tetramethyl-benzidine (TMB). Ready for use.

<b>STOP</b>
-------------

<b>STOP Solution</b>	1 x 15 mL	2–8 °C until expiry date stated on the vial
----------------------	-----------	---

Contains 0.12 M hydrochloric acid. Ready for use.

<b>WASHBUF</b>	<b>25X</b>
----------------	------------

<b>Wash Concentrate</b>	1 x 50 mL	2–8 °C until expiry date stated on the bottle
-------------------------	-----------	---

A Tris-HCl buffered salt solution with Tween 20. Contains Germall II as preservative. To be diluted with water 25 times before use.

### Indications of instability

The TMB HRP-Substrate should be colourless or slightly bluish. A blue colour indicates that the reagent has been contaminated and should be discarded.

### WARNINGS AND PRECAUTIONS

#### For in vitro diagnostic use.

- For professional use only.
- Please refer to the US Department of Health and Human Services (Bethesda, Md., US) publication No. (CDC) 88-8395 on laboratory safety or any other local or national regulation.
- Handle all patient specimens as potentially infectious.
- Follow local guidelines for disposal of all waste material.

### CLP (1272/2008) HAZARD CLASSIFICATION

Information about CLP (1272/2008) HAZARD CLASSIFICATION can be found at the end of this document.

# Protocol Sheet

## CanAg SCC EIA REF 800-10

Prepare the components directly before use. Use wash and incubation conditions according to the Instructions.

Step	Vial/Plate	Procedure																					
1. Prepare SCC Calibrators	<table border="1"><tr><td>CAL</td><td>SCC</td></tr></table> A, B, C, D, E	CAL	SCC	Add 0.75 mL of distilled water to each vial and mix gently. Allow to stand for at least 15 minutes. <b>NOTE:</b> The exact concentration of each calibrator is stated on the label. This value of the calibrators should be used for calculations.																			
CAL	SCC																						
2. Prepare Wash Solution	<table border="1"><tr><td>WASHBUF</td><td>25X</td></tr></table>	WASHBUF	25X	Dilute 50 mL of Wash Concentrate with 1200 mL of distilled water or deionized water.																			
WASHBUF	25X																						
3. Prepare Antibody Solution	<table border="1"><tr><td>CONJ</td><td>Anti-SCC</td></tr><tr><td>BIOTIN</td><td>Anti-SCC</td></tr></table>	CONJ	Anti-SCC	BIOTIN	Anti-SCC	Mix 50 $\mu$ L of Tracer, HRP Anti-SCC with 1 mL of Biotin Anti-SCC per strip:																	
CONJ	Anti-SCC																						
BIOTIN	Anti-SCC																						
	<table border="1"><thead><tr><th>No. of Strips</th><th>HRP Anti-SCC (<math>\mu</math>L)</th><th>Biotin Anti-SCC (mL)</th></tr></thead><tbody><tr><td>1</td><td>50</td><td>1</td></tr><tr><td>2</td><td>100</td><td>2</td></tr><tr><td>3</td><td>150</td><td>3</td></tr><tr><td>4</td><td>200</td><td>4</td></tr><tr><td>5</td><td>250</td><td>5</td></tr><tr><td>6</td><td>300</td><td>6</td></tr></tbody></table>	No. of Strips	HRP Anti-SCC ( $\mu$ L)	Biotin Anti-SCC (mL)	1	50	1	2	100	2	3	150	3	4	200	4	5	250	5	6	300	6	
No. of Strips	HRP Anti-SCC ( $\mu$ L)	Biotin Anti-SCC (mL)																					
1	50	1																					
2	100	2																					
3	150	3																					
4	200	4																					
5	250	5																					
6	300	6																					

				7	350	7
				8	400	8
				9	450	9
				10	500	10
				11	550	11
				12	600	12
4. Wash	<b>MICROPLA</b>			Wash each well once with Wash Solution		
5. Add calibrators and samples	<b>CAL</b>	<b>SCC</b>		25 µL in each well		
	A, B, C, D, E					
6. Add Antibody Solution	<b>ANTIBODY SOLUTION</b>			100 µL in each well		
7. Incubate	<b>MICROPLA</b>			1 hour shaking at room temperature		
8. Wash	<b>MICROPLA</b>			Wash each well six times with Wash Solution		
9. Add TMB HRP-Substrate	<b>SUBS</b>	<b>TMB</b>		100 µL in each well		
10. Incubate	<b>MICROPLA</b>			30 min shaking at room temperature		
11. Read absorbance	<b>MICROPLA</b>			620 nm		
Alt.11 Add Stop Solution	<b>STOP</b>			100 µL in each well		
Alt.12 Incubate	<b>MICROPLA</b>			1 min shaking at room temperature		
Alt.13 Read absorbance	<b>MICROPLA</b>			Read at 405 nm within 5 min		

## SPECIMEN COLLECTION AND HANDLING

The CanAg SCC EIA is intended for use with serum. Collect blood by venipuncture and separate the serum according to common procedures. Samples can be stored at 2–8 °C for 1 day. For longer periods it is recommended to store the samples at –70 °C or below. Avoid repeated freezing and thawing of the samples. Allow frozen samples to thaw slowly, preferably at 2–8 °C over night and then bring the samples to room temperature before analysis.

## PROCEDURE

### Materials required but not supplied with the kit

**1. Microplate shaker**

Shaking should be medium to vigorous, approximately 700-1100 oscillations/min.

**2. Microplate wash device**

Automatic plate washer capable of performing 1 and 6 washing cycles with a minimal fill volume of 350 µL/well/washcycle.

An 8-channel pipette with disposable plastic tips for delivery of 350 µL is recommended if an automatic microplate washer is not used.

**3. Microplate spectrophotometer**

With a wavelength of 620 nm and/or 405 nm and an absorbance range of 0 to 3.0.

**4. Precision pipettes**

With disposable plastic tips to deliver microlitre and millilitre volumes. An 8-channel pipette or dispenser pipette with disposable plastic tips for delivery of 100 µL is useful but not essential.

**5. Distilled or deionized water**

For reconstitution of SCC Calibrators and for preparation of Wash Solution.

### Procedural notes

1. A thorough understanding of this package insert is necessary to ensure proper use of the CanAg SCC EIA kit. The reagents supplied with the kit are intended for use as an integral unit. Do not mix identical reagents from kits having different lot numbers. Do not use the kit reagents after the expiry date printed on the outside of the kit box.
2. Reagents should be allowed to reach room temperature (20–25 °C) prior to use. The assay should only be performed at temperatures between 20–25 °C to

obtain accurate results. Frozen specimens should be brought to room temperature slowly and must be gently but thoroughly mixed after thawing.

3. Before starting to pipette calibrators and patient specimens it is advisable to mark the strips to be able to clearly identify the samples during and after the assay.
4. The requirement for efficient and thorough washing for separation of bound and unbound antigen and reagents from the solid-phase bound antibody-antigen complexes is one of the most important steps in an EIA. In order to ensure efficient washing make sure that all wells are completely filled to the top edge with wash solution during each wash cycle, that wash solution is dispensed at a good flow rate, that the aspiration of the wells between and after the wash cycles is complete and that the wells are empty. If there is liquid left, invert the plate and tap it carefully against absorbent paper.
  - Automatic strip washer: Follow the manufacturer's instructions for cleaning and maintenance diligently and wash the required number of wash cycles prior to and after each incubation step. It's highly recommended to use *strip* process mode and *overflow* wash mode with a dispensing volume of 800  $\mu$ L. The aspiration/wash device should not be left standing with the Wash Solution for long periods, as the needles may get clogged resulting in poor liquid delivery and aspiration.
5. The TMB HRP-Substrate is very sensitive for contamination. For optimal stability of the TMB HRP-Substrate, pour the required amount from the vial to a carefully cleaned reservoir or preferably a disposable plastic tray to avoid contamination of the reagent. Be sure to use clean disposable plastic pipette tips (or respenser pipette tip).
6. Be sure to use clean disposable plastic pipette tips and a proper pipetting technique when handling samples and reagents. Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or surface of the liquid. A proper pipetting technique is of particular importance when handling the TMB HRP-Substrate solution.

Preparation of reagents	Stability of prepared reagent
-------------------------	-------------------------------

<b>SCC Calibrators</b>	4 weeks at 2–8 °C 3 months at –20 °C or below
------------------------	--

Add exactly 0.75 mL of distilled water to each vial and mix gently. Allow to stand for at least 15 minutes to reconstitute. **NOTE:** The concentration of the calibrators is stated on the labels and should be used for calculation of results.

<b>Wash Solution</b>	2 weeks at 2–25 °C in a sealed container
----------------------	--

Pour the 50 mL Wash Concentrate into a clean container and dilute 25- fold by adding 1200 mL of distilled or deionized water to give a buffered Wash Solution.

<b>Antibody Solution</b>	3 weeks at 2–8 °C
--------------------------	-------------------

Prepare the required quantity of Antibody Solution by mixing 50 µL of Tracer, HRP Anti-SCC with 1 mL of Biotin Anti-SCC per strip (see table below and the Protocol Sheet):

No. of Strips	Tracer, HRP Anti-SCC (µL)	Biotin Anti-SCC (mL)
1	50	1
2	100	2
3	150	3
4	200	4
5	250	5
6	300	6
7	350	7
8	400	8
9	450	9
10	500	10
11	550	11
12	600	12

Be sure to use a clean plastic or glass bottle for preparation of the Antibody Solution.

**Alternative:** Pour the content of the Tracer, HRP Anti-SCC into the vial of Biotin Anti-SCC and mix gently. Make sure that all of the Tracer, HRP Anti-SCC is transferred to the vial of Biotin Anti-SCC.

**NOTE:** The Antibody Solution is stable for 3 weeks at 2–8 °C. Do not prepare more Antibody Solution than will be used within this period and make sure that it is stored properly.

### Assay procedure

Perform each determination in duplicate for calibrators and patient samples. A calibration curve should be run with each assay. All reagents and samples must be brought to room temperature (20–25 °C) before use.

1. Start to prepare SCC Calibrators, Wash Solution and Antibody Solution. It is important to use clean containers. Follow the instructions carefully.
2. Transfer the required number of microplate strips to a strip frame. (Immediately return the remaining strips to the aluminium pouch containing a desiccant and reseal carefully). Wash each strip once with the Wash Solution. Do not wash more strips than can be handled within 30 min.
3. Pipette 25 µL of the SCC Calibrators (CAL A, B, C, D, E) and patient samples (unknowns-Unk) into the strip wells according to the following scheme:

	1	2	3	4	5	6	7 etc
A	Cal A	Cal E					
B	Cal A	Cal E					
C	Cal B	Unk1					
D	Cal B	Unk1					
E	Cal C	Unk2					
F	Cal C	Unk2					
G	Cal D	Etc.					
H	Cal D						

4. Add 100 µL of Antibody Solution to each well using a 100 µL precision pipette (or an 8-channel 100 µL precision pipette). Avoid carry-over by holding the pipette tip slightly above the top of the well and avoid touching the plastic strip or the surface of the liquid.

5. Incubate the frame containing the strips for 1 hour ( $\pm 5$  min) at room temperature (20–25 °C) with constant shaking of the plate using a microplate shaker.
6. Wash each strip 6 times, using the wash procedure described in Procedural notes item 4.
7. Add 100  $\mu$ L of TMB HRP-Substrate to each well using the same pipetting procedure as in item 4. The TMB HRP-Substrate should be added to the wells as quickly as possible and the time between the addition to the first and last well should not exceed 5 min.
8. Incubate for 30 min ( $\pm 5$  min) at room temperature with constant shaking. Avoid direct sunlight.
9. Immediately read the absorbance at 620 nm in a microplate spectrophotometer.

### Option

If the laboratory does not have access to a microplate spectrophotometer capable of reading at 620 nm, the absorbance can be determined as follows:

**Alt. 9.** Add 100  $\mu$ L of Stop Solution. Mix and read absorbance at 405 nm in a microplate spectrophotometer within 5 min after addition of Stop Solution.

### Measurement range

The CanAg SCC EIA measures concentrations between 0.3 and 50  $\mu$ g/L. If SCC concentrations above the measuring range are to be expected, it is recommended to dilute samples with normal human serum prior to analysis. **NOTE:** The serum used for dilution should also be measured in order to determine the endogenous SCC concentration (see “Calculation of results”).

### Quality control

CanChek Tumor Marker Control Sera Levels 1 and 2 (available separately, REF 107-20) are recommended for validation of the assay series. If values outside of the specified range are obtained, a complete check of reagents and reader performance should be made and the analysis repeated.

### Reference material

Since no common reference material is available for SCC antigen, CanAg SCC Calibrator values are assigned against a set of in-house reference standards.

## CALCULATION OF RESULTS

If a microplate spectrophotometer reader with built-in data calculation program is used, refer to the manual for the plate reader and create a program using the concentration stated on the labels of each of the SCC Calibrators.

For automatic calculation of SCC results it is recommended to use either of the following methods:

- Cubic spline curve fit method. Calibrator 0 should be included in the curve with the value 0 µg/L.
- Spline smoothed curve fit method. Calibrator 0 should be used as plate blank.
- Interpolation with point-to-point evaluation. Calibrator 0 should be included in the curve with the value 0 µg/L.
- Quadratic curve fit method. Calibrator 0 should be included in the curve with the value 0 µg/L.

**Note:** 4-parametric or linear regression should not be used.

For manual evaluation, a calibration curve is constructed by plotting the absorbance (A) values obtained for each SCC calibrator against the corresponding SCC concentration (in µg/L), see figure below. The unknown SCC concentrations can then be read from the calibration curve using the mean absorbance value of each patient specimen.


If samples in an initial analysis give SCC levels higher than 50 µg/L the samples should be diluted 1/10 with normal human serum and reanalyzed to obtain the accurate SCC concentration. **NOTE:** The sample used for dilution should also be measured in order to determine the endogenous SCC concentration.

The SCC concentration of the undiluted sample is calculated as:

$$\text{Dilution 1/10: } 10 \times ([\text{SCC}]_{\text{Diluted sample}} - (0.9 \times [\text{SCC}]_{\text{Normal serum}}))$$

## Example of results

Specimen	Calibrator values	Mean abs value (A)	SCC ( $\mu\text{g/L}$ )	
CAL	SCC	A	0 $\mu\text{g/L}$	0.043
CAL	SCC	B	1 $\mu\text{g/L}$	0.119
CAL	SCC	C	5 $\mu\text{g/L}$	0.437
CAL	SCC	D	24 $\mu\text{g/L}$	1.715
CAL	SCC	E	50 $\mu\text{g/L}$	2.818
Specimen A		0.245	2.6	
Specimen B		1.363	18.3	

Example (do not use this curve or table above to determine actual assay results).

## LIMITATIONS OF THE PROCEDURE

SCC antigen is present in normal squamous cell epithelia and elevated SCC antigen may be found in skin disorders involving hyperkeratinization eg. psoriasis and eczema. Elevated levels also occur in benign conditions such as inflammatory lung disease and liver or renal insufficiency (4, 9).

The level of SCC thus cannot be used as absolute evidence for the presence or absence of malignant disease, and the SCC test should not be used in cancer screening. The results of the test should be interpreted only in conjunction with other investigations and procedures in the diagnosis of disease and the management of patients, and the SCC test should not replace any established clinical examination.

SCC antigen is present in skin, sweat and saliva, and is easily distributed in aerosols (e.g. as a result of sneezing). In order to avoid false elevated values due to contamination, gloves should be used throughout the test procedure when handling reagent vials, microplate, pipette tips etc. In addition, all elevated values should be confirmed by repeat testing.

Anti-reagent antibodies (human anti-mouse antibody (HAMA) or heterophilic antibodies) in the patient sample may occasionally interfere with the assay, even though specific blocking agents are included in the buffer.

Biotin may interfere with the assay giving false low results. This should be taken into consideration for patients taking dietary supplements or receiving therapy containing high (> 5mg/day) or extremely high (300 mg/day) biotin doses. Peak serum levels have been reported to occur 1-3h post ingestion and a physiological half-life of 8-16h depending on renal function (13,14). In a study by Grimsey et al. (13) a specimen concentration of 30 ng/mL was reached 8h following intake of 10 mg of biotin. For patients on very high doses of biotin it is recommended to stop taking biotin for at least 2 days before blood draw (15).

## EXPECTED VALUES

CanAg SCC EIA was used to measure SCC antigen in 175 healthy blood donors. The lower and upper extremes of the normal range were examined using IFCC recommended non-parametric statistical treatment. The reference interval contains the central 95% fraction of the reference distribution. The upper reference limit was accordingly estimated as the 97.5% upper fractile.

	Mean ( $\mu\text{g/L}$ )	SD ( $\mu\text{g/L}$ )	Median ( $\mu\text{g/L}$ )	Range ( $\mu\text{g/L}$ )	Upper reference limit ( $\mu\text{g/L}$ )
Healthy blood donors n=175	0.58	0.24	0.54	0.16–1.5	1.2 $\mu\text{g/L}$

It is recommended that each laboratory establish their own normal range to account for such local environmental factors as diet, climate, living conditions, patient selection, etc. It should also be borne in mind that the individual patient's own baseline results provides the most important reference point for interpretation of marker results.

## PERFORMANCE CHARACTERISTICS

### Precision

Total precision was calculated according to NCCLS guideline EP5-A (10) using four levels of frozen pooled human serum containing added human SCC antigen and 18 different CanAg SCC EIA reagent combinations. Each sample was randomly pipetted (n=2/analysis) and analysed twice each day over 20 days.

Sample	Replicates	Mean ( $\mu\text{g/L}$ )	Within-run SD ( $\mu\text{g/L}$ )	Within-run CV %	Between-day SD ( $\mu\text{g/L}$ )	Between-day CV %
SCC 1	80	2.62	0.05	1.9	0.04	1.3
SCC 2	80	7.77	0.16	2.0	0.15	1.9
SCC 3	80	17.7	0.34	1.9	0.20	1.1
SCC 4	80	30.2	0.71	2.4	0.38	1.3

### Detection limit

The detection limit of the CanAg SCC EIA is  $\leq 0.3 \mu\text{g/L}$  defined as the concentration corresponding to the mean of the absorbance values of the SCC calibrator A plus 2 standard deviations according to formula:

$$\frac{2 \times \text{SD CAL A}}{\text{OD CAL B} - \text{OD CAL A}} \times [\text{CAL B}] \mu\text{g/L}$$

### Recovery

Spiked serum samples were prepared by adding human SCC antigen to normal serum samples. The recovery of the added antigen was in the range 90–110 %.

### Hook effect

No hook effect has been noticed with samples up to 50 000  $\mu\text{g/L}$ . **NOTE:** In very high samples the colour of the substrate will change from blue to greenish (and eventually yellow in extremely high samples). This will lead to a falsely low absorbance at 620 nm, and in extreme cases the absorbance may fall within the calibration curve range and noticed as a hook.

### Linearity

Patient samples were serially diluted with normal human serum and analysed. The obtained values were within 90–110% of the expected values.

### Specificity

The CanAg SCC EIA is based on two mouse monoclonal antibodies, the catching MAb SCC 140 and the detecting MAb SCC 107 (11). The NCCLS guideline EP7-P (12) was followed to determine possible sources of interference. The following substances and concentrations were tested and found not to interfere with the test.

	Concentration with no significant ( $\pm 10\%$ ) interference
Lipemia (Intralipid®)	10 mg/mL
Bilirubin, unconjugated	0.6 mg/mL
Hemoglobin	5 mg/mL

## Biotin interference

A study was conducted to evaluate biotin interference. Low and high serum control samples were spiked to final biotin concentrations of 15, 30, 60, 200, and 600 ng/mL. The mean SCC concentration was determined for each sample and the percent recovery for each biotin concentration was calculated using the formula: Recovery (%) = 100x (Mean SCC concentration w. biotin added/Mean SCC concentration w. diluent only added).

SCC analyte level	Biotin test conc. (ng/mL)	Expected SCC conc. (µg/L)	Observed SCC conc. (µg/L)	Recovery (%)
Low	15	1,85	1,75	95
High 1	15	17,4	16,8	97
High 2	15	10,6	10,2	96
Low	30	1,85	1,76	95
High 1	30	17,9	16,9	94
High 2	30	10,5	9,91	95
Low	60	1,82	1,63	89
High 1	60	18,0	15,7	87
High 2	60	10,4	9,3	89
Low	200	1,83	0,40	22
High 1	200	17,6	4,00	23
High 2	200	10,6	3,03	29
High 1	600	17,8	1,16	6
High 2	600	11,1	0,88	8

Based on linear regression analysis, the lowest concentration of biotin found to influence test results ( $\geq 10\%$ ) was 41 ng/mL.

## Method comparison

The CanAg SCC EIA was compared to the Imx SCC MEIA.

For 72 human samples ranging in values from 0–4 µg/L, linear regression analyses of the results yielded:

$$\text{CanAg SCC} = 1.02 \times \text{Imx SCC} + 0.03 \quad r=0.86$$

For 138 human samples ranging in values from 0–50 µg/L, linear regression analyses of the results yielded:

$$\text{CanAg SCC} = 0.82 \times \text{Imx SCC} + 0.06 \quad r=0.98$$

## CLP (1272/2008) HAZARD CLASSIFICATION

The following warnings and precautions apply to

**SUBS** **TMB**

### Hazard pictograms:

### Signal word:

Danger

### Hazard Statement:

Repr. 1B: H360D May damage the unborn child.

### Prevention statement:

P202 Do not handle until all safety precautions have been read and understood.

### Prevention:

P280 Wear protective gloves / protective clothing / eye protection / face protection.

### Precautionary statement response:

P308+P313 IF exposed or concerned get medical advice/attention.

### Precautionary statement disposal:

P501 Dispose of contents / container to an approved hazardous / special waste disposal facility in accordance with local and national regulations.

### Restricted to professional users.

**Hazardous substances:** 2- Pyrrolidone

### Other hazards

None of the mixtures in the kit contains any substances considered to meet the criteria classifying them as PBT and/or vPvB.

## WARRANTY

The performance data presented here were obtained using the assay procedure indicated. Any change or modification of the procedure not recommended by Fujirebio Diagnostics may affect the results, in which event Fujirebio Diagnostics disclaims all warranties expressed, implied or statutory including the implied warranty of merchantability and fitness for use.

## REFERENCES

1. Suminami Y., Kishi F., Sekiguchi K., Kato H. (1991) Squamous cell carcinoma antigen is a new member of the serine protease inhibitors. *Biochem Biophys Res Commun* 181, 51-58.
2. Kato H. and Torigoe T. (1977) Radioimmunoassay for tumor antigen of human cervical squamous cell carcinoma. *Cancer* 40, 1621-1628.
3. Schneider S.S., Schick C., Fish K.E., Miller E., Pena J.C., Treter S.D., Hui S.M., Silverman G.A. (1995) A serine protease inhibitor locus at 18q21.3 contains a tandem duplication of the human squamous cell carcinoma antigen gene. *Proc Natl Acad Sci USA*, 92, 3147-3151.
4. de Bruijn H.W.A., Duk J.M., van der Zee A.G.J., Pras E., Willemse P.H.B., Boonstra H., Hollema H., Mourits M.J.E., de Vries E.G.E., Aalders J.G. (1998) The Clinical Value of Squamous cell Carcinoma Antigen in Cancer of the Uterine Cervix. *Tumor Biol* 19, 505-516.
5. Vassiliakopoulos T., Troupis T., Sotiropoulou C., Zacharatos P., Katsaounou P., Parthenis D., Noussia O., Troupis G., Papiris S., Kittas C., Roussos C., Zakynthinos S., Gorgoulis V. (2001) Diagnostic and prognostic significance of squamous cell carcinoma antigen in non-small cell lung cancer. *Lung Cancer* 32, 137-144.
6. Snyderman C.H., Dámico F., Wagner R., Eibling D.E. (1995) A Reappraisal of the Squamous Cell Carcinoma Antigen as a Tumor Marker in Head and Neck Cancer. *Arch Otolaryngol Head Neck Surg* 121, 1294-1297.
7. Duk J.M., Groenier K.H., de Bruijn H.W.A., Hollema H., ten Hoor K.A., van der Zee A.G.J., Aalders J.G. (1996) Pretreatment Serum Squamous cell Carcinoma Antigen: A Newly Identified Prognostic Factor in Early-Stage Cervical Carcinoma. *J Clin Oncol* 14, 111-118. 1997
8. Fleisher M. et al. (2002) Practice guidelines and recommendations for use of tumor markers in the clinic. *National Academy of Clinical Biochemistry* 15: p.19.
9. Yuyama N., et al., (2002) Analysis of novel disease-related genes in bronchial asthma. *Cytokine* 19(6), 287-296.
10. National Committee for Clinical Laboratory Standards, Evaluation of Precision Performance of Clinical Chemistry Devices. Approved Guideline EP5-A (1999).
11. Röijer E., Nilsson K., Oskarsson M., Dahlén U., Andersson I., Nilsson O. (2003) Development of Monoclonal Antibodies and Immunoassays against different forms of Squamous Cell Carcinoma Antigens (SCCA). *Tumor Biol* 24, p. 83.
12. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation protocol Number 7, Vol. 6, No 13, August (1986).

13. Grimsey P et al. Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. *J Pharmacokinet.* 2017; 2(4); 247-256.
14. Jenkins Colon P, Greene D.N. Biotin Interference in Clinical Immunoassays. *JALM* 2018; 2(6); 941-951.
15. Chun KY. Biotin Interference in Diagnostic Tests. *Clinical Chemistry* 2017; 63; 619-620

---

**CanAg<sup>®</sup>** is a registered trademark of Fujirebio Diagnostics AB

**Fujirebio Diagnostics AB**

**Elof Lindälvs gata 13**

**SE-414 58 Göteborg**

**Sweden**

**Phone + 46 31 85 70 30**

**Fax + 46 31 85 70 40**

**info@fdab.com**

**www.fdab.com**