

## Amplite™ Colorimetric Glucose-6-Phosphate Assay Kit

| Ordering Information               | Storage Conditions                         | Instrument Platform           |
|------------------------------------|--|-------------------------------|
| Product Number: 13805 (200 assays) | Keep in freezer<br>Avoid exposure to light | Absorbance microplate readers |

### Introduction

Glucose-6-phosphate (G6P) is a key intermediate for glucose transport into cells. G6P may also be converted to glycogen or starch for storage in the liver and muscles. G6P is utilized by glucose-6-phosphate dehydrogenase (G6PD) to generate the reducing equivalents in the form of NADPH. This is particularly important in red blood cells where a G6PDH deficiency leads to hemolytic anemia. AAT Bioquest's Amplite™ Colorimetric Glucose-6-Phosphate Assay Kit provides a simple, sensitive and rapid method for detecting G6P in biological samples such as serum, plasma, urine, as well as in cell culture samples. In the coupled enzyme assay, the G6P concentration is proportionally related to NADPH that is specifically monitored by a chromogenic NADPH sensor. The absorption signal can be read by an absorption microplate reader at the absorbance ratio of  $A_{575nm}/A_{605nm}$ . With the Amplite™ G6P Assay Kit, we were able to detect as little as 1  $\mu$ M G6P in a 100  $\mu$ L reaction volume.

### Kit Components

| Components                                | Amount                         |
|---|--------------------------------|
| Component A: Enzyme Probe                 | 2 bottles (lyophilized powder) |
| Component B: Assay Buffer                 | 1 bottle (10 mL)               |
| Component C: NADP                         | 1 vial                         |
| Component D: Glucose-6-Phosphate Standard | 1 vial (3.04 mg)               |

### Assay Protocol for One 96-Well Plate

#### Brief Summary

**Prepare G6P assay mixture (50  $\mu$ L) → Add G6P standards or test samples (50  $\mu$ L) → Incubate at room temperature for 30 minutes ~ 2 hours → Monitor absorbance ratio increase at  $A_{575nm}/A_{605nm}$**

*Note: Thaw kit components at room temperature before starting the experiment.*

#### 1. Prepare NADP stock solution (100X):

Add 100  $\mu$ L of H<sub>2</sub>O into the vial of NADP (Component C) to make 100X NADP stock solution.

#### 2. Prepare G6P stock solution:

Add 100  $\mu$ L of H<sub>2</sub>O or 1xPBS buffer into the vial of G6P Standard (Component D) to make 100 mM G6P standard solution.

*Note: The unused G6P standard stock solution should be divided into single use aliquots and stored in a freezer.*

#### 3. Prepare serial dilutions of G6P standard (0 to 100 $\mu$ M):

3.1 Add 10  $\mu$ L of G6P stock solution (from Step 2) into 990  $\mu$ L 1x PBS buffer to generate 1 mM G6P standard solution.

And then add 10  $\mu$ L of 1mM G6P stock solution into 990  $\mu$ L 1x PBS buffer to generate 100  $\mu$ M G6P standard solution.

*Note: Diluted G6P standard solution is unstable, and should be used within 4 hours.*

3.2 Take 200  $\mu$ L of 1 mM G6P standard solution to perform 1:3 serial dilutions to get 30, 10, 3, 1, 0.3, 0.1 and 0  $\mu$ M serial dilutions of G6P standard.

3.3 Add serial dilutions of G6P standard and G6P containing test samples into a white clear bottom 96-well microplate as described in Tables 1 and 2.

**Table 1:** Layout of G6P standards and test samples in a white clear bottom 96-well microplate

|       |       |      |      |      |      |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
|-------|-------|------|------|------|------|--|--|--|--|--|--|--|--|--|--|--|--|--|--|
| BL    | BL    | TS   | TS   | .... | .... |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| G6P 1 | G6P 1 | .... | .... | .... | .... |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| G6P 2 | G6P 2 |      |      |      |      |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| G6P 3 | G6P 3 |      |      |      |      |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| G6P 4 | G6P 4 |      |      |      |      |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| G6P 5 | G6P 5 |      |      |      |      |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| G6P 6 | G6P 6 |      |      |      |      |  |  |  |  |  |  |  |  |  |  |  |  |  |  |
| G6P 7 | G6P 7 |      |      |      |      |  |  |  |  |  |  |  |  |  |  |  |  |  |  |

Note: G6P=G6P Standards, BL=Blank Control, TS=Test Samples.

**Table 2** Reagent composition for each well

| G6P Standard                  | Blank Control                | Test Sample |
|-------------------------------|------------------------------|-------------|
| Serial Dilutions*: 50 $\mu$ L | Dilution Buffer : 50 $\mu$ L | 50 $\mu$ L  |

Note: Add the serially diluted G6P standards from 0.1  $\mu$ M to 100  $\mu$ M into wells from G6P1 to G6P7 in duplicate.

#### 4. Prepare G6P assay mixture:

4.1 Add 5 mL of Assay Buffer (Component B) into one bottle of Enzyme Mixture (Component A).

4.2 Add 50  $\mu$ L NADP stock solution (100X, from Step 1) into the bottle of Component A (from Step 4.1), and mix well.

Note: This G6P assay mixture is enough for one 96-well plate. It is unstable at room temperature, and should be used promptly within 2 hours and avoid exposure to light.

#### 5. Run G6P assay:

5.1 Add 50  $\mu$ L of G6P assay mixture (from Step 4.2) to each well of G6P standard, blank control, and test samples (see Step 3.3) to make the total assay volume of 100  $\mu$ L/well.

Note: For a 384-well plate, add 25  $\mu$ L of sample and 25  $\mu$ L assay mixture into each well.


5.2 Incubate the reaction at room temperature for 30 minutes to 2 hours, protected from light.

5.3 Monitor the absorbance ratio increase with an absorbance plate reader at  $A_{575\text{nm}}/A_{605\text{nm}}$ .

### Data Analysis

The absorbance in blank wells (with the dilution buffer only) is used as a control, and is subtracted from the values for those wells with the G6P reactions. A typical G6P standard curve is shown in Figure 1.

Note: The absorbance background increases with time, thus it is important to subtract the absorbance value of the blank wells for each data point.



**Figure 1.** G6P dose response was measured with Amplitude™ Colorimetric G6P Assay Kit in a 96-well white clear bottom plate using a SpectraMax Plus (Molecular Devices) microplate reader. As low as 1 $\mu$ M G6P in 100  $\mu$ L volume can be detected with 1 hour incubation.

### References

1. Corpas FJ, Barroso JB, Sandalio LM, Distefano S, Palma JM, Lupiáñez JA, Del Río LA (March 1998). "A dehydrogenase-mediated recycling system of NADPH in plant peroxisomes". *Biochem. J.* **330** (Pt 2): 777–84.
2. Gaskin RS, Estwick D, Peddi R (2001). "G6PD deficiency: its role in the high prevalence of hypertension and diabetes mellitus". *Ethnicity & disease* **11** (4): 749–54.
3. Wajcman H, Galactéros F (2004). "[Glucose 6-phosphate dehydrogenase deficiency: a protection against malaria and a risk for hemolytic accidents]". *C. R. Biol.* **327** (8): 711–20.
4. Frank JE (October 2005). "Diagnosis and management of G6PD deficiency". *Am Fam Physician* **72** (7): 1277–82.
5. Aster J, Kumar V, Robbins SL, Abbas AK, Fausto N, Cotran RS (2010). *Robbins and Cotran pathologic basis of disease*

**Warning: This kit is only sold to end users. Neither resale nor transfer to a third party is allowed without written permission from AAT Bioquest®, Inc. Chemical analysis of the kit components is strictly prohibited. Please call us at 408-733-1055 or e-mail us at [info@aatbio.com](mailto:info@aatbio.com) if you have any questions.**