

## Amplite™ Colorimetric Glucose Quantitation Kit

### Ordering Information

### Storage Conditions

### Instrument Platform

Product Number: 40004 (500 assays)

Keep at -20 °C and protect from light

Absorbance microplate readers

### Introduction

Glucose, a monosaccharide, is the most important carbohydrate in biology. It is a source of energy and metabolic intermediate for cell growth. As one of the main products of photosynthesis, glucose starts cellular respiration in both prokaryotes and eukaryotes. Glucose level is a key diagnostic parameter for many metabolic disorders, e.g., diabetes.

This Amplite™ Colorimetric Glucose Quantitation Kit provides a quick and sensitive method for the measurement of glucose. It uses glucose oxidase-based enzyme coupled reactions to detect glucose through the production of hydrogen peroxide, which is monitored by our Amplite™ Red peroxidase substrate. Amplite™ Red peroxidase substrate can be read by an absorbance microplate reader at ~570 nm. The assay is robust, and can be readily adapted for a wide variety of applications that require the measurement of glucose. The assay has very low background since it is run in the red visible range that significantly reduces the interference from biological samples. With the Amplite™ Colorimetric Glucose Quantitation Kit, we can detect as little as 3 µM D-glucose.

### Kit Key Features

<b>Sensitive:</b>	Detect as low as 3 µM D-glucose in solution.
<b>Continuous:</b>	Easily adapted to automation without a separation step.
<b>Convenient:</b>	Formulated to have minimal hands-on time. No wash is required.
<b>Non-Radioactive:</b>	No special requirements for waste treatment.

### Kit Components

Components	Amount
Component A: Amplite™ Red (light-sensitive)	1 vial
Component B: Assay Buffer	1 bottle (50 mL)
Component C: Horseradish Peroxidase (HRP)	1 vial (10 units)
Component D: Glucose Oxidase	1 vial (100 units)
Component E: DMSO	1 vial (200 µL)
Component F: Glucose	1 vial (144 mg)

### Assay Protocol for One 96-Well Plate

#### Brief Summary

**Prepare assay reaction mixture (50 µL) → Add Glucose standards or test samples (50 µL) → Incubate at 37 °C for 10-30 minutes → Monitor absorbance increase at OD of 570±5 nm**

*Note: Thaw all the kit components to room temperature before starting the experiment.*

#### 1. Prepare stock solutions:

- 1.1 **250X Amplite™ Red stock solution:** Add 100 µL of DMSO (Component E) into the vial of Amplite™ Red substrate (Component A). The stock solution should be used promptly. Any remaining solution should be aliquoted and refrozen at -20 °C.


*Note 1: Avoid repeated freeze-thaw cycles.*

*Note 2: The Amplite™ Red substrate is unstable in the presence of thiols such as dithiothreitol (DTT) and 2-mercaptoethanol. The final concentration of DTT or 2-mercaptoethanol in the reaction should be no higher than 10 µM. The Amplite™ Red substrate is also unstable at high pH (> 8.5). Therefore, the reaction should be performed at pH 7–8. The provided assay buffer (pH 7.4) is recommended.*



## Data Analysis

The absorbance reading in blank wells (with the assay buffer only) is used as a control, and is subtracted from the values for those wells with glucose standards and test samples. The standard curve of glucose is shown in Figure 1. *Note: The absorbance background increases with time, thus it is important to subtract the absorbance of the blank wells for each data point*



**Figure 1.** Glucose dose response was measured with Amplite™ Colorimetric Glucose Quantitation Kit (Cat #40004) on a 96-well clear bottom plate using a SpectraMax microplate reader (Molecular Devices) with path check on. As low as 3 µM glucose was detected with 30 minutes incubation (n=3).

## References

1. Delva P, Degan M, Trettene M, Lechi A. (2006) Insulin and glucose mediate opposite intracellular ionized magnesium variations in human lymphocytes. *J Endocrinol*, 190, 711.
2. Delva P, Degan M, Pastori C, Faccini G, Lechi A. (2002) Glucose-induced alterations of intracellular ionized magnesium in human lymphocytes. *Life Sci*, 71, 2119.
3. Wang XT, Au SW, Lam VM, Engel PC. (2002) Recombinant human glucose-6-phosphate dehydrogenase. Evidence for a rapid-equilibrium random-order mechanism. *Eur J Biochem*, 269, 3417.
4. Leira F, Louzao MC, Vieites JM, Botana LM, Vieytes MR. (2002) Fluorescent microplate cell assay to measure uptake and metabolism of glucose in normal human lung fibroblasts. *Toxicol In Vitro*, 16, 267.

**Warning: This kit is only sold to end users. Neither resale nor transfer to a third party is allowed without written permission from AAT Bioquest. Chemical analysis of the kit components is strictly prohibited. Please call us at 408-733-1055 or e-mail us at [info@aatbio.com](mailto:info@aatbio.com) if you have any questions.**