

Cell Meter™ Live Cell Caspase 8 Binding Assay Kit

Red Fluorescence

Catalog number: 20116
Unit size: 25 Tests

Component	Storage	Amount
Component A: iFluor 647-LETD-FMK	Freeze (< -15 °C), Minimize light exposure	1 vial
Component B: Washing Buffer	Freeze (< -15 °C), Minimize light exposure	1 bottle (100 mL)
Component C: 500X Nuclear Green™ DCS1	Freeze (< -15 °C), Minimize light exposure	1 vial (100 µL)
Component D: 500X Hoechst	Freeze (< -15 °C), Minimize light exposure	1 vial (100 µL)

OVERVIEW

Our Cell Meter™ live cell caspases activity assay kits are based on fluorescent FMK inhibitors of caspases. These inhibitors are cell permeable and non-cytotoxic. Once inside the cell, the caspase inhibitors bind covalently to the active caspases. This Cell Meter™ Live Cell Caspase 8 Activity Assay Kit is designed to detect cell apoptosis by measuring caspase 8 activation in live cells. It is used for the quantification of activated caspase 8 activities in apoptotic cells, or for screening caspase 8 inhibitors. iFluor 647-LETD-FMK, the red label reagent, allows for direct detection of activated caspase 8 in apoptotic cells by fluorescence microscopy, flow cytometer, or fluorescent microplate reader. The kit provides all the essential components with an optimized assay protocol.

AT A GLANCE

1. Prepare cells with test compounds at a density of 5×10^5 to 2×10^6 cells/mL
2. Add iFluor 647-LETD-FMK into cell solution at 1:150 ratio
3. Incubate at room temperature for 1 hour
4. Pellet the cells, wash and resuspend the cells with buffer or growth medium
5. Analyze the cells with flow cytometry using 660/20 nm filter (APC channel)

Important

Thaw all the components at room temperature before use.

KEY PARAMETERS

Flow cytometer

Excitation	640 nm laser
Emission	660/20 nm filter
Instrument specification(s)	APC channel

Fluorescence microscope

Excitation	Cy5 filter set
Emission	Cy5 filter set
Recommended plate	Black wall/clear bottom
Instrument specification(s)	Black wall/clear bottom

Fluorescence microplate reader

Excitation	640 nm
Emission	660 nm
Cutoff	665 nm
Recommended plate	Black wall/clear bottom
Instrument specification(s)	Bottom read mode

PREPARATION OF STOCK SOLUTIONS

Unless otherwise noted, all unused stock solutions should be divided into single-use aliquots and stored at -20 °C after preparation. Avoid repeated freeze-thaw cycles.

iFluor 647-LETD-FMK stock solution (150X)

Make 150 X iFluor 647-LETD-FMK stock solutions by adding 50 µL of DMSO to the vial of iFluor 647-LETD-FMK (Component A).

Note Store the unused iFluor 647-LETD-FMK stock solution at -20°C in single use aliquots in dark place

SAMPLE EXPERIMENTAL PROTOCOL

1. Culture cells to a density optimal for apoptosis induction according to your specific induction protocol, but not to exceed 2×10^6 cells/mL. At the same time, culture a non-induced negative control cell population at the same density as the induced population for every labeling condition. Here are a few examples for inducing apoptosis in suspension culture:
 - a. Treating Jurkat cells with 2 µg/mL camptothecin for 3 hours.
 - b. Treating Jurkat cells with 1 µM staurosporine for 3 hours.
 - c. Treating HL-60 cells with 4 µg/mL camptothecin for 4 hours.
 - d. Treating HL-60 cells with 1 µM staurosporine for 4 hours.

Note Each cell line should be evaluated on an individual basis to determine the optimal cell density for apoptosis induction.

2. Add iFluor 647-LETD-FMK stock solution (150X) into the cell solution at a 1:150 ratio, and incubate the cells in a 37 °C, 5% CO₂ incubator or room temperature for 1 hour.

Note The cells can be concentrated up to $\sim 5 \times 10^6$ cells/mL for iFluor 647-LETD-FMK labeling.

Note For adherent cells, gently lift the cells with 0.5 mM EDTA to keep the cells intact, and wash the cells once with serum-containing media prior to incubation with iFluor 647-LETD-FMK.

Note The appropriate incubation time depends on the individual cell type and cell concentration used. Optimize the incubation time for each experiment.

3. Spin down the cells at ~200 g for 5 minutes, and wash cells with 1 mL wash buffer (Component B) twice. Resuspend the cells in desired amount of washing buffer.

Note iFluor 647-LETD-FMK is fluorescent, thus it is important to wash out any unbound reagent to eliminate the background.

4. If desired, label the cells with a DNA stain such as Nuclear Green DCS1 (Ex/Em = 490/520 nm) for dead cells, or Hoechst (Ex/Em = 350/450 nm) for whole population of the cell nucleus stain.
5. Monitor the fluorescence intensity by appropriate application.
 - a. For flow cytometry, monitor the fluorescence intensity using 660/20 nm filter (APC channel). Gate on the cells of

interest, excluding debris.

- b. For fluorescence microscopy and fluorescent microplate reader. Place 100 μ L of the cell suspensions into each of wells of a 96-well black microtiter plate. Observe cells under a fluorescence microscope using Cy5 filter (FITC channel for Nuclear Green DCS1 staining, DAPI channel for Hoechst staining). Or monitor the fluorescence intensity using Ex/Em = 640/680 nm (cut off at 665 nm) bottom read mode with a fluorescent microplate reader.

Note If it is necessary to equilibrate the cell concentrations, adjust the suspension volume for the induced cells to approximate the cell density of the non-induced population. This adjustment step is optional if your cell treatment does not result in a dramatic loss in stimulated cell population numbers.

EXAMPLE DATA ANALYSIS AND FIGURES




Figure 1. Flow cytometric analysis of active caspase 8 using Cell Meter™ Live Cells Caspase 8 detection kit in Jurkat cells. The cells were treated with 1 μ M staurosporine for 5 hours (Green) while untreated cells were used as a control (Red). The staurosporine response was inhibited by Z-VAD-FMK (caspase inhibitor) shown as blue. Cells were incubated with iFluor 647-LETD-FMK for 1 hour at RT. The fluorescent intensity was measured using NovoCyte flow cytometer with 660/20 nm filter (APC channel).

DISCLAIMER

AAT Bioquest provides high-quality reagents and materials for research use only. For proper handling of potentially hazardous chemicals, please consult the Safety Data Sheet (SDS) provided for the product. Chemical analysis and/or reverse engineering of any kit or its components is strictly prohibited without written permission from AAT Bioquest. Please call 408-733-1055 or email info@aatbio.com if you have any questions.