

**CytoWatch™ Wash-free fluorescence cell imaging buffer \*10X\***

 Catalog number: 20045, 20046  
 Unit size: 10 mL, 100 mL

Component	Storage	Amount (Cat No. 20045)	Amount (Cat No. 20046)
CytoWatch™ Wash-free fluorescence cell imaging buffer *10X*	Freeze (< -15 °C), Minimize light exposure	10 mL	100 mL

**OVERVIEW**

CytoWatch™ Wash-free fluorescence cell imaging buffer \*10X\* is a ready-to-use buffer optimized for fluorescence cell imaging. In some cases, this buffer significantly enhances the imaging signal. It is used in wash steps when performing immunohistochemistry (IHC) or immuno-labeling with tissue or 3D cell culture. CytoWatch™ Wash-free fluorescence cell imaging buffer is 10X concentrated and should be diluted to 1X with PBS before use.

**SAMPLE EXPERIMENTAL PROTOCOL**
**Typical Assay Protocol (for one 96-well plate)**

1. Thaw 10X CytoWatch™ Wash-free fluorescence cell imaging buffer at room temperature before use.  
  
**Note:** It is OK to use if the buffer has precipitates.
2. Prepare a 1X CytoWatch™ Wash-free fluorescence cell imaging buffer by adding 1 mL of 10X CytoWatch™ Wash-free fluorescence cell imaging buffer to 9 mL of HHBS (1X Hank's with 20 mM Hepes buffer, pH 7.0, cat#20001) or a buffer of your choice, and mix well.  
  
**Note:** 10 mL of 1X staining buffer is enough for one plate. The buffer is stable at room temperature. It is recommended to aliquot and store any unused 10X assay buffer at < -20 °C. Protect from light. Avoid repeated freeze-thaw cycles.
3. Add the cell staining dye stock solution (generally, a concentrated DMSO solution) into 1X CytoWatch™ Wash-free fluorescence cell imaging buffer (from Step 2) to make the final well concentration 2X of the desired concentration.
4. Add the 2X Assay Solution (from Step 3) to the microplate well, making sure it's the same volume as the cell culture medium (e.g., 100 µL/well/96-well or 25 µL/well/384-well).
5. Incubate the cells in a 37 °C, 5% CO2 incubator, or as desired.

**Note:** The staining dye has the potential to disrupt the effectiveness of the 1X CytoWatch™ Wash-free fluorescence cell imaging buffer. If this occurs, it is advisable to utilize a preferred cell staining method and swap out the cell staining solution with either the cell growth medium or HHBS. Following this, add 100 µL/well/96-well (25 µL/well/ 384-well) of 1X CytoWatch™ Wash-free fluorescence cell imaging buffer into each respective well.

6. Observe the cells with a fluorescence microscope or a plate reader as required.

**DISCLAIMER**

AAT Bioquest provides high-quality reagents and materials for research use only. For proper handling of potentially hazardous chemicals, please consult the Safety Data Sheet (SDS) provided for the product. Chemical analysis and/or reverse engineering of any kit or its components is strictly prohibited without written permission from AAT Bioquest. Please call 408-733-1055 or email [info@aatbio.com](mailto:info@aatbio.com) if you have any questions.