

HiScreen MabSelect Xtra

Product Information

Cat#No# Hi-011P

Product Overview

HiScreen MabSelect Xtra is an excellent choice for method optimization and parameter screening for capture of monoclonal antibodies (MAbs).

Description

MabSelect Xtra is a recombinant protein A-based affinity medium engineered to give an exceptionally high dynamic binding capacity for monoclonal antibodies. In addition, MabSelect Xtra is optimized for Fc-fusion proteins. Higher capacity translates directly to lower cost of production.

Characteristic

High dynamic binding capacity and suitable for high expression feedstocks.
Improved process economics through reduced raw materials costs and/or reduced number of cycles.
Very low, unspecific binding due to high ligand selectivity and matrix hydrophilicity.
High capacity for many Fc-fusion proteins.

Maximum operating pressure

3 bar [0.3 MPa] (44 psi)

Average particle size

75 µm

Ligand

Recombinant protein A (E. coli).

Coupling chemistry

Epoxy

Dynamic binding capacity

~ 40 mg human IgG/ml medium

Recommended flow rate

HiScreen MabSelect Xtra

100–300 cm/h

Recommended column height

100 mm

Chemical stability

Stable in all aqueous buffers commonly used in protein A chromatography: 0.1 M 20% sodium citrate/HCl (pH 3), 6 M Gua-HCl, 6 M urea, 20% ethanol, 2% benzyl alcohol

pH working range

3–10

CIP stability

2–12

Temperature stability

2°C–40°C

Storage

2°C to 8°C in 20% ethanol.

Shipping

20% ethanol

Elution buffer

0.1 M sodium citrate, pH 3.0 to 3.6.

Cleaning-in-place

Precipitated or denatured substances: 1. Wash with 2 column volumes (CV) of 6 M guanidine hydrochloride, contact time at least 10 min. 2. Wash immediately with at least 5 CV filtered start buffer. or 1. Wash with 2 CV 50 mM NaOH of 1.0 M NaCl or 50 mM NaOH in 0.5 M Na₂SO₄, contact time ~ 10 min. 2. Wash immediately with at least 5 CV filtered start buffer.

Hydrophobically bound substances: 1. Wash with 2 CV 50 mM NaOH 1.0 M NaCl or 50 mM NaOH in 0.5 M Na₂SO₄, contact time ~ 10 min. Wash with 2 CV nonionic detergent (e.g., conc. 0.1%), contact time ~ 10

HiScreen MabSelect Xtra

min. 2. Wash immediately with at least 5 CV start buffer. or 1. Wash with 3 to 4 CV 70% ethanol or 30% isopropanol, contact time ~ 10 min. Increasing gradients may be applied to avoid air bubble formation when using high concentrations of organic solvents. 2. Wash immediately with at least 5 CV start buffer.

Sanitization

1. Wash the column with 0.1 M acetic acid in 20% ethanol. 2. Leave the column in contact with the solution for 1 hour. 3. Wash with at least 5 CV sterile start buffer at pH 7 to 8. or 1. Wash the column with 70% ethanol. 2. Leave the column in contact with the solution for 12 hours. 3. Wash with at least 5 CV sterile start buffer at pH 7 to 8.

Pack size

1 × 4.7 mL

Maximum flow velocity

300 cm/h

Dimensions

7.7 × 100 mm

Column volume

4.7 ml

Column i.d.

7.7 mm

Column hardware pressure limit

0.8 MPa (8 bar)
