

Superose 6 Increase columns

Product Information

Cat#No# Su-323P

Product Overview

Superose 6 Increase SEC columns are designed for high-resolution, small-scale preparative purification, and analysis of large globular proteins, protein complexes, and other macromolecules in the molecular weight range $M_r \sim 5K$ to $\sim 5M$.

Description

Superose 6 Increase prepacked columns are designed for rapid separation and analysis of proteins and other biomolecules by size exclusion chromatography (SEC, also known as gel filtration), in the molecular weight (M_r) range for globular proteins between 5000 and 5 000 000. This very wide fractionation range of the chromatographic resin makes it suitable for purification of protein complexes, membrane proteins, and other macromolecules. The columns are also useful as a screening tool to explore the molecular-weight distribution of unknown samples. This new generation SEC resin replaces the well-known Superose 6 columns.

Characteristic

Wide fractionation range: optimized to purify large biomolecules, including membrane proteins, protein complexes, and DNA fragments. Very high resolution: small bead size and a narrow particle size distribution provide high resolution, for high protein purity.

High flow rates: rigid beads give excellent pressure/flow properties.

Enhanced performance: improved resolution and runtime compared to predecessor.

Matrix

Composite of cross-linked agarose.

Average particle size

$\sim 8.6 \mu m$

Recommended column height

300 to 310 mm

Chemical stability

Superose 6 Increase columns

Resistant to most solutions used in liquid chromatography except hydrocarbons, aromatic solvents and chlorinated hydrocarbons.

pH working range

3 to 12

CIP stability

1 to 14

Storage

4°C to 30°C, 20% ethanol

Maximum flow velocity

1.50 mL/min, water at room temperature; 0.75 mL/min, 20% ethanol at room temperature; 0.75 mL/min water at low temperature; 0.35 mL/min 20% ethanol at low temperature.

Dimensions

10 × 300 mm

Column volume

~ 24 mL

Column i.d.

10 mm

Bed support

Filter

Material (bed support)

Polyethylene (PE) filter

Material (column tube)

Borosilicate glass

Material (column hardware)



Superose 6 Increase columns

Borosilicate glass, polyetheretherketone (PEEK), polypropylene (PP), ethylene propylene diene monomer (EPDM), and perfluoro-rubber (PFR).
