

# Flexible Plasmid Extraction

## — All-in-One Solution from Mini to Maxi

- Plasmid Mini Kit: Quickly extract 20–50 µg of high-purity plasmid DNA from 1–10 mL *E. coli* LB culture. Residual endotoxin is extremely low (<0.1 EU/µg DNA).
- Plasmid Maxi Kit: Extract 0.5–2 mg of high-copy plasmid DNA from 150–300 mL *E. coli* LB culture with yields up to 80–90%. Residual endotoxin remains ultra-low (<0.1 EU/µg DNA).



Feature	Traditional Plasmid Prep	Low-Endotoxin Plasmid Prep
Endotoxin level	Often >5 EU/µg	Ultra-low, safe for sensitive cells
Workflow	Complex, time-consuming	Streamlined, ~30 minutes
Cell impact	Can cause apoptosis, stress responses, false positives	Reliable transfection and stable cell performance
Applications	General cloning	Stem cell transfection, viral vector packaging, gene therapy

### Features

- Flexible Options** Choose Mini, Midi, or Maxi kits to match your workflow.
- Wide Sample Range** Supports small- to large-scale bacterial cultures (2–300 mL).
- High DNA Purity** Free from protein and RNA; ready for digestion, transformation, sequencing, and sensitive applications.
- High Yield** Mini kits up to 50 µg; Midi/Maxi kits up to 4 mg plasmid DNA.
- Ultra-Low Endotoxin** <0.01 EU/µg for mid-to-high copy plasmids, much lower than conventional methods.
- Easy Operation** Push-filter design simplifies handling; compatible with various rotor types.

### Performance

#### 1. Plasmid Miniprep Kit

##### (A). High Yields

Sample No.	Kit	Concentration (ng/µL, Qubit)	Concentration (ng/µL, Nanodrop)	A <sub>260/280</sub>	A <sub>260/230</sub>
1	Supplier F*	29.6	72.836	1.948	2.335
2		32.4	81.289	1.941	2.349
3	Yeasen-19001	74.8	184.149	1.938	2.335
4		78.2	171.941	1.925	2.318
5	Supplier M*	69.6	141.349	1.877	2.32
6		52.8	134.178	1.882	2.325




Figure 1. Comparison of Plasmid DNA Yields from Equal Volumes of Bacterial Culture  
Plasmid yield: Yeasen-19001 > Supplier M\* > Supplier F\*.

## Yeasen Biotechnology

### Leading Innovation in Molecular Enzymes and Reagents

[www.yeasenbio.com](http://www.yeasenbio.com)

US site: 209 Perry Pkwy, Suite 13, Gaithersburg, MD 20877, USA  
Singapore Site: 55 AYER RAJAH CRESCENT #07-05 SINGAPORE 139949

+1 (240)-472-6069 (U.S.)    info@yeasenbio.com

Y0030-V1



Free Sample

## (B). High Quality with Ultra-low Endotoxin




Figure 2. Transfection efficiency of plasmid DNA extracted using Yeasen 19021, Supplier O\*, and Supplier T\* kits.

HEK293 cells were transfected with pEGFP-C1, and bright-field and green fluorescence images were captured 72 hours post-transfection. Fluorescence intensity indicates transfection efficiency in the order: Yeasen 19021 > Supplier T\* > Supplier O\*.

## 2. Plasmid Maxiprep Kit


### (A). Higher Plasmid Yield and Purity

Kit	Concentration (ng/ $\mu$ L)	First Elution Volume (mL)	Yield 1st Elution ( $\mu$ g)	Yield 2nd Elution ( $\mu$ g)	1st Elution Recovery Rate	A <sub>260/280</sub>	A <sub>260/230</sub>
Supplier A*	853.0	1.03	879	1068	82.30%	1.81	2.29
Supplier A*	631.6	1.00	630	739	85.24%	1.81	2.26
19038ES	942.9	1.07	1005	1189	84.53%	1.84	2.41
19038ES	928.8	1.06	983	1161	84.66%	1.84	2.41

Result: The Hieff™ Pro Endo-free Plasmid Maxi Kit (19038ES) achieved 1.26x higher plasmid yield compared with Supplier A\*, while maintaining comparable purity.

### (B). Lower Endotoxin Residual

Kit	Endotoxin Residual	Kit	Endotoxin Residual
Supplier A*	~0.1 EU/ $\mu$ g	19038ES	<0.01 EU/ $\mu$ g
Supplier A*	~0.1 EU/ $\mu$ g	19038ES	<0.01 EU/ $\mu$ g
19037ES	~0.5 EU/ $\mu$ g		
19037ES	~0.5 EU/ $\mu$ g		



Result: The 19038ES kit consistently delivered plasmid DNA with endotoxin levels below 0.01 EU/ $\mu$ g, significantly lower than Supplier A\*.

## Customer Feedback

From 200 mL bacterial culture, up to 1.8 mg plasmid DNA was recovered with excellent purity and ultra-low endotoxin, fully meeting daily transfection needs.

User-friendly workflow – the push-column filter allows easy removal of bacterial debris with a single pull.

Kit	Concentration (ng/ $\mu$ L)	Yield ( $\mu$ g)	A <sub>260/280</sub>	A <sub>260/230</sub>	Endotoxin Residual
Supplier B*	943.5	1415.3	1.83	2.21	~0.5 EU/ $\mu$ g
Supplier B*	852.5	1278.8	1.84	2.23	~0.5 EU/ $\mu$ g
19038ES	1221.4	1832.0	1.82	2.33	<0.1 EU/ $\mu$ g
19038ES	1264.1	1896.2	1.82	2.30	<0.1 EU/ $\mu$ g

Result: Compared to Supplier B\*, the 19038ES kit yielded higher DNA concentrations and lower endotoxin levels (<0.1 EU/ $\mu$ g vs ~0.5 EU/ $\mu$ g), ensuring superior plasmid quality for sensitive applications such as stem cell transfection and viral packaging.

## Ordering Information

Category	Name	Cat. No.	Size
Plasmid miniprep	Hieff™ Plasmid Mini Kit	19001ES50/70	50 T/200 T
	Hieff™ Endo-free Plasmid Mini Kit	19021ES50/70	50 T/200 T
Plasmid midiprep	Hieff™ Endo-free Plasmid Mini-Midi Kit	19022ES08/50/70	8 T/50 T/200 T
	Hieff™ Endo-free Plasmid Midi Kit	19023ES02/20	2 T/20 T
Plasmid maxiprep	Hieff™ Endo-free Plasmid Maxi Kit V2	19037ES02/10	2 T/10 T
	Hieff™ Pro Endo-free Plasmid Maxi Kit	19038ES02/10	2 T/10 T

## Citations & References

[1] Zhang C, Liu Y, Zhang T, Lv C, Zang J, Zhao G. Structural comparison between the DNA-protective ability of scallop and shrimp ferritin from iron-induced oxidative damage. *Food Chem.* 2022;386:132827. doi:10.1016/j.foodchem.2022.132827 (IF:7.514)

[2] Hei H, Gao J, Dong J, et al. BK Knockout by TALEN-Mediated Gene Targeting in Osteoblasts: KCNMA1 Determines the Proliferation and Differentiation of Osteoblasts. *Mol Cells.* 2016;39(7):530-535. doi:10.14348/molcells.2016.0033 (IF:2.670)

[3] Han P, Cao P, Yue J, et al. Knockdown of hnRNPA1 Promotes NSCLC Metastasis and EMT by Regulating Alternative Splicing of LAS1L exon 9. *Front Oncol.* 2022;12:837248.

Published 2022 Jun 23. doi:10.3389/fonc.2022.837248 (IF:5.738)

[4] Zhang J, Wang W, Feng N, Jiang X, Zhu S, Chen YQ. Ndufa6 regulates adipogenic differentiation via Scd1. *Adipocyte.* 2021;10(1):646-657. doi:10.1080/21623945.2021.2007590 (IF:4.534)

[5] Zhu S, Zhang J, Wang W, Jiang X, Chen YQ. Blockage of NDUFB9-SCD1 pathway inhibits adipogenesis: Blockage of NDUFB9-SCD1 pathway inhibits adipogenesis. *J Physiol Biochem.* 2022;78(2):377-388. doi:10.1007/s13105-022-00876-7 (IF:4.158)