


# WSL-1850 *CytoWatcher II*

*Compact Live Cell Imaging System*



***CytoWatcher II*** is a digital microscope that allows you to image living cells and tissues over long periods of time without damaging them. The space-saving and power-saving design allows it to fit compactly inside a CO<sub>2</sub> incubator, minimizing the impact on the environmental temperature. It is compatible with not only bright field imaging but also fluorescence imaging, and the color imaging with excellent reproducibility is suitable for observing tissue sections and biological tissues.






# Digital microscope ideal for live cell imaging



1

## 5 million pixel color CMOS sensor that captures cell dynamics

The CMOS image sensor used in CytoWatcher II has improved performance compared to conventional sensors, and is now able to take clearer images than ever before. The sensor's light-receiving area uses a back-illuminated structure and a light-shielding structure, and the pixel size (area) has been increased by 1.6 times, improving sensitivity and saturation characteristics and suppressing noise. In addition, it uses a global shutter method (a sensor method that captures the entire subject at the same time and then outputs it), so even if you take an image of a fast-moving subject, the image will not be distorted and high-precision images that are accurate to the actual situation can be obtained. The resolution remains at 5 million pixels, but the increased pixel size makes it possible to capture a wider area of images. This makes it easier to find and focus on your subject. Operability has been improved and clearer images can now be obtained.



2

## Imaging while installed inside the CO<sub>2</sub> incubator




CytoWatcher II has a moisture-resistant structure, so you can perform time-lapse imaging for a long period time even if it is installed in a CO<sub>2</sub> incubator, which has a high humidity environment of over 95%. It uses a sensor that does not require cooling even inside the incubator, and has a power-saving design that consumes almost no power by switching to standby mode when not capturing. Only when capturing is the device activated, the light turns on, the camera turns on, and the image is taken. It does not become a source of unnecessary heat and minimizes the impact on the environmental temperature. Multiple devices can be installed within the incubator. Of course, even with the device inside, it does not interfere with normal cell culture environment.

3

## Space-saving installation area of B6 size

CytoWatcher II is a space-saving B6 size instrument with an installation area equivalent to two 96-well plates. It is short and very compact, so it does not take up much space even when installed inside a CO<sub>2</sub> incubator. Although it has a compact body, the stage area is large, and 10 cm dishes, T-75 flasks, and 96-well plates can also be used. It is easy to carry and can be moved and used in various locations. Of course, it is also suitable for use in limited spaces such as clean benches.




# WSL-1850 *CytoWatcher II* WSL-1850-B *CytoWatcher II FL*

4

## Just connect to the USB terminal of your PC




When CytoWatcher II is connected to a Windows PC using a USB 3.0 cable, Power is supplied and operation can be controlled using the dedicated operation software "ImageSaver T". Fluorescent light sources can also be used by simply connecting them to the USB terminal. After adjusting the focus and determining the observation area on the Live imaging screen, preparations are complete. All you have to do is click the start button to start imaging. It is possible to connect two CytoWatcher II to one PC and take images at the same timing under the same imaging conditions. This is a convenient function when performing experiments with multiple culture conditions or comparisons with controls.



5

## Compatible with fluorescence imaging inside living cells



CytoWatcher II FL is equipped with a blue LED light source (465 nm) for excitation, a short pass filter (480 nm) for the excitation light source, and a band pass filter (525/45 nm) for fluorescent imaging. Green fluorescence such as GFP in living cells can be taken images over time while they are still alive. Time-lapse imaging of bright field and fluorescence images of the same cell is also possible.



6

## Color imaging with a wide view and high reproducibility

CytoWatcher II has a wide view of 1.720 mm x 1.439 mm, making it suitable for observing tissues. Not only does it have a large field of view thanks to its low magnification lens, but it also has higher resolution than a stereo microscope, allowing you to capture precise images. The high-performance color CMOS sensor enables color photography with high reproducibility similar to the colors of actual tissues, such as animal tissue sections and plant tissues.



# Easy time-lapse imaging with intuitive operation



WSL-1850

## CytoWatcher II

### Imaging window

There's an objective lens under the window

### white LED light source

(for bright field photography)  
Turning on/off/light intensity controlled by PC



LED light source viewed from below

### Focus adjustment handle

Adjust the focus by manually turning the handle, allowing for precise adjustments.

### Control unit

Control unit for excitation LED

### 5 mega pixel color CMOS

4 x objective lenses (fixed)  
Digital zoom (~16x)

### Main body: Moisture-proof structure

Can be used in high humidity environments

WSL-1850-B

## CytoWatcher II FL

Model compatible with fluorescence imaging

### Filter holder

Equipped with an orange filter, which can be pulled out and used to reduce background during fluorescence photography.



### FLuorescent unit

Moisture-proof structure with built-in excitation LED light source, excitation filter, mirror, and fluorescence filter

Excitation light source: 465 nm


Excitation filter: 480 nm

Fluorescence filter: 525/45 nm



## Display images clearly with various display settings

"Color/monochrome", "grid", "saturation detection", and "inversion" displays make it easier to see captured images. From the settings screen, you can adjust "color/white balance", "gain", "sharpness", etc.



Switch between color/monochrome.

Display a grid on the image.

Displays areas where the brightness is saturated in red.

Inverts the tone.



## Easy-to-understand settings with a simple operation screen

CytoWatcher II has two imaging modes: “Live imaging” and “Interval imaging”. After setting the cells, simply open the dedicated software “ImageSaver T” and start imaging.

### Live imaging

Captured images are displayed in real time

### Setting sample

Place the cells and adjust the focus while watching the live recording

### Live imaging conditions

Lighting and the intensity adjustment  
 Brightfield: 1-63  
 Fluorescent: 1-127  
 Exposure time  
 Range: 1 to 10,000 ms

### Interval imaging condition

- ① Exposure time: 1 to 10,000 ms
- ② Shooting interval
- ③ Shooting period
- ④ Light source White LED: Bright field photography  
 Fluorescent LED: Fluorescent photography  
 (Simultaneous bright field and fluorescence imaging is also possible)


### Interval imaging

Automatically capture and save images at regular intervals

## Export data to various storage formats

The list of saved images will be displayed as thumbnails. You can adjust “contrast”, “scale bar”, etc. and save it under a different name. There are three storage formats: 8bit TIFF/JPEG/BMP.

### Image list




### Contrast




### Save as

contrast adjustment  
 Scale bar display



### Save image



contrast adjustment  
 Scale bar display

# Imaging for various purposes




## Observation of fibroblasts

The figure below shows an image of NIH3T3 cells taken with CytoWatcher II in live imaging mode. The area enclosed in a square will be enlarged and displayed using digital zoom.

### Condition

Exposure time: 10 ms  
Lighting: White LED lighting  
Container:  $\phi$ 35 mm dish




## Fluorescence imaging of GFP-expressing cells


The figure below shows images taken with CytoWatcher II FL of NIH3T3 cells 20 hours after transfection of the EGFP expression vector. The area enclosed in a square will be enlarged using digital zoom. The same cells were also imaged using interval photography 18 to 30 hours after gene transfection. If you are interested in the video, please access it using the QR code.

### Condition

Exposure time: 1 s  
Interval: 20 min  
Period: 2 day  
Lighting: Exciting LED lighting  
Container:  $\phi$ 35 mm dish



## Wound healing assay



A scratch (wound) was created on confluent NIH3T3 cells using a chip, and the healing process was taken images at 10-minute intervals using CytoWatcher II. The image on the left shows images every 2 hours. You can observe the dynamic movement of cells. Please access the video using the QR code.



### Condition

Exposure time: 10 ms  
Interval: 10 min  
Period: 1 day  
Lighting: White LED lighting  
Container:  $\phi$ 35 mm dish




## Tissue section stained image


The figure below shows color images of various tissue sections taken using CytoWatcher II in live imaging mode. Color imaging with excellent reproducibility is also suitable for observing tissue sections.

### Condition


Exposure time: 10 ms  
Lighting: White LED lighting  
Slide glass




Spleen (Masson trichrome stain)




Spleen (EvG staining)




Heart (Masson trichrome stain)



Liver (Masson trichrome stain)




Intestine (HE staining)




Basal ganglia (Masson trichrome staining)

**This is different from conventional CytoWatcher**  
**Expands the imaging area while maintaining the same resolution** 


### CytoWatcher II






### Fluorescent image



### Bright field image



### CytoWatcher



The figure shows images of GFP-expressing NIH3T3 cells taken with CytoWatcher and CytoWatcher II (left: 4x) and an enlarged image of the boxed area (center). Bright field images are magnified images of different fields of the same cell, similar to fluorescent images (right). CytoWatcher II has an expanded shooting area of 1.6 times compared to the previous CytoWatcher. The image resolution remains the same as the previous model, 5 million pixels, but the increased pixel size allows for better sensitivity.

\*Enlarged bright field image and fluorescent image are different fields of view.

# Specification

Bright field imaging  
WSL-1850 CytoWatcher II



Fluorescent/Bright field imaging  
WSL-1850-B CytoWatcher II FL



Fluorescent  
imaging unit

| Model               | WSL-1850 CytoWatcher II                                                                                                           | WSL-1850-B CytoWatcher II FL                                                 |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------|
| Camera              | 5 mega pixels / Color CMOS camera (Global shutter method)                                                                         |                                                                              |
| Resolution          | 2448 x 2048 pixels                                                                                                                |                                                                              |
| Magnification       | 4x (Digital zoom : ~ 16x)                                                                                                         |                                                                              |
| Field of view size  | 1.720 mm x 1.439 mm                                                                                                               |                                                                              |
| Focus               | Manual (coarse/fine movement)                                                                                                     |                                                                              |
| Light source        | White LED (trans-illumination)                                                                                                    | White LED (trans-illumination)<br>Side illuminating epi-blue LED (465 nm)    |
| Filter              | —                                                                                                                                 | Excitation: 480 nm shortpass filter<br>Emission: 525 / 45 nm bandpass filter |
| Moisture resistance | Available at 95%RH of humidity (Usable in CO2 incubator.)                                                                         |                                                                              |
| Software            | ImageSaverT/ImageSaverT for Windows<br>Capture mode: Live / Time-lapse<br>Save format: 8bit TIFF / BMP / JPEG, Color / Monochrome |                                                                              |
| Connection to PC    | USB 3.0 x 1                                                                                                                       | USB 3.0 Type-A x1 and USB 2.0 (or 3.0) x1                                    |
| OS                  | Windows 10 / 11 (64/32bit)                                                                                                        |                                                                              |
| Language            | Japanese/English Switchable with ImageSaverT software                                                                             |                                                                              |
| Size, Weight        | 130(W) x 180(D) x 190(H) mm, 2.5 kg                                                                                               | 130(W) x 180(D) x 190(H) mm, 2.9 kg                                          |
| Power source        | USB feeding (bus power), DC5V, 7W                                                                                                 |                                                                              |

## ATTO Corporation

Head Office:  
3-2-2 Moto-asakusa, Taito-ku, Tokyo  
111-0041, JAPAN  
URL: <https://www.attoeng.site/>

