



# Condalab

Inspired by knowledge


## **MICROBIOLOGICAL ANALYSIS IN THE FOOD INDUSTRY**

PROCEDURE IN ACCORDANCE WITH ISO STANDARDS



**Edition No. 3**

# Index




# Who are we?

European leaders  
in the manufacture  
of culture media.

Founded in 1960, we are one of the leading manufacturers of dehydrated culture media for microbiology and molecular biology in Europe and have positioned ourselves as a leading private company in the international market.

From our factory located in Madrid, Spain we export to more than 130 countries worldwide directly or through an extensive network of authorized distributors.

The key to success is our distribution channel, together with a professional team and the wide range of products we offer.



**Experts.**  
Mastery of culture media.

**Flexible and reliable.**  
Tailored to the needs of our customers.


**Innovative.**  
Inspiring the future.



## What do we do?

Develop, manufacture and distribute high-quality culture media for microbiology and molecular biology.

The experience gained in the design and manufacture of culture media has made us specialists. We develop, produce and distribute culture media of the highest quality for microbiology and molecular biology with the design of more than 700 dehydrated media. Condalab is also known for providing key ingredients such as agar, peptones and agarose, among others. Our catalog also includes media for molecular biology.



### Microbiology.

Dehydrated culture media  
Prepared culture media  
Supplements  
Microbial sensitivity tests  
Colorants  
Condagene®



### Molecular biology.

Dehydrated culture media  
Agaroses  
Dyes for molecular biology



### Bioingredients.

Agaroses  
Peptones  
Carbohydrates

# Who do we do it for?

We have an extensive network of clients as a result of active listening and the search for optimal solutions.

Condalab products target the following market niches:



## Quality control.

Food and beverage industry, breweries, pharmaceutical industry and cosmetics industry.



## Clinical analysis.

Hospitals, veterinary clinics, clinical and food control laboratories.



## Production processes.

Fermentation processes, vaccines, probiotics and culture media manufacturers.



## R+D.

Laboratories, research centers and universities.




# How do we do it?

## We opt for quality.

We continue to improve and increase our production to achieve the highest quality standards. We have ISO 9001:2015, ISO 13485:2018 and the CE mark for Invitro medical devices.

Our formulations meet the international standards of European Pharmacopoeia, FDA, APHA, USP and AOAC. We follow strict controls throughout production before, during and after each manufacturing process, to ensure quality and consistency from batch to batch.





# THE ROLE OF STANDARD METHODS

How important are  
standard methods  
in food safety?

# Foods that are not microbiologically tested can become a health risk.

MICROORGANISMS RESPONSIBLE FOR MOST CASES AND NOTIFICATIONS IN THE EU


*Campylobacter\**    *Salmonella*    *Yersinia*    STEC    *Listeria*

*\*Most commonly reported since 2005*

MAIN MICROORGANISMS ASSOCIATED WITH HIGH RATES OF DEATH OR SERIOUS DISEASES

*Salmonella*    STEC    *Listeria*





Foods that are not microbiologically tested can become a health risk because of the possible diseases they can cause. Diarrhoeal diseases, for example, are the leading cause of death in children and the second in adults, and in many cases are related to the ingestion of contaminated food.

**Microbiological criteria must therefore be established to** ensure food safety and protect consumer health. This implies, among other things, setting out the **analytical method** to use to ensure compliance with these criteria.

In addition, the use of **internationally recognized standardized methods** for microbiological analysis contributes to the generation of comparable results and increases confidence in the results obtained, which ensures the quality and safety of the food being tested.

That is why more and more protocols developed in the **ISO and CEN Committees** are indicated as **standard methods** both by accreditation bodies and to establish microbiological criteria (EU Regulations 2073/2005 and 1441/2007).

Currently, the strategy for the development and revision of **ISO methods** includes the simplification and incorporation of new technologies for optimization and better application in food control laboratories.

From Condalab we place at your disposal the whole range of culture media under **ISO formulations**, as well as the various **methods of analysis** to facilitate implementation in the various food quality control laboratories.

HOW TO READ A

# WORKFLOW



## Heading

## Type of method

# DETECTION AND COUNT OF LISTERIA MONOCYTOGENES AND LISTERIA SPP.

PART 1: DETECTION METHOD  
PROCEDURE IN ACCORDANCE WITH ISO 11290-1:2017

ISO standard, parts,  
amendments (Amd.)  
and year of publication

## Enrichment.

### (A) Primary culture

25 g/ml sample + 225 ml Listeria 1/2 Fraser Broth

30°C | 25 h ± 1 h



### (B) Second enrichment

0.1 ml of primary culture + 10 ml Listeria Fraser Broth

37°C | 24 h ± 2 h

Culture medium and  
proportions

Incubation conditions

## Presumed isolation.

0.1 ml of A and B separately in ALOA Agar

37°C | 24 - 48 h ± 2 h

**Listeria spp.**  
Blue-green colonies  
**L. monocytogenes**  
Blue-green colonies and presence of opaque halo



0.1 ml of A and B separately in Oxford Agar or Palcam Agar

35°C ± 2°C | 24 - 48 h ± 2 h

Consult TDS for the medium used for the  
identification of suspicious colonies

2 mandatory  
steps to perform

Colonies reading

## Isolation suspect colonies.

Once suspected colonies have been identified, seed  
maximum 5 colonies/plate in a non-selective agar:  
Nutritive Agar, Blood N°2 Agar or TSYE Agar

37°C | 24 - 48 h ± 2 h

## Confirmation.

### Listeria spp.

Microscopy  
**Short bacilli or cocobacilli**  
Catalase activity  
**Positive**

### Listeria monocytogenes

Microscopy  
**Short bacilli or cocobacilli**  
Hemolytic activity  
**Positive**  
Fermentation of carbohydrates  
• Rhamnose: **Positive** • Xylose: **Negative**

Tests and results for  
each

Culture medium and  
reagents

The image features a background of numerous spherical cells, likely red blood cells, viewed under a microscope. The cells are arranged in a somewhat regular pattern, with some showing a distinct biconcave disc shape. The overall color palette is a range of reds, from deep maroon to bright pink. A dark red, semi-transparent vertical bar runs down the center of the image, serving as a backdrop for the text. The text is centered within this bar and is rendered in a clean, white, sans-serif font.

**Sample  
preparation**

## Rules for consultation.

---

**ISO 6887-2:2017** Preparation of meat and meat products.

**ISO 6887-3:2017** Preparation of fish and fishery products.

**ISO 6887-4:2017** Preparation of miscellaneous products.

**ISO 6887-5:2021** Preparation of milk and milk products.

**ISO 6887-6:2013** Preparation of samples taken in the primary production stage.

## Initial dilution.

---


$\geq 10$  g/ml sample + (9\*g/ml of the sample) ml of diluent.

  $\leq 27^{\circ}\text{C}$  |  $\leq 45$  min

## Additional dilutions.

---

1 ml of initial dilution + 9 ml diluent.

  $\leq 27^{\circ}\text{C}$  | 5 s

Repeat this step with the new dilution  $10^{-2}$  the times necessary (n) to obtain the desired number of dilutions  $10^{-n}$ .

## Alternative dilutions.

---

Required in special cases, for example 1:2 or 1:5 are prepared identically respecting the ratio between the initial suspension and the diluent.

A vertical strip of a scanning electron micrograph (SEM) showing numerous rod-shaped bacteria with long, thin flagella. The image is divided into three vertical sections: a blue-tinted left section, a central red-tinted section, and a blue-tinted right section. The text 'Pathogen detection' is centered in the red section.

# Pathogen detection

# Pathogen detection

**DETECTION AND ENUMERATION OF *CAMPYLOBACTER* SPP.**

**DETECTION OF *CRONOBACTER* SPP.**

**DETECTION OF *ESCHERICHIA COLI* O157**

**DETECTION AND ENUMERATION OF *LISTERIA MONOCYTOGENES* AND *LISTERIA* SPP.**

**DETECTION, ENUMERATION AND SEROTYPING OF *SALMONELLA***

**DETERMINATION OF *VIBRIO* SPP.**

## Initial dilution.

---

See ISO 6887-1 and any international standards appropriate for the product concerned.

## Presumed isolation.

---

0.1 ml of sample or initial dilution in MYP Agar.



30°C | 18 h - 24 h\*

\*Possibility of incubating additional 24 hours.

Count for all colonies suspected of being *Bacillus cereus*.



### ***Bacillus cereus***

Pink colonies, generally surrounded by a precipitation zone.

## Isolation of suspected colonies.

---

Once suspected colonies are identified, streak at least 5 colonies in Blood N°2 Agar.



30°C | 24 h ± 2 h

## Confirmation.

---

Select well-isolated colonies from the selective medium plates and perform the following test:

Hemolytic activity: **Positive**

## Enrichment.

### A

- ↓ N° *Campylobacter* spp.
- ↓ n° flora accompanying

10 ml/g in 90 ml  
of Bolton Broth.



37°C | 4 h - 6 h

+

41,5°C | 44 ± 4 h

### B

- ↓ *Campylobacter* spp.
- ↑ n° flora accompanying

10 ml/g in 90 ml  
of Preston Broth.



41,5°C | 24 ± 2 h

### C.

- ↑ N° *Campylobacter* spp.

Streak directly into  
the CCDA Agar.

## Presumed isolation.

### CCDA Agar



#### *Campylobacter* spp.

Typical gray colonies, often with metallic shine, flat and wet, with a tendency to spread.

### Secondary\*: E.g. Preston Agar

\*Optional in cases A) and C)



41,5°C | 44 ± 2 h

## Isolation of suspected colonies.

Non-selective blood agar (eg: Columbia Agar).



41,5°C | 24 h - 48 h

NOTE: Cultivate in all cases in microaerobic atmosphere.

## Confirmation.

Morphology and mobility: **Small mobile curved bacilli**

Microaerobic growth at 25°C: **Negative**

Aerobic growth at 41,5°C: **Negative**

Oxidase: Positive

## Initial dilution.

---

Prepare the samples and dilutions, from the sample if liquid or from the initial suspension in other products, in accordance with ISO 6887.

## Presumed isolation.

---

0.1 ml in CCDA Agar.



41,5°C | 44 ± 2 h

\*Possibility of incubating for additional 24 hours



***Campylobacter* spp.**

Typical gray colonies, often with metallic shine, flat and wet, with a tendency to spread.

## Isolation of suspected colonies.

---

Once identified, select plates with < 150 suspected colonies and take 5 colonies for subculture and confirmation.

Non-selective blood agar (e.g.: Columbia Agar).



41,5°C | 24h - 48 h

NOTE: Cultivate in all cases in microaerobic atmosphere.

## Confirmation.

---

Morphology and mobility: **Small mobile curved bacilli**

Microaerobic growth at 25°C: **Negative**

Oxidase: **Positive**

## Enrichment.

See the ISO 6887 series of standards and any international standards appropriate for the product concerned.

### (A) Primary culture

10 g/ml sample +  
90 ml Buffered Peptone Water



34°C - 38°C | 18 ± 2 h



### (B) Second enrichment

0.1 ml of primary culture +  
10 ml CSB



41,5°C | 18 ± 2 h

## Presumed isolation.

0, 1 ml of **B** in CCI Agar.



41,5°C | 24 h ± 2 h



### ***Cronobacter* spp.**

Blue to blue-green colonies.

### **Others in CCI agar**

White colonies without/with green, gray or black center  
Yellow or red colonies.

## Isolation of suspected colonies.

Once suspected colonies have been identified, streak 5 colonies in a non-selective agar: TSA Agar.



34°C - 38°C | 21 h ± 3 h

## Confirmation.

Select well-isolated colonies from selective media plates to perform the following biochemical confirmation tests:

Oxidase: **Negative**

Hydrolysis of a  $\alpha$ -D PNP substrate glucopyranoside:  
**Positive**

L-lysine decarboxylase: **Negative**

L-ornithine decarboxylase: **Negative**

Fermentation of carbohydrates:

- Sucrose: **Positive**

OPTIONAL

Methyl red: **Strain dependent**

Voges-Proskauer: **Negative**

## Enrichment.

X g/ml sample + 9 g/ml mTSB + N.


 41,5°C | 18 h - 24 h


## Presumed isolation.

Concentration by immunomagnetic separation (IMS) at 6h and 12h - 18h of incubation.

NOTE: At 6h it can give a presumed positive result that can turn negative after incubating for 12h - 18h.

50 µl of magnetic particles in CT-SMAC Agar.

 37°C | 18 h - 24 h

 ***E. coli* O157:**  
Transparent colonies with a pale brown-yellowish appearance.

50 µl of magnetic particles in CondaChrome® *E. coli* O157:H7 Agar.

 37°C | 18 h - 24 h

 ***E. coli* O157:**  
Colonies with a pale pink appearance.

## Isolation of suspected colonies.

Once the suspected colonies have been identified in each isolation medium, streak at least 5 colonies in Nutritive Agar.

 37°C | 18 - 24 h

## Confirmation.

Select well-isolated colonies from selective medium plates to perform the following tests:

### BIOCHEMISTRY

Indol: **Positive**

### SEROLOGICAL

Only for positive indol colonies

### OPTIONAL

Flagellar antigens

Pathogenic characteristics

NOTE: Negative indol mutations have been found

# DETECTION AND ENUMERATION OF *LISTERIA MONOCYTOGENES* AND *LISTERIA* SPP.

PART 1: DETECTION METHOD  
PROCEDURE IN ACCORDANCE WITH ISO 11290-1:2017

## Enrichment.

### (A) Primary culture

25 g/ml sample + 225 ml Listeria 1/2 Fraser Broth.

 30°C | 25 h ± 1 h

### (B) Second enrichment


0.1 ml of primary culture + 10 ml Listeria Fraser Broth.


 37°C | 24 h ± 2 h



## Presumed isolation.


0.1 ml of A and B separately in Agar ALOA.


 37°C | 24 - 48 h ± 2 h

 **Listeria spp.**  
Blue-green colonies.  
**L. monocytogenes**  
Blue-green colonies and presence of  
opaque halo.




0.1 ml of A and B separately in Oxford Algar  
or Palcam Agar.

 35°C ± 2°C | 24 - 48 h ± 2 h

 Consult TDS for the medium used for the  
identification of suspicious colonies.

## Isolation of suspected colonies.

Once suspected colonies have been identified, streak maximum 5 colonies/plate in a non-selective agar: Nutritive Agar, Blood Agar N°2 or Agar TSYE.

 37°C | 24 - 48 h ± 2 h

## Confirmation.

### *Listeria* spp.

Microscopy  
**Short bacilli or cocobacilli**

Catalase activity  
**Positive**

### *Listeria monocytogenes*

Microscopy  
**Short bacilli or cocobacilli**

Hemolytic activity  
**Positive**

Fermentation of carbohydrates

• Rhamnose: **Positive** • Xylose: **Negative**

# DETECTION AND ENUMERATION OF *LISTERIA MONOCYTOGENES* AND *LISTERIA* SPP.

PART 2: ENUMERATION METHOD  
PROCEDURE IN ACCORDANCE WITH ISO 11290-2:2017

## Initial dilution.

See the ISO 6887 series of standards and any international standards appropriate for the product concerned.

## Presumed isolation.

0.1 ml of primary culture/dilutions in ALOA Agar.



37°C | 24 h ± 2\*

\*Possible additional incubation under the same conditions.

Count for all colonies suspected of being *L. monocytogenes* and *Listeria* spp.



### **Listeria** spp.

Blue-green colonies.

### **L. monocytogenes**

Blue-green colonies and presence of opaque halo.

## Isolation of suspected colonies.

Once suspected colonies have been identified, streak 5 colonies in a non-selective agar: Nutritive Agar, Blood N° 2 Agar or TSYE Agar.



Consult TDS for the medium used to apply the incubation conditions.

## Confirmation.

### **Listeria** spp.

Microscopy: **Short bacilli or cocobacilli**

Catalase activity: **Positive**

### OPTIONAL

Voges-Proskauer Test: **Positive**

Motility: **Positive**

### **Listeria monocytogenes**

Hemolytic activity: **Positive**

Fermentation of carbohydrates:

- Rhamnose: **Positive**
- Xylose: **Negative**

### OPTIONAL

Microscopy: **Short bacilli or cocobacilli**

Catalase activity: **Positive**

Motility: **Positive**


CAMP test: **Positive**

# DETECTION, ENUMERATION AND SEROTYPING OF SALMONELLA

PART 1: DETECTION OF *SALMONELLA* SPP.  
PROCEDURE IN ACCORDANCE WITH  
6579-1:2017 | AMD 1:2020

## Pre-enrichment.

25 g in 225 ml\* in Buffered Peptone Water.

 34°C - 38°C | 18 ± 2 h

\*For specific products consult ISO 6887.


## Selective enrichment.

0.1 ml in 10 ml of RVS Broth or MSRV Agar.

 41,5°C | 24 ± 3 h



1 ml in 10 ml of MKTTn Broth.


 34°C - 38°C | 24 ± 3 h

NOTE 1: For some products it is necessary to incubate for a further 24 ± 3 h.

NOTE 2: MSRV Agar only for mobile strains of *Salmonella* spp.: turbid gray area is observed around the inoculum.


## Presumed isolation.

### XLD Agar

 34°C - 38°C | 24 ± 3 h


*Salmonella* spp. Typical black colony with light reddish halo.  
For specific phenotypes consult TDS of the medium.

### Secondary: SS, Hektoen, Chromogenic or BGA Agar

 See TDS of the medium used for incubation conditions and identification of suspicious colonies.

## Isolation of suspected colonies.

Non-selective agar: Nutritive Agar or Sodium Chloride Enriched Nutritive Agar.

 34°C - 38°C | 24 ± 3 h

## Confirmation.

TSI Agar  
Urea Agar  
LDC medium  
Latex test

OPTIONAL  
Detection of β-galactosidase  
Indol reaction

## Initial dilution.

See ISO 6887 or ISO 8261 for the product in question.

## Enrichment.

Dilution 1:10 + Shigella Broth.



41,5°C ± 1 °C | 16 h - 20 h in anaerobiosis.

## Isolation.

Three selective agar plates are inoculated.

MacConkey Agar  
(Low selectivity)



***Shigella sonnei***

Colorless to pale pink colonies.

***Shigella spp.***

Colorless/translucent colonies  
(lactose negative).

XLD Agar  
(Intermediate selectivity)



37°C ± 1°C | 20 h - 24 h



***Shigella sonnei***

Translucent colonies with  
red/cherry center.

***Shigella spp.***

Translucent colonies with  
red/cherry center.

Hektoen Agar  
(High selectivity)



***Shigella sonnei***

Raised green wet colonies.

***Shigella spp.***

Green wet colonies.

## Isolation of suspected colonies.

Once suspected colonies have been identified, streak maximum 5 colonies/plate in a non-selective agar: Nutritive agar.



37°C ± 1°C | 20 h - 24 h

## Confirmation.

Select well-isolated colonies from the selective medium plates and perform the following biochemical and serological confirmation tests:

Triple sugar iron Agar (TSI inclined)

Semi-solid nutritive agar: **Negative mobility**

Urea Agar: **Negative**

L-lysine decarboxylation: **Negative**

Decarboxylation of L-ornithine: **Positive**

(*S. sonnei*), **Negative** (*Shigella spp.*)

Detection of Indol formation:

**Negative** (*S. sonnei*), **variable** (*Shigella spp.*)

Detection of β-galactosidase

Use of sugars

Sodium acetate (complementary):  
**no growth or very poor growth**

Antigen differentiation

Agglutination tests

# DETERMINATION OF

# VIBRIO SPP.

PART 1: DETECTION OF POTENTIALLY ENTEROPATHOGENIC *VIBRIO PARAHAEMOLYTICUS*, *VIBRIO CHOLERAЕ* AND *VIBRIO VULNIFICUS*  
PROCEDURE IN ACCORDANCE WITH ISO 21872-1:2017

## Enrichment.

Dilutions: See ISO 6887 or any other international standard specific to the product in question.

### A | Initial suspension

25 g/ml sample +  
225 ml Alkaline Peptone Water



41,5°C ± 1°C | 6 h ± 1 h  
37°C ± 1°C | 6 h ± 1 h



### B | Selective enrichment

1 ml of initial suspension +  
10 ml Alkaline Peptone Water



41,5°C ± 1°C | 6 h ± 1 h  
37°C ± 1°C | 6 h ± 1 h

NOTE: The incubation conditions shall be determined by the target species and the treatment of the product.

## Presumed isolation.

0.1 µl of **A and B** in TCBS Agar.



37°C ± 1°C | 24 h ± 3 h



***V. parahaemolyticus* and *V. vulnificus***  
greenish colonies of 2 - 3 mm  
***V. cholerae***  
yellowish colonies of 1 - 2 mm

0.1 µl **A and B** in CondaChrome® Vibrio Agar.



35°C ± 1°C | 24 h - 48 h



***V. parahaemolyticus***  
blue-green colonies  
***V. alginolyticus***  
colorless colonies  
***V. cholerae* and *V. vulnificus***  
pink colonies

## Isolation of suspected colonies.

Once suspected colonies have been identified, streak at least 5 colonies in a non-selective agar: Saline Nutritive Agar (SNA).



37°C ± 1°C | 24 h ± 3 h

## Confirmation.

Select well-isolated colonies from the selective media plates and perform the following biochemical confirmation tests:

Oxidase: **Positive**

Lysine decarboxylase in saline medium: **Positive**

Arginine dihydrolase in saline medium: **Negative**

Detection of β-galactosidase: **Species dependent**

Indol: **Positive**

Halotolerance: **Species dependent and NaCl []**

### OPTIONAL

Microscopic examination

Motility

# DETECTION OF PATHOGENIC *YERSINIA* *ENTEROCOLITICA*

PROCEDURE IN ACCORDANCE WITH ISO 10273:2017

## Enrichment.

See the ISO 6887 series of standards and any international standards appropriate for the product concerned.

### (A) Initial suspension

25 g/ml sample + 225 ml PSB Broth

 25°C | 44 h ± 4 h


### (B) Selective enrichment

10 ml of initial suspension + 90 ml ITC Broth

 25°C | 44 h ± 4 h


## Isolation.

0.1 ml of **A** directly without enrichment in CIN Agar.


 30°C | 24 h ± 2 h

### TREATMENT BY KOH<sup>a</sup>

0.5 ml separately enriched A and B + 4.5 ml of KOH solution.


 20 s ± 5 s

After treatment, **A and B<sup>a</sup>** separately in CIN Agar.

 **Pathogenic *Yersinia enterocolitica***  
Small colonies (≤1 mm) with transparent edges and dark red center.

## Isolation of suspected colonies.

Once suspected colonies have been identified, streak at least 5 colonies in a non-selective agar: Nutritive Agar, Nutritive Blood Agar or TSA Agar.

 Consult TDS for the medium used to apply the incubation conditions.

## Confirmation.

Select well-isolated colonies from selective media plates to perform the following biochemical confirmation tests:

Urea: **Positive**

Lysine decarboxylase: **Negative**

Arginine dihydrolase: **Negative**

Phenylalanine deaminase: **Negative**

Fermentation of carbohydrates:

- Sucrose: **Positive**
- Sorbitol: **Positive**

Fermentation of carbohydrates:

- Rhamnose: **Negative**
- Melibiose: **Negative**
- Citrate: **Negative**

**OPTIONAL** (BIOLOGICAL TYPING)

Aesculin


Xylose

Pirazinamidase

Tween-esterase/lipase

Trehalose

Indol

A microscopic view of various microorganisms, including rod-shaped and spherical forms, set against a dark blue background. The image is overlaid with a teal-to-purple gradient. The text "Indicator and spoilage microorganisms" is centered in white.

**Indicator and spoilage  
microorganisms**

# Indicator and spoilage microorganisms

**ENUMERATION OF MESOPHILIC LACTIC ACID BACTERIA**

**DETECTION AND ENUMERATION OF *CLOSTRIDIUM* SPP.**

31

**DETECTION AND ENUMERATION OF *ENTEROBACTERIACEAE***

**ENUMERATION OF B-GLUCURONIDASE-POSITIVE *ESCHERICHIA COLI***

## ENUMERATION OF COAGULASE-POSITIVE STAPHYLOCOCCI

39

40

41

42

## ENUMERATION OF YEASTS AND MOULDS

43

44

45

# ENUMERATION OF MESOPHILIC LACTIC ACID BACTERIA

COLONY-COUNT TECHNIQUE AT 30°C  
PROCEDURE IN ACCORDANCE WITH ISO 15214:1998

## Initial suspension.

---

See ISO 6887 series for the product concerned.

## Presumed isolation.

---

1 ml of sample or initial suspension in MRS Agar.


 30°C | 72 h ± 3 h

Plate count with 300 > 15 colonies from at least two successive dilutions.

## Confirmation.


---

Due to the possible growth in the MRS medium of species other than lactic acid ones, it may be necessary to confirm the colonies by:

- Gram staining
- Catalase test

## Initial dilution.


See the ISO 6887 series and the ISO standard specific to the product in question. Series of decimal dilutions to obtain between 10 and 150 colonies per 90 mm plate or 10 and 360 colonies per plate 140 mm plate.

 80°C | 10±1' \*Optional: Heat treatment for the selection of spores.

## Presumed isolation.

1 ml of the suspension or each dilution on an empty plate + Iron Sulfite Agar at 44°C-47°C: 90 mm plate - 12 to 15 ml or 140 mm plate - 45 to 50 ml.

Once solidified, add a second layer of Iron Sulfite Agar: 90 mm plate - 10 ml or 140 mm plate - 20 ml, and allow its solidification.

 37°C | 48 ± 2h in anaerobiosis\*



### **Sulfite-reducing bacteria**

Black-grayish to yellow-brown colonies.

## Isolation of suspected colonies.

Once suspected colonies are identified, streak 5 colonies/plate in a non-selective agar (Columbia, TSA or BHI Agar) in duplicate, to incubate under aerobiotic and anaerobiotic conditions.

 37°C | 20 ± 2

## Confirmation.

If the growth of a typical colony occurs in the incubated plate under anaerobiotic conditions and not aerobiosis, the colony belongs to the genus *Clostridia* and is counted as sulfite-reducing *Clostridium* spp.

If the growth of a typical colony occurs under both conditions, the colony does not belong to the genus *Clostridia*, therefore they cannot be counted as sulfite-reducing *Clostridium* spp.

## Initial suspension.

See the ISO 6887 series and the ISO standard specific to the product in question. Series of decimal dilutions to obtain between 10 and 150 colonies per 90 mm plate or 10 and 360 colonies per plate 140 mm plate.

 80°C | 10±1'


\*Optional: heat treatment for spore selection.

## Presumed isolation.

1 ml of the sample (liquid sample) or 10<sup>-1</sup> dilution (other products) in an empty plate + Tryptose Sulfite Cycloserine Agar at 44°C-47°C: 90 mm plate - 12 to 15 ml or 140 mm plate - 30 to 35 ml.

Once solidified, add a second Tryptose Sulfite Cycloserine Agar layer: 90 mm plate - 5 ml or 140 mm plate - 12 ml and allow to solidify.

With a sterile pipette, the process is repeated with 10<sup>-1</sup> dilutions (liquid sample) or 10<sup>-2</sup> dilutions (other products). If necessary, repeat with more dilutions, at least two consecutive decimal dilutions.

 37°C | 20 h ± 2 h in anaerobiosis



***Clostridium perfringens*:**


Black or gray to yellowish-brown colonies.

Count *C. perfringens* characteristic colonies in plates with less than 150 colonies (90 mm) or less than 365 colonies (140 mm)\*.

\*Do not leave more than 30' out of the anaerobic atmosphere before counting.

## Isolation of suspected colonies.

Once suspected colonies are identified, streak 5 colonies/plate in a non-selective agar (Columbia Agar, TSA or BHI Agar) in duplicate to incubate under anaerobic conditions.

 37°C | 20 h ± 2 h in anaerobiosis

## Confirmation.

Select well-isolated colonies to perform for the following tests:

### Acid Phosphatase Method

Spread the colonies on filter paper and add 2 to 3 drops of the acid phosphatase reagent.



***Clostridium perfringens*:**

Purplish color in 3-4'

### OPTIONAL

Differentiation of pathogenic and non-pathogenic strains for humans.

### Motility Test Method in Nitrate Motility Medium

Add colonies into Nitrate Motility Medium tubes



37°C | 20 h ± 2 h in anaerobiosis



***Clostridium perfringens*:**

Sulphite production: **positive**

Motility: **negative**

Indole production: **negative**

## Initial dilution.

---

See ISO 6887 or ISO 8261 for the product in question.

## Isolation.

---

1 ml of initial suspension in VRBL Agar medium, inoculate a plate per dilution.

 30°C or 37°C | 24 h ± 2 h

## Isolation of suspected colonies.


---

Perform the count of suspected colonies on plates with less than 150 colonies.

## Confirmation.

---

Inoculate 5 of these atypical colonies in Brilliant Green Bile Broth.

 30°C or 37°C | 24 h ± 2 h



**Coliforms:**

Pink, reddish or purple colonies, sometimes with a precipitation halo.

## Enrichment.

Prepare the samples and dilutions, depending on the type of sample, in accordance with ISO 6887, ISO 7218, ISO 8261. Dilution in Buffered Peptone Water 0.1% is recommended.

 37°C | 18 h ± 2 h

## Presumed isolation.

Inoculate a VRBG Agar plate with the enriched primary culture.

 37°C | 18 h ± 2 h




### Enterobacteria

Pink, reddish or purple colonies with and without precipitation halos.

## Isolation of suspected colonies.

Once suspected colonies have been identified, streak in a non-selective agar: Nutritive Agar.

 Consult TDS for the medium used to apply the incubation conditions.

## Confirmation.

Select well-isolated colonies from the selective media plates and perform the following biochemical confirmation tests:

Oxidase: **Negative**

Fermentation of carbohydrates

- Glucose: **Positive**

## Presumed isolation.

---

Inoculate a VRBG Agar plate with the initial suspension and/or dilutions.

 37°C | 18 h ± 2 h




### Enterobacteria

Pink, reddish or purple colonies with and without precipitation halos.

## Isolation of suspected colonies.

---

Perform the count of suspected colonies and streak 5 of these colonies/plate in a non-selective agar: Nutritive Agar.

 Consult TDS for the medium used to apply the incubation conditions.

## Confirmation.

---

Select well-isolated colonies from the selective media plates and perform the following biochemical confirmation tests:

Oxidase: **Negative**

Fermentation of carbohydrates

- Glucose: **Positive**

# DETECTION AND ENUMERATION OF ESCHERICHIA COLI

MOST PROBABLE NUMBER TECHNIQUE  
PROCEDURE IN ACCORDANCE WITH ISO 7251:2005

## Initial dilution.

---

See ISO 6887 series or ISO 8261 for the product concerned.  
For **enumeration**, use 3 tubes per dilution; in some cases up to 5 tubes are required.

## Enrichment.

---

1 ml of the initial suspension + 9 ml Lauryl Sulfate Broth or 10 ml of the initial suspension + 10 ml Lauryl Sulfate broth x2 (double concentration).



37°C | 24 h ± 2 h\*

\*Possibility of incubating for an additional 24 h if no opacity or gas production is observed.

## Isolation.

---

If opacity or gas is observed, inoculate in EC Medium.



44°C | 24 h ± 2 h

\*Possibility of incubating for an additional 24 h if no opacity or gas production is observed.

## Identification of suspected colonies

---



If visible gas is observed, inoculate in pre-heated Peptone Water at 44°C.



44°C | 48 h ± 2 h

## Confirmation.

---

Add 0.5 ml of indol reagent to the Peptone Water tubes, mix and examine after 1 min.

Presence of gas in EC Medium: **Positive**

Indol: **Positive**

# ENUMERATION OF B-GLUCURONIDASE-POSITIVE ESCHERICHIA COLI

PART 1: COLONY-COUNT TECHNIQUE  
AT 44°C USING MEMBRANES  
PROCEDURE IN ACCORDANCE WITH ISO 16649-1:2018

## Initial dilution.

---

In accordance with ISO 6887 and the standard for the specific product to be analyzed.

## Enrichment.

---

1 ml of initial suspension on MMBG\* Broth.



37°C | 4 ± 0.25 h

\*If only the initial suspension is used, prepare in duplicate.

## Presumed isolation.

---

Transfer membrane to TBX Agar.



44°C | 20 h - 24 h

## Reading of results.

---

CFU count in plates with typical colonies < 150 and < 300 (typical and non-typical).



**B-glucuronidase-positive *E. coli***

Typical colonies of blue or blue-green color.

# ENUMERATION OF B-GLUCURONIDASE-POSITIVE ESCHERICHIA COLI

PART 2: COLONY-COUNT TECHNIQUE AT 44°C  
PROCEDURE IN ACCORDANCE WITH ISO 16649-2:2001

## Initial suspension.

---

In accordance with ISO 6887 and the specific standard for the specific product to be analyzed.

## Presumed isolation.

---

1 ml of sample or initial suspension ( $10^{-1}$ ) on an empty Petri plate.

Add 15 ml of TBX Agar at a  $T^{\circ}$  of 44 - 47°C and mix carefully.



44 ± 1°C | 18 h - 24 h

Note 1: If you suspect the presence of stressed cells, pre-incubate for 4 h at 37°C, and then raise the temperature to 44°C for 18 - 24 h.

## Reading of results.

---

CFU count for plates with typical colonies < 150 and < 300 (typical and non-typical).



**B-glucuronidase-positive *E. coli***

Typical colonies of blue or blue-green color.

# ENUMERATION OF B-GLUCURONIDASE-POSITIVE ESCHERICHIA COLI

PART 3: DETECTION AND MOST  
PROBABLE NUMBER TECHNIQUE  
PROCEDURE IN ACCORDANCE WITH ISO 16649-3:2015

## Initial suspension.

---

In accordance with ISO 6887 and the specific standard for the specific product to be analyzed.

## Enrichment.

---

Inoculate\* with sample or initial suspension (10-1) in 3-5 tubes of MMBG Broth concentration 2x and 3-5 tubes of concentration 1x.



$37 \pm 1^\circ\text{C}$  |  $24 \pm 2$  h

\*Large sample volumes (10 ml) are added to equal volumes of medium 2x, while the volumes of 1 ml (or the dilution thereof) are added to 5 ml of medium 1x.

## Presumed isolation.

---

Streak from tubes with change of coloration by acid production on a TBX Agar plate.



$44 \pm 1^\circ\text{C}$  |  $22 \pm 2$  h

## Reading of results.

---

Count plates with typical colonies < 150 cfu and < 300 cfu (typical and non-typical).



**B-glucuronidase-positive *E. coli***

Typical colonies of blue or blue-green color.

# ENUMERATION OF COAGULASE-POSITIVE STAPHYLOCOCCI

PART 1: METHOD USING  
BAIRD-PARKER AGAR MEDIUM  
PROCEDURE IN ACCORDANCE WITH ISO 6888-1:2021

## Initial dilution.

See the ISO 6887 series of standards and any international standards appropriate for the product concerned.

## Presumed isolation.

0.1 ml of initial suspension in Baird Parker Agar.

See ISO 7218 for the number of plates to use in accordance with the dilutions in your analysis.



34°C - 38°C | 48 h ± 4 h



### Coagulase-positive staphylococci

Black or gray colonies surrounded by a transparent area with an opalescent ring.

After 48 h ± 4 h the colonies can lose their typical appearance, so the count may be incorrect.

## Confirmation.

To confirm the purity of the selected colony, streak the suspension in a non-selective agar: Nutritive Agar or Blood Agar.

Transfer an inoculum (suspected colony) to BHI Broth.

For the coagulase test, add 0.1 - 0.3 ml of rabbit plasma.



34°C - 38°C | 24 h ± 2 h

Coagulase test: **Positive**

# ENUMERATION OF COAGULASE-POSITIVE STAPHYLOCOCCI

PART 2: METHOD USING RABBIT PLASMA  
FIBRINOGEN AGAR MEDIUM  
PROCEDURE IN ACCORDANCE WITH ISO 6888-2:2021

## Initial dilution.

See the ISO 6887 series of standards and any international standards appropriate for the product concerned.

## Presumed isolation.

1 ml of initial suspension in 18-20 ml of RPF Medium (Baird Parker Agar + RPF supplement), freshly prepared (3 mm). Mix and let solidify, and then place in the incubator under the specified conditions.



34°C - 38°C | 24 - 48 h ± 4 h



### Coagulase-positive staphylococci

\*Black, gray or white colonies surrounded by an opaque halo of precipitation.

\*After 48 h ± 4 h the colonies may lose their typical appearance, so the count may be incorrect.

## Confirmation.

No confirmation is required, as coagulase activity is detected with the RPF Medium.

# ENUMERATION OF COAGULASE-POSITIVE STAPHYLOCOCCI

PART 3: DETECTION AND MPN  
TECHNIQUE FOR LOW NUMBERS  
PROCEDURE IN ACCORDANCE WITH ISO 6888-3:2003

## Enrichment.

See the ISO 6887 series of standards and any international standards appropriate for the product concerned.

### DETECTION

#### A | Initial suspension

1 ml initial dilution + 9 ml GC Broth (Giolitti-Cantoni)  
or 10 ml initial suspension + 10 ml GC Broth [2x]  
(double concentration).

#### B | Selective enrichment

A + 90 ml deaerated GC  
broth supplemented with  
potassium tellurite.


E.g. 5 ml/g of sample to 45 ml of GC Broth in 450 ml of supplemented GC broth.

### ENUMERATION

#### B | Selective enrichment

10 ml or 1 ml sample/dilution [2x] + 3 tubes per medium: GC Broth and  
GC Broth 2x, deaerated and supplemented with potassium tellurite.


Pour on a layer of agar or paraffin to form a seal, then allow to solidify.

 37°C | 24 h ± 2 h\*

\*Possibility of incubating a further 24 h ± 2 h if no blackening is observed.

## Presumed isolation.

With the aid of a wire loop, seed **B** in Baird Parker Agar or Medium RPF.

  $t_1$  37°C | 24 h ± 2 h  
 $t_2$  37°C | 48 h ± 2 h




Consult TDS for the medium used to identify the appearance  
of typical colonies of coagulase-positive staphylococci.

## Confirmation.

#### Baird Parker Agar

Transfer an inoculum (suspected colony)  
to BHI Broth.

Coagulase test:  
Add 0.1-0.3 ml of rabbit plasma.

 37°C | 24 h ± 2 h      Coagulase test: **Positive**

#### RPF medium

No confirmation is required, as  
coagulase activity is detected  
with the RPF Medium.

## Initial dilution.

---

Prepare the samples and dilutions, depending on the type of sample, in accordance with ISO 6887.

## Isolation.

---

1 ml of sample (liquid) / 1 ml of initial dilution ( $10^{-1}$ ) (other products) on empty Petri plate.



1 ml of  $10^{-1}$  (liquids) / 1 ml of  $10^{-2}$  (other products) on empty Petri plate.

Add 12-15 ml of PCA Agar at a T° of 44°C to 47°C and mix carefully.



30°C ± 1°C | 72h ± 3 h

Note 1: If you suspect microorganism overgrowth on the surface, once cooled, pour 4 ml more of PCA Agar over the inoculated medium as a coating layer.

## Reading of results.

---



Count for plates with number of colonies between 15 and 300.

## Initial dilution.

---

Prepare the samples and dilutions, depending on the type of sample, in accordance with ISO 6887.

## Isolation.

---

0.1 ml of sample (liquid) / 0.1 ml  
of initial suspension (other products)  
in PCA Agar.



0.1 ml of 10<sup>-1</sup> (liquid products) / 0.1 ml  
of 10<sup>-2</sup> (other products) in PCA Agar.



30°C ± 1°C | 72h ± 3 h

## Reading of results.

---



Count for plates with number of colonies between 15 and 300.

Note 1: If extended colonies are expected, the plates are examined after 24-48 h to indicate visible colonies and compare after incubation is completed.

## Initial dilution.

Prepare the samples and dilutions, depending on the type of sample, in accordance with ISO 6887, ISO 7218, ISO 8261. Dilution in Buffered Peptonate Water 0.1% is recommended.

## Isolation.

0.1 ml of sample / 0.1 of initial  
dilution in DRBC Agar.



0.1 ml of  $10^{-1}$  (liquids) / 0.1 ml of  $10^{-2}$   
(other products) in DRBC Agar.



$25^{\circ}\text{C} \pm 1^{\circ}\text{C}$  | 5-7 days

Note 1: In case of a low number of molds and yeasts, volumes of up to 0.3 ml of  $10^{-1}$  may be added.

## Reading of results.



Count for plates with number of colonies/propagules  $< 150$ .

Note 2: If necessary, count yeast colonies and molds separately.

Note 3: In the event of fast-growing molds, perform the count at 2 days and then at 5-7 days.

## Initial dilution.

Prepare the samples and dilutions, depending on the type of sample, in accordance with ISO 6887, ISO 7218, ISO 8261. Dilution in Buffered Peptonate Water 0.1% is recommended.

## Isolation.

0.1 ml of sample / 0.1 of initial  
dilution in DG18 Agar.



0.1 ml of  $10^{-1}$  (liquids) / 0.1 ml of  $10^{-2}$   
(other products) in DG18 Agar.



$25^{\circ}\text{C} \pm 1^{\circ}\text{C}$  | 5-7 days

Note 1: In case of a low number of molds and yeasts, volumes of up to 0.3 ml of  $10^{-1}$  may be added.

Note 2: If you suspect *Xeromyces bisporus*, incubate for 10 days.

## Reading of results.



Count for plates with number of colonies/propagules < 150.

Note 3: If necessary, count yeast colonies and molds separately.

Note 4: In the event of fast-growing molds, perform the count at 2 days and then at 5-7 days.

# ENUMERATION OF PRESUMPTIVE *PSEUDOMONAS SPP.*

METHOD FOR THE ENUMERATION  
PRESENT IN MEAT AND MEAT PRODUCTS  
PROCEDURE IN ACCORDANCE WITH ISO 13720:2010

## Initial dilution.


---

See the ISO 6887 series of standards (ISO 6887-1 and ISO 6887-2) and any international standards appropriate for the product in question.

## Presumed isolation.

---

0.1 ml of the initial suspension in CFC Agar.

 25°C ± 1°C | 44 h ± 4 h

 *Pseudomonas spp.*: Growth in CFC Agar.

Count colonies and retain plates with < 150 colonies.

## Confirmation.

---

Select 5 colonies randomly from each of the retained plates:

\*Oxidase: **Positive**

\*Use a plastic or platinum-iridium alloy loop.

The image shows a cover page for a 'Product list'. The background is a solid orange color. Overlaid on this are several thin, light blue lines that form overlapping circles and arcs, creating a geometric pattern. The text 'Product list' is centered in the middle of the page in a white, bold, sans-serif font.

# Product list

CAT	PRODUCT	FORMAT
1402	Buffered Peptone Water	500 g
6702	Buffered Peptone Water	2x5L
6705	Buffered Peptone Water	3x3L
6707	Buffered Peptone Water	5x2L
5020	Buffered Peptone Water	10 x 225 ml
4250	Buffered Peptone Water	20 tubes
1405	Saline Peptone Water	500 g
6710	Saline Peptone Water (Maximum Recovery Diluent - MRD)	3x3L
6711	Saline Peptone Water (Maximum Recovery Diluent - MRD)	2x5L
4101	Ringer Solution 1/4	20 tubes
1406	Buffered Saline Peptone Water	500 g
1405	Saline Peptone Water	500 g
4044	Saline Peptone Water	20 tubes
5182	Saline Peptone Water	10 x 90 ml
5153	Buffered Peptone Water	10 x 100 ml
5171	Buffered Peptone Water	10 x 90 ml

CAT	PRODUCT	FORMAT
1343	Selective Agar Base for Bacillus Cereus (MYP)	500 g
6021	Supplement for Bacillus Cereus	10 vials
5152	Egg yolk emulsion	100 ml
945	Selective Agar for Bacillus Cereus (MYP)	20 plates
1328	Blood Agar No. 2 Base	500 g
912	Blood Agar No. 2	20 plates

CAT	PRODUCT	FORMAT
1441	Bolton Broth ase for selective enrichment	500 g
4658	Bolton broth base	10 x 225 ml
4659	Bolton broth base	10 x 250 ml
4675	Bolton broth base	10 x 90 ml
6070	Bolton Selective Supplement	10 vials
2166	Preston Broth base for Campylobacter	500 g
6081	Preston Campylobacter supplement	10 vials
1129	Blood-free Campylobacter Agar Base (CCDA)	500 g
6053	CCDA Supplement (Blood-free Campylobacter)	10 vials
878	Campylobacter CCDA Agar	20 plates
1104	Columbia Agar Base	500 g
931	Columbia Agar + 5% Lamb Blood	20 plates
1131	Campylobacter Agar Base (Preston)	500 g

CAT	PRODUCT	FORMAT
1402	Buffered Peptone Water	500 g
6702	Buffered Peptone Water	2x5L
6705	Buffered Peptone Water	3x3L
6707	Buffered Peptone Water	5x2L
5020	Buffered Peptone Water	10 x 225 ml
4250	Buffered Peptone Water	20 tubes
5153	Buffered Peptone Water	10 x 100 ml
5171	Buffered Peptone Water	10 x 90 ml
2143	Cronobacter Selective Broth (CSB)	500 g
1446	Chromogenic agar for Cronobacter isolation (CCI)	500 g
866	Chromogenic Agar for Cronobacter CCI	20 plates
1068	Soy and Trypticasein Agar (TSA)	500 g
904	Soy and Trypticasein Agar (TSA)	20 plates
4003	Soy and Trypticasein Agar (TSA)	20 tubes
5000	Soy and Trypticasein Agar (TSA)	10 x 100 ml
5157	Soy and Trypticasein Agar (TSA)	10 x 200 ml
2149	Means for decarboxylation of L-Ornithine	500 g
1208	lysine Decarboxylase Broth	500 g

CAT	PRODUCT	FORMAT
1292	Trypticasein Soy Broth Modified with Novobiocin (mTSB)	500 g
1099	MacConkey Agar with Sorbitol (CT-SMAC)	500 g
6064	Cefixime Tellurite (CT) Supplement	10 vials
1588	Chromogenic Agar Base E. coli O157:H7	500 g
6064	Cefixime Tellurite (CT) Supplement	10 vials
1060	Nutritive Agar	500 g
902	Nutritive Agar	20 plates
1237	Tryptophane Culture Broth	500 g
4027	Tryptophane Broth	20 tubes
5205	Kovac reagent	100 ml

CAT	PRODUCT	FORMAT
1182	Listeria Fraser Broth Base	500 g
1183	1/2 Fraser Broth Base	500 g
4053	Listeria 1/2 Fraser Broth Enrichment	20 tubes
5022	Listeria 1/2 Fraser Broth Enrichment	10 x 225 ml
6703	1/2 Fraser Broth	2:x5L
6706	1/2 Fraser Broth	3 x x3L
6708	1/2 Fraser Broth	5x2L
6001	Selective supplement for Listeria Fraser	10 vials
6002	Selective supplement for Listeria 1/2 Fraser	10 vials

CAT	PRODUCT	FORMAT
1345	Chromogenic Agar Base for Listeria in accordance with Ottaviani And Agosti (ALOA)	500 g
1141	Listeria Palcam Agar Base	500 g
1133	Listeria Oxford Agar Base	500 g
884	Chromogenic agar for Listeria in Accordance with to Ottaviani and Agosti (ALOA)	120 plates
939	Chromogenic agar for Listeria in Accordance with to Ottaviani and Agosti (ALOA)	20 plates
955	Listeria Palcam Agar	20 plates
4024	Listeria Fraser Enrichment Broth	20 tubes
4053	Listeria 1/2 Fraser Broth Enrichment	20 tubes
5022	Listeria 1/2 Fraser Broth Enrichment	10 x 225 ml
5023	Listeria Fraser Broth Enrichment	10 x 225 ml
6003	Listeria Oxford Selective Supplement	10 vials
6004	Listeria Palcam Selective Supplement	10 vials
6031	Listeria lipasa C. Supplement	10 vials
6040	Listeria Selective Chromogenic Supplement	10 vials
1060	Nutritive Agar	500 g
902	Nutritive Agar	20 plates
1328	Blood No. 2 Agar Base	500 g
912	Blood Agar No. 2	20 plates
1398	TSYEA Agar(Tryptone Soy Yeast Extract Agar)	500 g
1342	Base of broth for the use of carbohydrates	500 g
6001	Selective supplement for Listeria Fraser	10 vials
6002	Selective supplement for Listeria 1/2 Fraser	10 vials
1345	Chromogenic Agar base for Listeria de in accordance with Ottaviani And Agosti (ALOA)	500 g
884	Chromogenic agar for Listeria	120 plates
939	Listeria chromogenic Agar	20 plates
6031	Listeria lipasa C. Supplement	10 vials
6040	Listeria Selective Chromogenic Supplement	10 vials
1060	Nutritive Agar	500 g
902	Nutritive Agar	20 plates
1328	Blood Agar No. 2 Base	500 g

CAT	PRODUCT	FORMAT
1376	Modified semi-solid medium Rappaport Vassiliadis (MSRV)	500 g
1174	Rappaport Soy Broth (Vassiliadis)	500 g
4245	Rappaport Soy Broth (Vassiliadis)	20 tubes
1173	Muller-Kauffmann Broth Base with Bright Green And Novobiacin (MKTTN)	500 g
4667	Muller Kauffmann Tetrionate-Novobiocin (MKTTN)	10 x 100 ml
4021	Muller Kauffman Tetrionate Broth Base	20 tubes
1274	XLD Agar (Xylose Lysine Seoxycholate Agar)	500 g
847	XLD Agar (Xylose Lysine Seoxycholate Agar)	20 plates
1172	Triple Sugar and Iron Agar(TSI)	500 g
4204	Triple Sugar and Iron Agar(TSI)	20 tubes
2180	Urea Agar Base (Christensen)	500 g
5100	Solution Urea 40%	100 ml
1110	Urea Agar Base (Christensen)	500 g

CAT	PRODUCT	FORMAT
1176	Decarboxylation of Lysine Medium	500 g
1208	Lysine Decarboxylase Broth	500 g
4229	Lysine Decarboxylase Broth	20 tubes
1078	Brilliant Green Agar (BGA)	500 g
915	Brilliant Green Agar (BGA)	20 plates
1064	Salmonella Shigella Agar (SS Agar)	500 g
909	Salmonella Shigella Agar (SS Agar)	20 plates
1122	Salmonella Chromogenic Agar	500 g
1030	Hektoen Enteric Agar	500 g
918	Hektoen Enteric Agar	20 plates
6141	Salmonella Latex Test	50 test

CAT	PRODUCT	FORMAT
2078	Shigella Broth	500 g
1052	MacConkey Agar	500 g
900	MacConkey Agar	20 plates
1274	XLD Agar (Xylose Lysine Seoxycholate Agar)	500 g
847	XLD Agar (Xylose Lysine Seoxycholate Agar)	20 plates
1030	Enteric Hektoen Agar	500 g
918	Enteric Hektoen Agar	20 plates
1172	Triple Sugar and Iron Agar (TSI)	500 g
4204	Triple Sugar and Iron Agar (TSI)	20 tubes
2046	Semi-solid nutritive Agar	500 g
1176	Decarboxylation of Lysine Medium	500 g
4229	Lysine Decarboxylase Broth	20 tubes
2180	Urea Agar Base (Christensen)	500 g
1237	Tryptophane Culture Broth	500 g
1192	Differential Acetate Agar	500 g

CAT	PRODUCT	FORMAT
2155	Saline Alkaline Peptone Water	500 g
4018	Alkaline Peptone Water	20 tubes
4685	Alkaline Peptone Water	10 x 225 ml
1074	TCBS Agar	500 g
957	TCBS AgarD	20 plates
1237	Tryptophane Culture Broth	500 g
4027	Tryptophane Broth	20 tubes
5205	Kovac reagent	100 ml

CAT	PRODUCT	FORMAT
1298	Peptone Sorbitol and Bile Salts Broth (PSB)	500 g
4683	Peptona Bile and Sorbitol Broth (PBS)	10 x 100 ml
1361	Irgasan Ticarcillin and Potassium Chlorate (TTI) Broth	500 g
6051	ITC supplement	10 vials
1126	Selective Yersinia Agar Base (CIN)	500 g
6033	Yersinia Selective Supplement (CIN)	10 vials
933	Yersinia CIN Agar	20 plates
1060	Nutritive Agar	500 g
902	Nutritive Agar	20 plates
1328	Blood Agar No. 2 Base	500 g
912	Blood agar No. 2	20 plates
1068	Soy and Trypticasein Agar (TSA)	500 g
904	Soy and Trypticasein Agar (TSA)	20 plates
4003	Soy and Trypticasein Agar (TSA)	20 tubes
5000	Soy and Trypticasein Agar (TSA)	10 x 100 ml
5157	Soy and Trypticasein Agar (TSA)	10 x 200 ml
1208	Lysine Decarboxylase Broth	500 g
1176	Decarboxylation of Lysine Medium	500 g
1040	Phenylalanine Agar	500 g
1014	Simmons Citrate Agar	500 g
4001	Simmons Citrate Agar	20 tubes
1031	Bile Esculina Agar	500 g
1227	Urea Indol Broth	500 g
4212	Urea Indol Broth	20 tubes
1364	Kligler Iron Agar	500 g

CAT	PRODUCT	FORMAT
4684	MRS low pH Agar	10 x 200 ml
1433	MRS Agar Low pH	500 g

CAT	PRODUCT	FORMAT
1029	TSC Agar Base (Tryptose Sulfite Cycloserine)	500 g
4660	TSC Agar (Tryptose Sulfite Cycloserine)	10 x 100 ml
4661	TSC Agar (Tryptose Sulfite Cycloserine)	10 x 200 ml
6020	Clostridium Perfringens Supplement (TSC)	10 vials
4728	TSC Agar (Tryptose Sulfite Cycloserine)	30 plates
1565	Nitrate Motility Medium Base	500 g
1559	Iron Sulfite Agar	500 g

CAT	PRODUCT	FORMAT
1093	Violet Red Bile Lactose Agar (VRBL)	500 g
910	Violet Red Bile Lactose Agar (VRBL)	20 plates
4532	Violet Red Bile Lactose Agar (VRBL)	30 plates
4682	Violet Red Bile Lactose Agar (VRBL)	10 x 200 ml
5161	Violet Red Bile Lactose Agar (VRBL)	10 x 100 ml
1228	Brilliant Green Bile Broth 2%	500 g
4213	Brilliant Green Bile Broth 2%	20 tubes

CAT	PRODUCT	FORMAT
1092	Violet Red Bile Glucose Agar (VRBG)	500 g
911	Violet Red Bile Glucose Agar (VRBG)	20 plates
4524	Violet Red Bile Glucose Agar (VRBG)	30 plates
4670	Violet Red Bile Glucose Agar (VRBG)	10 x 200 ml
5158	Violet Red Bile Glucose Agar (VRBG)	10 x 100 ml
1355	Nutritive Agar enriched with Sodium Chloride	500 g
2150	Glucose OF Medium	500 g
4026	Glucose OF Medium	20 tubes

CAT	PRODUCT	FORMAT
1310	Lauryl Sulfate Broth (Lauryl Tryptose Broth – LTB)	500 g
4039	Lauryl Sulfate Broth (Lauryl Tryptose Broth – LTB)	20 tubes
4040	Lauryl Sulfate Broth (Lauryl Tryptose Broth – LTB) C+	20 tubes
1522	EC Medium	500 g
4265	EC Medium	20 tubes
1403	Peptone Water (Tryptone Water)	500 g
4201	Peptone Water (Tryptone Water)	20 tubes
5160	Peptone Water (Tryptone Water)	10 x 100 ml

CAT	PRODUCT	FORMAT
1151	TBX Chromogenic Agar (Tryptone Bile X-Glucuronide)	500 g
982	TBX Chromogenic Agar (Tryptone Bile X-Glucuronide)	20 plates
982	TBX Chromogenic Agar (Tryptone Bile X-Glucuronide)	20 plates
982	TBX Chromogenic Agar (Tryptone Bile X-Glucuronide)	20 plates
1365	Minerals Modified Glutamate Broth (MMBG)	500 g

CAT	PRODUCT	FORMAT
1100	Baird Parker Agar Base	500 g
5129	Egg yolk Emulsion + Tellurite	100 ml
1331	Brain Heart Infusion Broth	500 g
4202	Brain Heart Infusion Broth	20 tubes
1060	Nutritive Agar	500 g
902	Nutritive Agar	20 plates
1328	Blood Agar No. 2 Base	500 g
912	Blood agar No. 2	20 plates
1100	Baird Parker Agar Base	500 g
6024	RPF supplement	10 vials
984	Baird-Parker Agar RPF	20 plates
4669	Baird Parker RPF Agar Base	10 x 90 ml
1287	Giolitti-Cantoni Broth	500 g
4228	Giolitti-Cantoni Broth	20 tubes
4261	Giolitti-Cantoni Broth (double concentration)	20 tubes
1100	Baird Parker Agar Base	500 g
5129	Egg yolk Emulsion + Tellurite	100 ml
6024	RPF supplement	10 vials
1331	Brain Heart Infusion Broth	500 g
4202	Brain Heart Infusion Broth	20 tubes

CAT	PRODUCT	FORMAT
1056	Standard Methods Agar (PCA)	500 g
903	Standard Methods Agar (PCA)	20 plates
4706	Standard Methods Agar (PCA)	30 plates
5112	Standard Methods Agar (PCA)	10 x 200 ml
5115	Standard Methods Agar (PCA)	10 x 100 ml
4105	Plate Count Agar (PCA)	20 tubes
4681	Standard Methods Agar (PCA) with Milk	10 x 200 ml
5172	Standard Methods Agar (PCA) with Milk	10 x 100 ml

CAT	PRODUCT	FORMAT
1160	Rose Bengal Agar + Chloramphenicol + Dicloran (DRBC Agar)	500 g
833	Rose Bengal Agar + Chloramphenicol + Dicloran (DRBC Agar)	20 plates
1161	Dichloran Glycerol Agar (DG 18)	500 g

CAT	PRODUCT	FORMAT
1356	Pseudomonas Agar Base CFC	500 g
6036	CFC supplement	10 vials
829	Pseudomonas Agar CFC	20 plates

The background consists of a light green field with diagonal stripes. A darker green vertical band runs through the center, containing the text. The text is white and centered within this band.

**Annexes:  
Rapid alternative methods**

## Perform your analyses faster and easier using Condagene®

Until now, traditional culture methods have been our best ally when considering reliability in our microbiological analyses. But thanks to Condagene® you can perform your analyses **more quickly and easily**, providing an ideal complement to your habitual techniques.



Time




Specificity



Sensitivity

### Salmonella

#### ISO 6579-1 method vs Condagene®



## Why choose PCR for foodborne pathogen detection?



Complex and multi-step methods



Powerful screening tool



High number of samples to be analyzed

### Pathogenic bacteria

*Salmonella* spp.; *L. monocytogenes*; *Campylobacter*; *Cronobacter*; STEC; *Vibrio* spp.

### Allergens

Almond ; Cashew; Celery; Hazelnut; Peanut; Lupin; Mustard; Pecan; Sesame; Soy bean, Walnut.

### Learn more about CondaChrome®

CondaChrome® media contain in their composition a colorless chromogenic substrate, which thanks to the specific enzymatic activities of each microorganism, degrades and results in the release of an element known as a chromophore, which gives the colony an intense, specific color that allows **identification** of the bacteria at **first sight**.



Quick results



Time/space saving




Easy interpretation



Minimum investment

### CondaChrome® Agar Salmonella




- For isolation of *Salmonella* spp.

Reading of colonies

- Magenta: *S. interiditis*, *S. typhimurium*, *S. typhi*, *S. dysenteriae*
- Colorless: *Proteus vulgaris*
- Blue green: *Escherichia coli*

### CondaChrome® Agar E. coli O157:H7




- For detection of *E. coli* O157:H7

Reading of colonies

- Pale Rose: *Escherichia coli* O157:H7
- Total inhibition: *Enterobacter*, *Salmonella*, *Escherichia coli*, *Enterococcus* and *Staphylococcus*

### CondaChrome® Agar Vibrio



- For isolation and detection of *Vibrio* spp.

Reading of colonies

- Rose: *V. cholerae*, *V. vulnificus*
- Blue green: *V. parahaemolyticus*
- Colorless: *V. alginolyticus*

## Extraction kits

CAT	PRODUCT	FORMAT
6500	Extraction Condagene® Complex	100 rxn
6504	Condagene® Extraction Rapid	100 rxn
6505	Condagene® Extraction Rapid	250 rxn

CAT	PRODUCT	FORMAT
6506	Condagene® Extraction Column	50 rxn
6507	Condagene® Extraction Column	250 rxn

## Pathogens

CAT	PRODUCT	FORMAT
6518	Condagene® Salmonella spp.	100 rxn
6519	Condagene® Listeria monocytogenes	100 rxn
6520	Condagene® Campylobacter jejuni	100 rxn
6521	Condagene® Cronobacter spp.	100 rxn
6522	Condagene® Vibrio spp.	100 rxn
6523	Condagene® STEC Screening	50+50+50 rxn
6524	Condagene® STEC ID	50x5 rxns
6541	Condagene® E. coli O157:H7 / O157	100 rxn

## Allergens

CAT	PRODUCT	FORMAT
6530	Condagene® Allergens Almond	100 rxn
6531	Condagene® Allergens Cashew	100 rxn
6532	Condagene® Allergens Celery	100 rxn
6533	Condagene® Allergens Hazelnut	100 rxn
6534	Condagene® Allergens Lupin	100 rxn
6535	Condagene® Allergens Mustard	100 rxn
6536	Condagene® Allergens Peanut	100 rxn
6537	Condagene® Allergens Pecan (Walnut)	100 rxn
6538	Condagene® Allergens Sesame	100 rxn
6539	Condagene® Allergens Soy bean	100 rxn
6540	Condagene® Allergens Walnut	100 rxn

## Equipment

CAT	PRODUCT
CDL96	CDL-96 RealTime System

CAT	PRODUCT
6503	Condagene® B-Light

## CondaChrome®

CAT	PRODUCT	FORMAT
1122	CondaChrome® Salmonella Agar	500 g
1180	CondaChrome® EC with MUG Fluorogenic Agar	500 g
1285	CondaChrome® EC with MUG Fluorogenic Broth	500 g
1340	CondaChrome® E. coli-Coliforms Chromogenic Medium	500 g
1465	CondaChrome® Lauryl Sulfate Chromogenic Broth	500 g
1585	CondaChrome® Standard Method Chromogenic Agar (PCA)	500 g
1588	CondaChrome® E. coli O157:H7 Chromogenic Agar Base	500 g
2054	CondaChrome® Vibrio Agar	500 g
2076	CondaChrome® Staphylococcus Agar	500 g

# ADDITIONAL RESOURCES



## Product list

Conda<sup>ene</sup>  
*Detection by qPCR*

Conda<sup>Chrome</sup>  
Chromogenic culture med

**Analysis in the  
dairy industry**

**Discover  
our webinars**

**Subscribe to  
our newsletter**

**Connect  
with us**





# Condalab

MICROBIOLOGICAL ANALYSIS  
IN THE FOOD INDUSTRY

[www.condalab.com](http://www.condalab.com)