

## iFluor™ 350 Anti-human CD138 Antibody \*MI15\*

Catalog number: 11380010, 11380011 Unit size: 100 tests, 500 tests

**Product Details** 

Storage Conditions 2-8°C with minimized light exposure. Do not freeze.

Expiration Date 12 months upon receiving

Concentration 0.1 mg/mL

Formulation Phosphate-buffered saline (PBS, pH 7.2), 0.09% sodium azide, 0.2% (w/v) BSA

**Antibody Properties** 

Species Reactivity Human

Class Primary

Clonality Monoclonal

Host Mouse

Isotype Mouse IgG1

Immunogen CD138 (Syndecan-1)

Clone MI15

Conjugate iFluor™ 350

**Biological Properties** 

Appearance Off-white liquid

Preparation Antibody purified by affinity chromatography and then conjugated with iFluor™ 350 under

optimal conditions

Application Flow Cytometry (FACS), Fluorescence Imaging

**Spectral Properties** 

Conjugate iFluor™ 350

Excitation Wavelength 345 nm

Emission Wavelength 450 nm

**Applications** 

The MI15 monoclonal antibody reacts with human CD138, a 100 - 200 kD transmembrane protein typically expressed on the surface of b cells, plasma cells, epithelial cells and endothelial cells. CD138 is involved with important cellular pathways, in particular, the cytokine-mediated signaling pathway and canonical Wnt signaling pathway. Also, in certain organisms, it enhances extracellular exosome assembly and is a

promoter of exosomal secretion. From a research standpoint, it is of biological interest due to its association with critical macromolecules/ligands such as Fibronectin. CD138 is a fairly uncommon antibody target, with a little more than 7000 publications in the last decade. Even still, CD138 is frequently used in flow cytometry applications as a phenotypic marker for differentiation of cell types, specifically in the study of cell motility/cytoskeleton/structure, cell biology and synaptic biology. This antibody was purified through affinity chromatography and conjugated to iFluor™ 350 (ex/em = 345/450 nm). It is compatible with the 355 nm laser and 450/50 nm bandpass filter (for example, as in the BD Special Order LSRFortessa™ Cell Analyzer).