

## Product Information

### Single Test Kit

**CAT#:** EKIT-009

Single Test Kit is designed to be a qualitative In Vitro end-point endotoxin test for human and animal parenteral drugs, biological products, and medical devices. Limulus Amebocyte Lysate (LAL) as supplied is to be reconstituted with sample or control directly. There is no need to reconstitute the system with LAL Reagent Water before use. After incubation for 1 hour, and in the presence of endotoxin, gelation occurs; in the absence of endotoxin, gelation does not occur. Our supply a series of Single Test Kits with different sensitivity (0.015 EU/ml, 0.03 EU/ml, 0.06 EU/ml, 0.125 EU/ml and 0.25 EU/ml). When using the Single Test Kit, you just need to add 200µl of sample or control to LAL, then, wait for the required incubation time to allow for gelation to occur.

CAT#	PRODUCT	Package	Storage
EKIT-009	Single Test Kit	40 assay	
EKIT-009	Single Test Kit	20 assay	

### Product Information

Product Sensitivity	0.03 EU/ml-0.25 EU/ml
Highlights	High reproducibility Competitive price Broad detection range Easy to use
Usage	For Research Use Only! Not For Use in Humans.