Catalog # SPN-Y169



Source

Anti-SARS-CoV-2 Spike Antibody, Mouse IgG1 (2C9C9) (XBB.1.5/Omicron Specific) is isolated from a Spike protein infected Mouse and is recombinantly produced from Hybridoma

Clone
2C9C9
Species
Mouse
Isotype
Mouse IgG1 Mouse Kappa
Conjugate
Unconjugated
Antibody Type
Hybridoma Monoclonal
Reactivity
Virus
Immunogen

Recombinant SARS-CoV-2 Spike Trimer Protein (XBB.1.5/Omicron) erived from human 293 cells.

Specificity

This product is a specific antibody specifically reacts with SARS-CoV-2 Spike Trimer Protein, His Tag (XBB/Omicron) (MALS verified) (Cat. No. SPN-C5248), SARS-CoV-2 Spike Trimer Protein, His Tag (XBB.1/Omicron) (MALS verified) (Cat. No. SPN-C522t) and SARS-CoV-2 Spike Trimer Protein, His Tag (XBB.1.5/Omicron) (MALS verified) (Cat. No. SPN-C524i). No cross-reactivity is detected with Protein of WT (Cat. No. SPN-C52H9), Alpha (Cat. No. SPN-C52H6), Beta, Gamma (Cat. No. SPN-C52Hg), Delta (Cat. No. SPN-C52H6), Beta, Gamma (Cat. No. SPN-C52Hg), Delta (Cat. No. SPN-C5223), BA.3/Omicron (Cat. No. SPN-C522b), BA.2/Omicron (Cat. No. SPN-C5229), BA.5/Omicron (Cat. No. SPN-C522e), BA.2.12.1/Omicron (Cat. No. SPN-C522d), BQ.1.1/Omicron (Cat. No. SPN-C522s), BA.2.75/Omicron (Cat. No. SPN-C522f), BA.4.6/Omicron (Cat. No. SPN-C522m), BF.7/Omicron (Cat. No. SPN-C522q).

Purity

>95% as determined by SDS-PAGE.

Purification

Protein A purified/ Protein G purified

Formulation

Lyophilized from 0.22 μm filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage

For long term storage, the product should be stored at lyophilized state at -20° C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

Application

Application

Recommended Usage



>>> www.acrobiosystems.com

7/24/2024

Anti-SARS-CoV-2 Spike Antibody, Mouse IgG1 (2C9C9) (XBB.1.5/Omicron Specific)



Catalog # SPN-Y169

ApplicationRecommended UsageELISA0.8-200 ng/mL

SDS-PAGE



Anti-SARS-CoV-2 Spike Antibody, Mouse IgG1 (2C9C9) (XBB.1.5/Omicron Specific) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 95% (With <u>Star Ribbon Pre-stained Protein Marker</u>).

Bioactivity-ELISA



Immobilized SARS-CoV-2 Spike Trimer Protein, His Tag (XBB/Omicron) (MALS verified) (Cat. No. SPN-C5248), SARS-CoV-2 Spike Trimer Protein, His Tag (XBB.1/Omicron) (MALS verified) (Cat. No. SPN-C522t) and SARS-CoV-2 Spike Trimer Protein, His Tag (XBB.1.5/Omicron) (MALS verified) (Cat. No. SPN-C524i) can bind Anti-SARS-CoV-2 Spike Antibody, Mouse IgG1 (2C9C9) (XBB.1.5/Omicron Specific) (Cat. No. SPN-Y169) with a linear range of 0.391-25 ng/mL (QC tested). The antibody does not bind Spike Protein of WT (Cat. No. SPN-C52H9), Alpha (Cat. No. SPN-C52H6), Beta, Gamma (Cat. No. SPN-C52Hg), Delta (Cat. No. SPN-C52H6), Beta, B.1.1.529/Omicron (Cat. No. SPN-C52Hz), BA.2/Omicron (Cat. No. SPN-C5223), BA.3/Omicron (Cat. No. SPN-C5225), BA.4/Omicron (Cat. No. SPN-C5229), BA.5/Omicron (Cat. No. SPN-C522e), BA.2.12.1/Omicron (Cat. No. SPN-C522d), BQ.1.1/Omicron (Cat. No. SPN-C522s), BA.2.75/Omicron (Cat. No. SPN-C522f), BA.4.6/Omicron (Cat. No. SPN-C522m), BF.7/Omicron (Cat. No. SPN-C522q) (QC tested).



>>> www.acrobiosystems.com

7/24/2024

Catalog # SPN-Y169



Background

It's been reported that SARS-CoV-2 can infect the human respiratory epithelial cells through interaction with the human ACE2 receptor. The spike protein is a large type I transmembrane protein containing two subunits, S1 and S2. S1 mainly contains a receptor binding domain (RBD), which is responsible for recognizing the cell surface receptor. S2 contains basic elements needed for the membrane fusion. The S protein plays key parts in the induction of neutralizing-antibody and T-cell responses, as well as protective immunity.

Clinical and Translational Updates



>>> www.acrobiosystems.com

