

DiagSupport[™] Chloromethyl Polystyrene Resin, 200-400 Mesh, 0.8-1.4 mmol/g

Cat.No: SPS-RA23-029

DESCRIPTION

Description The chloromethyl polystyrene resin (Merrifield resin) has in the past been the

standard support for the synthesis of peptide acids by Boc strategy. Originally, the cesium salt of a protected amino acid was anchored to the chloromethyl support via nucleophilic displacement of chlorine. Although, Me₄N salts, sodium salts in THF with Bu₄NF catalysis and more recently zinc salts in EtOH have also been used. Cleavage is normally effected by treatment with HF or trifluoromethanesulfonic acid or by hydrogenolysis, alcohols can be released by using diisobutylaluminium hydride

or LiBH₄. For cleavage scavengers may be required.

APPLICATION

Application Notes Standard support for the synthesis of peptide acids by Boc strategy.

PRODUCT INFORMATION

Particle Size 200-400 mesh

Functional Group Chloromethyl

Capacity 0.8-1.4 mmol/g

45-1 Ramsey Road, Shirley, NY 11967, USA

Email: info@cd-bioparticles.com

Tel: 1-631-633-6938 Fax: 1-631-938-8221