

Veterinary Rapid Test

PROVIDE UNIQUE SOLUTION TO WORLD VETERINARY DIAGNOSIS

QUALITY

INTEGRITY

RESPONSIBILITY TO SERVING SOCIETY

HEALTHY WORLD



CPV Ag

Canine Parvo Virus Antigen Test



What is CPV ?

Canine parvovirus is a highly contagious disease caused by canine parvovirus (CPV). Dogs usually die within a short period of time.


Transmission route: fecal-oral route.

Infection symptom: vomiting, diarrhea and even bleeding. High morbidity, strong infectivity, and high mortality.

Hemorrhagic enteritis type: severe vomiting, hemorrhagic enteritis and significant decrease in white blood cells;

Myocarditis type: sudden death from myocarditis. vomiting, diarrhea and even bleeding.

Features






 Easy operation

 Fast result


 High sensitivity and accuracy

 Don't need the machine

General Information

Catalog no.	109103	Applicable location  Pet shop  Pet hospital  Pet beauty  Animal disease control center  Third-party testing institution
Product name	CPV Ag	
Purpose	Qualitative detection of Canine Parvo Virus Antigen in canine feces.	
Materials provided	*Test device *Swab *Buffer *Package insert *Workstation	
Animal	Canine	
Specimen	Feces	
Test time	10 minutes	
Specification	10pcs/box or 20pcs/box	
Storage	4-30°C	
Shelf life	24 months	

Test Procedure



CDV Ag

Canine Distemper Virus Antigen Test




What is CDV ?


Canine distemper virus is the oldest single stranded RNA virus infecting dogs. Canine distemper is a severe infectious disease caused by Canine distemper virus (CDV). The symptoms of simple viral infection are mild, but if the symptoms of secondary bacterial infection worsen, the mortality rate increases.

Transmission route: Through the digestive tract and respiratory tract.

Infection symptom: Early symptoms are similar to colds, followed by bronchitis, catarrhal pneumonia, gastroenteritis, and later spasms and convulsions.

Features






 Easy operation

 Fast result


 High sensitivity and accuracy

 Don't need the machine

General Information

Catalog no.	109101	Applicable location  Pet shop  Pet hospital  Pet beauty  Animal disease control center  Third-party testing institution
Product name	CDV Ag	
Purpose	Qualitative detection of Canine Distemper Virus Antigen in canine secretions.	
Materials provided	*Test device *Swab *Buffer *Package insert *Workstation	
Animal	Canine	
Specimen	Secretions	
Test time	10 minutes	
Specification	10pcs/box or 20pcs/box	
Storage	4-30 °C	
Shelf life	24 months	

Test Procedure



CCV Ag

Canine Coronavirus Antigen Test




What is CCV ?

Canine coronavirus often co-infects with canine parvovirus, rotavirus, and other gastrointestinal diseases. The mortality rate among puppies is high. It is a viral infectious disease that poses a serious threat to the dog breeding industry, livestock farming, and wildlife conservation.

Transmission route: Respiratory tract, digestive tract, feces, and contamination transmission.

Infection symptom: Canine coronavirus can cause gastroenteritis, with symptoms including vomiting and diarrhea in dogs. Seriously affected dogs may exhibit depression, excessive sleepiness, decreased or loss of appetite, and typically have no change in body temperature. This disease can occur throughout the year, but is more common in winter.

Features






 Easy operation

 Fast result


 High sensitivity and accuracy

 Don't need the machine

General Information

Catalog no.	109106	Applicable location  Pet shop  Pet hospital  Pet beauty  Animal disease control center  Third-party testing institution
Product name	CCV Ag	
Purpose	Qualitative detection of Canine Coronavirus Antigen in canine feces.	
Materials provided	*Test device *Swab *Buffer *Package insert *Workstation	
Animal	Canine	
Specimen	Feces	
Test time	10 minutes	
Specification	10pcs/box or 20pcs/box	
Storage	4-30°C	
Shelf life	24 months	

Test Procedure



FPV Ag

Feline Panleukopenia Virus Antigen Test




What is FPV ?

Feline panleukopenia is caused by feline distemper virus. It's a highly infectious and lethal feline viral disease, which spreads all over the world.

Transmission route: Transmitted through fecal orifice. Indirect contact is the most common route of infection, and FPV can be carried by pollutants (shoes, clothing), which means indoor cats also have a risk of infection.

Infection symptom: High fever, stubborn vomiting, diarrhea, dehydration, circulatory disorders, and rapid reduction of white blood cells. Cats that are not fully vaccinated or not vaccinated are more likely to develop feline fever, especially kittens aged 3-5 months.

Features






 Easy operation

 Fast result


 High sensitivity and accuracy

 Don't need the machine

General Information

Catalog no.	109125	Applicable location  Pet shop  Pet hospital  Pet beauty  Animal disease control center  Third-party testing institution
Product name	FPV Ag	
Purpose	Qualitative detection of Feline Panleukopenia Virus Antigen in feline feces.	
Materials provided	*Test device *Swab *Buffer *Package insert *Workstation	
Animal	Feline	
Specimen	Feces	
Test time	10 minutes	
Specification	10pcs/box or 20pcs/box	
Storage	4-30 °C	
Shelf life	24 months	

Test Procedure



FIPV Ab

Feline Infectious Peritonitis Virus Antibody Test



What is FIPV ?

Feline Infectious Peritonitis Virus (FIPV) is an important disease of domestic cats. It occurs worldwide in cats of all ages, but the disease is most common in young cats less than two years of age.


Transmission route: Direct contact, fecal-oral route.

Infection symptom:

In the wet (or effusive) form: Causes bloating and swelling in the abdomen (ascites) and may affect the heart and lungs. Cats with this form may pant and act sleepy or lethargic.

In the dry (or non-effusive) form: Usually affects the eyes and may have neurological symptoms such as trouble with balance and seizures.

Features






 Easy operation

 Fast result


 High sensitivity and accuracy

 Don't need the machine

General Information

Catalog no.	109126	Applicable location  Pet shop  Pet hospital  Pet beauty  Animal disease control center  Third-party testing institution
Product name	FIPV Ab	
Purpose	Qualitative detection of Feline Infectious Peritonitis Virus Antigen in feline feces.	
Materials provided	*Test device *Swab *Buffer *Package insert *Workstation	
Animal	Feline	
Specimen	Whole Blood/Serum/Plasma	
Test time	10 minutes	
Specification	10pcs/box or 20pcs/box	
Storage	4-30°C	
Shelf life	24 months	

Test Procedure



FIV Ab/FeLV Ag

Feline Immunodeficiency Virus Antibody/ Feline Leukemia Virus Antigen Combo Test



What are FIV and FeLV ?

FIV has highly variable clinical signs that are similar to FeLV and include fever, anemia, weight loss. FeLV breaks down a cat's immune system, leading to symptoms that include pale gums, yellow color of the mouth and whites of the eyes, weight loss, and poor coat condition.

Transmission route: Saliva, nasal secretions, urine, feces, and milk and bite wounds.

Infection symptom:

FeLV(May be none or non-specific):depression,weakness,enlarged lymph nodes,anemia, lymphoma (cancer of lymphocyte blood cells),blood in stool,chronic or recurring infections, weight loss,diarrhea,excessive drinking and urination.

FIV:may be none or signs of immune deficiency causing infections in any organ system, cancer,or neurologic issues .

Features

Easy operation

Fast result


High sensitivity and accuracy

Don't need the machine

General Information

Catalog no.	109155	Applicable location Pet shop Pet hospital Pet beauty Animal disease control center Third-party testing institution
Product name	FIV Ab/FeLV Ag	
Purpose	Qualitative detection of FIV Antibody/FeLV Antigen Combo Test in feline whole/serum/plasma.	
Materials provided	*Test device *Buffer *Package insert *Workstation	
Animal	Feline	
Specimen	Whole Blood/Serum/Plasma	
Test time	10 minutes	
Specification	10pcs/box or 20pcs/box	
Storage	4-30 C	
Shelf life	24 months	

Test Procedure





Product name	Catalog No.	Abbr.	Specimen	Format	Specification
Canine Distemper Virus Antigen Test	109101	★ CDV Ag	Secretions	Cassette	10T/20T
Canine Distemper Virus Antibody Test	109102	CDV Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Canine Parvo Virus Antigen Test	109103	★ CPV Ag	Feces	Cassette	10T/20T
Canine Parvo Virus Antibody Test	109104	CPV Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Canine Influenza Virus Antigen Test	109105	CIV Ag	Secretions	Cassette	10T/20T
Canine Coronavirus Antigen Test	109106	★ CCV Ag	Feces	Cassette	10T/20T
Canine Parainfluenza Virus Antigen Test	109107	CPIV Ag	Secretions	Cassette	10T/20T
Canine Adenovirus II Antigen Test	109108	CAV-II Ag	Secretions	Cassette	10T/20T
Canine Adenovirus I Antigen Test	109109	CAV-I Ag	Secretions	Cassette	10T/20T
Canine Toxoplasma Antibody Test	109111	★ TOXO Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Canine Heartworm Antigen Test	109112	★ CHW Ag	Whole Blood/Serum/Plasma	Cassette	10T/20T
Leishmania Canine Antibody Test	109113	LSH Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Canine Brucella Antibody Test	109114	C.Bru Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Ehrlichia Canine Antibody Test	109115	★ EHR Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Canine Leptospirosis Antibody Test	109116	★ Lepto Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Babesia gibsoni Antibody Test	109117	★ BG Ab	Secretions	Cassette	10T/20T
Rabies Virus Antigen Test	109118	★ Rabies Ag	Secretions	Cassette	10T/20T
Rabies Virus Antibody Test	109119	Rabies Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Lyme Disease Antibody Test	109120	Lyme Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Canis Pregnancy Relaxin Test	109121	★ RLN	Whole Blood/Serum/Plasma	Cassette	10T/20T
Canine Giardia Antigen Test	109122	C-GIA Ag	Feces	Cassette	10T/20T
Canine Parvo/Corona Virus Antigen Combo Test	109124	★ CPV/CCV Ag	Feces	Cassette	10T/20T
Canine Distemper/ Parainfluenza Virus Antigen combo Test	109123	CDV/CPIV Ag	Secretions	Cassette	10T/20T
Canine Anaplasma Antibody Test	109137	C.ANA Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Canine Rotavirus Antigen Test	109138	ROTA Ag	Feces	Cassette	10T/20T
Canine Parvo/Distemper Antibody Combo Test	109139	CPV/CDV Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Canine Distemper/Adeno Virus Antigen Combo Test	109140	CDV/CAV Ag	Secretions	Cassette	10T/20T
Canine Parvo/Corona/Rota Virus Antigen Combo Test	109141	CPV/CCV/Rota Ag	Feces	Cassette	10T/20T
Canine Parvo/Corona/Giardia Antigen Combo Test	109142	★ CPV/CCV/Giardia Ag	Feces	Cassette	10T/20T
Canine Distemper/Adeno/Influenza Virus Antigen Combo Test	109143	CDV/CAV/CIV Ag	Secretions	Cassette	10T/20T
Canine Ehrlichia/Anaplasma Antibody Combo Test	109145	EHR/ANA Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Ehrlichia/Lyme/Anaplasma Antibody Combo Test	109146	EHR/LYM/ANA Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Ehrlichia/Lyme/Anaplasma Antibody/Heartworm Antigen Combo Test	109147	EHR/LYM/ANA Ab/CHW Ag	Whole Blood/Serum/Plasma	Cassette	10T/20T
Ehrlichia/Babesia/Anaplasma Antibody Combo Test	109148	EHR/BAB/ANA Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Ehrlichia/Babesia/Anaplasma Antibody/Heartworm Antigen Combo Test	109149	EHR/BAB/ANA Ab/CHW Ag	Whole Blood/Serum/Plasma	Cassette	10T/20T

Feline

Product name	Catalog No.	Abbr.	Specimen	Format	Specification
Feline Panleukopenia Virus Antigen Test	109125	★ FPV Ag	Feces	Cassette	10T/20T
Feline Infectious Peritonitis Virus Antibody Test	109126	FIPV Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Feline Infectious Peritonitis Virus Antigen Test	109127	FIPV Ag	Feces	Cassette	10T/20T
Feline Coronavirus Antigen Test	109128	★ FCoV Ag	Feces	Cassette	10T/20T
Feline Leukemia Virus Antigen Test	109129	★ FeLV Ag	Whole Blood/Serum/Plasma	Cassette	10T/20T
Feline Immunodeficiency Virus Antibody Test	109130	★ FIV Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Feline Giardia Antigen Test	109131	F-GIA Ag	Feces	Cassette	10T/20T
Feline Anaplasma Antibody Test	109132	F.ANA Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Feline Toxoplasma Antibody Test	109133	F.TOXO Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Feline Calicivirus Virus Antigen Test	109134	FCV Ag	Secretions	Cassette	10T/20T
Feline Heartworm Antigen Test	109135	FHW Ag	Whole Blood/Serum/Plasma	Cassette	10T/20T
Feline Panleukopenia Virus Antibody Test	109152	★ FPV Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Feline Coronavirus Antibody Test	109153	★ FCoV Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Feline Herpes Virus Antigen Test	109154	★ FHV Ag	Secretions	Cassette	10T/20T
FIV Antibody/FeLV Antigen Combo Test	109155	FIV Ab/FeLV Ag	Whole Blood/Serum/Plasma	Cassette	10T/20T
Feline Herpes/ Calicivirus Antigen Combo Test	109156	FHV/FCV	Secretions	Cassette	10T/20T

Porcine

Product name	Catalog No.	Abbr.	Specimen	Format	Specification
African Swine Fever Virus Antibody Test	108926	★ ASFV Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
African Swine Fever Virus Antigen Test	108927	★ ASFV Ag	Whole Blood/Serum/Plasma	Cassette	10T/20T
Porcine Rotavirus Antigen Test	108901	PROTA Ag	Feces	Cassette	10T/20T
Swine Infectious Gastroenteritis Virus Antigen Test	108902	TGE Ag	Feces	Cassette	10T/20T
Porcine Epidemic Diarrhea anti-IgA Test	108903	PED IgA	Whole Blood/Serum/Plasma	Cassette	10T/20T
Porcine Circovirus Antibody Test	108904	PCV Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Porcine Trichinella spiralis Antibody Test	108905	PTS Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Classical Swine Fever Virus Antibody Test	108906	CSFV Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Porcine Pseudorabies -gE Antibody Test	108907	PRV gE Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Porcine Pseudorabies -gB Antibody Test	108908	PRV gB Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Porcine Reproductive and Respiratory Syndrome Antibody Test	108909	PRRS Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Swine Foot and Mouth Diseases Virus Serotype-O Antibody Test	108910	S.FMDV-O Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Swine Foot and Mouth Diseases Virus Serotype-A Antibody Test	108911	S.FMDV-A Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T

Bovine

Product name	Catalog No.	Abbr.	Specimen	Format	Specification
Bovine Foot and Mouth Diseases Virus Serotype-O Antibody Test	108917	B.FMDV-O Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Bovine Foot and Mouth Diseases Virus Serotype-A Antibody Test	108918	B.FMDV-A Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Bovine Brucella Antibody Test	108919	B.BRU Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Bovine Viral Diarrhea Antibody Test	108921	BVDV Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Bovine Infectious Rhinitis Antibody Test	108922	IBR Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T



Product name	Catalog No.	Abbr.	Specimen	Format	Specification
Newcastle Disease Virus Antigen Test	108912	NDV Ag	Cloacal secretions/Secretions	Cassette	10T/20T
Avian Influenza Virus Antigen Test	108913	★ AIV Ag	Cloacal secretions/Secretions	Cassette	10T/20T
Avian Influenza Virus H5 Antigen Test	108914	AIV H5 Ag	Cloacal secretions/Secretions	Cassette	10T/20T
Avian Influenza Virus H7 Antigen Test	108915	AIV H7 Ag	Cloacal secretions/Secretions	Cassette	10T/20T
Avian Influenza Virus H9 Antigen Test	108916	AIV H9 Ag	Cloacal secretions/Secretions	Cassette	10T/20T



Product name	Catalog No.	Abbr.	Specimen	Format	Specification
Sheep Brucella Antibody Test	108920	S.BRU Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Clostridium Perfringens Antibody Test	108923	CLP Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Clostridium Spoilage Antibody Test	108924	CLS Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Peste des Petits Ruminants Antibody Test	108925	PPR Ab	Whole Blood/Serum/Plasma	Cassette	10T/20T
Foot and Mouth Diseases Virus Non-structural Protein 3ABC Antibody Test	108928	FMDV (3ABC)	Whole Blood/Serum/Plasma	Cassette	10T/20T



Product name	Catalog No.	Abbr.	Specimen	Format	Specification
Vamber Canine Infectious Hepatitis/Parvovirus/Distemper Virus IgG Antibody Combo Test	109144	ICH/CPV/CDV	Whole Blood/Serum/Plasma	Elisa Key	10T/20T
Vamber Canine Lyme/Ehrlichia/Anaplasma/Leishmania/ Babesia IgG Antibody Combo Test	109160	LYM/EHR/ANA/LSH/BAB	Whole Blood/Serum/Plasma	Elisa Key	10T/20T
Vamber Canine Infectious Hepatitis/Parvovirus/Distemper Virus/Leptospira/Toxoplasma IgG Antibody Combo Test	109161	ICH/CPV/CDV/LEPTO/TOXO	Whole Blood/Serum/Plasma	Elisa Key	10T/20T
Vamber Canine Pancreatic Lipase Test	109162	CPL	Whole Blood/Serum/Plasma	Elisa Key	10T/20T
Vamber Feline Panleukopenia/Herpes/Calicivirus IgG Combo Test	109157	FPV/FHV/FCV	Whole Blood/Serum/Plasma	Elisa Key	10T/20T
Vamber Feline Panleukopenia/Herpes/Calicivirus/Coronavirus/Toxoplasma IgG Combo Test	109163	FPV/FHV/FCV/FCOV/TOXO	Whole Blood/Serum/Plasma	Elisa Key	10T/20T
Vamber Feline Pancreatic Lipase Test	109164	FPL	Whole Blood/Serum/Plasma	Elisa Key	10T/20T



Rapid Test Kit

Hangzhou Testsea Biotechnology Co., Ltd.

Add: No. 13-2 Guanshan Road, Yuhang District, 311115
Hangzhou, Zhejiang, China

Web: www.testsealabs.com

Email: sales@testsealabs.com

Tel: +860571-88600160

Rabies Ag

Rabies Virus Antigen Test




What is RBV ?

Rabies virus (RBV) belongs to the genus Rabies virus of Rhabdoviridae. It is bullet shaped and consists of a spiral symmetric nucleocapsid, envelope, and single stranded RNA. Mainly present in the salivary glands and saliva of diseased animals. Rabies is a zoonotic disease caused by Rabies virus, with a high mortality rate. Once rabies shows symptoms, the mortality rate can reach 100%.

Transmission route: It is usually transmitted from diseased animals to humans and animals through saliva by biting.

Infection symptom: Drooling, throat paralysis, hind limb paralysis, and inability to walk normally. Due to the frequent clinical manifestations of hydrophobia, it is also known as hydrophobia.

Features






 Easy operation

 Fast result


 High sensitivity and accuracy

 Don't need the machine

General Information

Catalog no.	109118	Applicable location  Pet shop  Pet hospital  Pet beauty  Animal disease control center  Third-party testing institution
Product name	Rabies Ag	
Purpose	Qualitative detection of Rabies Virus Antigen in canine secretions.	
Materials provided	*Test device *Swab *Buffer *Package insert *Workstation	
Animal	Canine	
Specimen	Secretions(saliva)	
Test time	10 minutes	
Specification	10pcs/box or 20pcs/box	
Storage	4-30 C	
Shelf life	24 months	

Test Procedure



FCoV Ag

Feline Coronavirus Antigen Test




What is FCoV Ag ?


Feline coronavirus (FCoV) is a positive-sense single-stranded RNA virus that is infectious to all species of feline animals worldwide. It exists in two different forms: Feline Enteric Coronavirus (FECV), which infects the intestines, and Feline Infectious Peritonitis Virus (FIPV), which causes Feline Infectious Peritonitis (FIP).

Transmission route: Fecal-oral Route


Infection symptom: Vomiting and moderate diarrhea may occur after infection, which can lead to dehydration if severe. The mortality rate is low, and the recovery rate is high. Feline enteric coronavirus only shows symptoms in the initial infection, but once infected, it carries the virus for life, with no apparent clinical manifestations afterwards.

Features






 Easy operation

 Fast result


 High sensitivity and accuracy

 Don't need the machine

General Information

Catalog no.	109128	Applicable location  Pet shop  Pet hospital  Pet beauty  Animal disease control center  Third-party testing institution
Product name	FCoV Ag	
Purpose	Qualitative detection of Feline Coronavirus Antibody in feline feces.	
Materials provided	*Test device *Swab *Buffer *Package insert *Workstation	
Animal	Feline	
Specimen	Feces	
Test time	10 minutes	
Specification	10pcs/box or 20pcs/box	
Storage	4-30°C	
Shelf life	24 months	

Test Procedure



FHV Ag

Feline Herpes Virus Antigen Test




What is FHV Ag ?


Feline Herpesvirus (FHV) is a virus that primarily causes acute upper respiratory infections in cats.

Transmission route: Direct contact (through contact with saliva, secretions from the eyes or nasal cavity); sneezing; sharing food bowls and contaminated environment (including shared items, beauty tools).

Infection symptom: Coughing, sneezing, runny nose with watery or purulent discharge, conjunctivitis, fever, reduced appetite, and swollen lymph nodes mainly affecting the respiratory tract.

Features






 Easy operation

 Fast result

 High sensitivity and accuracy

 Don't need the machine

General Information

Catalog no.	109154	Applicable location  Pet shop  Pet hospital  Pet beauty  Animal disease control center  Third-party testing institution
Product name	FHV Ag	
Purpose	Qualitative detection of Feline Panleukopenia Virus Antigen in feline feces.	
Materials provided	*Test device *Swab *Buffer *Package insert *Workstation	
Animal	Feline	
Specimen	Secretions	
Test time	10 minutes	
Specification	10pcs/box or 20pcs/box	
Storage	4-30°C	
Shelf life	24 months	

Test Procedure

