

Comprehensive Solution for Women's Fertility Health

Accompanying women at every stage

Each stage of a woman's life is unique and significant.
As a reliable partner, Testsealabs cares for women's health,
provides diagnostic products for women at every stage.



Hangzhou Testsea Biotechnology Co., Ltd.

sales@testsealabs.com

www.testsealabs.com

No. 13-2 Guanshan Road, Yuhang District, 311115 Hangzhou, Zhejiang, China



ADVANTAGE

- High sensitivity and high specificity
- User-friendly
- Various product forms are available
- Short inspection time

INTENDED USE




Prepare for pregnancy



Test for early pregnancy



Monitoring for Fetal health



Detection for menopause

Digital Ultra Early Pregnancy Test with Week Indicator

HCG Pregnancy Test

Clear, digital results: give you clear digital "Pregnant" or "Not Pregnant" result in Words
A unique flashing "stop light"



ADVANTAGE

- One step test
- Clear results in 4 min.
- High accuracy: > 99%
- Various product forms (strip/device/panel/cup/midstream/digital reader)
- Multiple sensitivity options


LH Ovulation Test

The Best Ovulation Tests We Tested at HOME

Pinpoint your fertile window and get pregnant faster with these easy-to-use ovulation tests

--- ADVANCE DIGITAL Ovulation Test --Identify Your Most Fertile Days, Find your 2 peak days, double your chances of getting pregnant

*The fertile days of the menstrual period start after this LH surge.



PRODUCT INFORMATION




| Combination set | |
|--------------------|---|
| Single Kit | 1 reader + HCG test strip |
| Ovulation Kit | 1 reader + 10 LH test strips |
| Combination | 1 reader + 10 LH test strips + 2 HCG test strips |
| Pregnancy Test Kit | 2 early pregnancy test sticks + 1 electronic pregnancy test stick |

Operation process

a) Either.

Put just the absorbent tip pointing downwards into your urine stream for 1-5 seconds. Take care not to get the reader wet.



Or

Collect a sample of your urine in a clean, dry container and put just the absorbent tip into the sample for 1-10 seconds.

b) Replace the cap, wipe off any excess urine. Keep the Reader pointing lying flat.

| | |
|-----------------------------|---|
| Not Pregnant X | Display "Not Pregnant" "X" for HCG negative result. |
| Pregnant 1-2 W | Display "Pregnant" "1-2" "W" for HCG positive result. |
| Pregnant 2-3 W | Display "Pregnant" "2-3" "W" for HCG positive result. |
| Pregnant 3+W | Display "Pregnant" "3+" "W" for HCG positive result. |
| No Ovulatory ↻ | Display "No Ovulatory" "↻" for LH negative result. |
| Ovulatory ☾ | Display "Ovulatory" "☾" for LH positive result. A rise in estrogen has been detected, increased chance of getting pregnant. |
| Ovulatory ☽ | Display "Ovulatory" "☽" for LH positive result: highest chance of getting pregnant. Peak ovulation usually lasts for 24-48 hours. |
| X | Device error for both HCG and LH |

| | Product name | Catalog No. | Specimen | Format | Specification | Certificate |
|--|--|-------------|-------------|-------------------|----------------|-------------|
| Convention Fertility Test Series | HCG Pregnancy Test (Urine) | 109702 | Urine | Strip | 50T | CE ISO |
| | | | | Cassette | 25T/40T | |
| | HCG Pregnancy Test (Serum) | 109702 | Serum | Midstream | 1T/20T | CE ISO |
| | | | | Cassette | 25T/40T | |
| | HCG Pregnancy Test (Serum/Urine) | 109702 | Serum/Urine | Cassette | 25T/40T | CE ISO |
| | | | | LH Ovulation Test | 109703 | Urine |
| Digital Fertility Test Series | FSH Menopause Test | 109704 | Urine | Strip | 50T | CE ISO |
| | | | | Cassette | 25T/40T | |
| | Digital HCG Pregnancy Test | 109702 | Urine | Midstream | 1T/2T/3T/5T | ISO |
| | Digital LH Ovulation Test | 109703 | Urine | Midstream | 7T/10T/15T/20T | ISO |
| | Digital Pregnancy & Ovulation Combination Test Set | 109700 | Urine | Midstream | customized | ISO |

Torch Test

Toxo, RV, CMV, HSV I / II

The "Five Prenatal Tests," also known as the "Five Teratogenic Infections," refer to five pathogens that, if a mother becomes infected, may pose significant risks to the fetus.


Pre-pregnancy Testing: Assess the immunity of women of childbearing age.

Prenatal Testing: Compare with pre-pregnancy status to determine infection conditions and evaluate the necessity of prenatal diagnosis for intrauterine infections.

Postnatal Testing: Conduct tests on the newborn to provide evidence of congenital infections after birth.

FEATURES & BENEFITS

- Specimen: Whole blood/Serum/Plasma
- Rapid results in 15 min.
- One step test





FSH Menopause Test

Qualitative detection of follicle-stimulating hormone (FSH) levels in human urine is primarily used clinically as an aid in the diagnosis of ovarian disorders.

WHEN TO USE FSH TEST?

- **Sample Type:** Urine
- **Formats Available:** Strip and Cassette
- **High Sensitivity:** Accurately detects ovarian decline
- **Rapid Results:** Observe results within 10 minutes
- **Wide Applications:** Suitable for promoting eugenics, preventing premature ovarian failure, addressing postnatal decline, and predicting menopause





TSH Thyroid Stimulating Hormone Test

The Thyroid-Stimulating Hormone (TSH) Test Kit is intended for the quantitative in vitro determination of TSH levels in human serum or plasma.

Sample Type: Serum/Plasma ●

Detection Time: Results readable within 5 minutes ●

Sensitivity: 98.72% ●

Specificity: 98.36% ●

HPV Test Seris

PEOPLE SUSCEPTIBLE TO HPV INFECTION

- Multiple sex partners
- Homosexuals
- Individuals with weakened immunity due to staying up late
- History of HPV infection



| | Product Name | Catalog No. | Specimen | Format | Specification | Certificate |
|----------------|---|-------------|-----------------|-----------|---------------|-------------|
| HPV Test Seris | Human Papillomavirus (HPV) Test Midstream | 101105 | Urine | Midstream | IT/20T | CE ISO |
| | HPV 16/18 E7 Triline Antigen Test Cassette | 101109 | Female Cervical | Cassette | IT/25T | CE ISO |
| | HPV 16+18 E7 Antigen Test | 101110 | Female Cervical | Cassette | IT/25T | ISO |
| | HPV L1+16/18 E7 Antigen Combo Test | 101112 | Female Cervical | Cassette | IT/25T | ISO |
| | HPV 16/18+L1 Combo Antigen Test Cassette | 101116 | Female Cervical | Cassette | IT/25T | CE ISO |
| | Medical anti-HPV Functional Protein Gynecological Gel | 101135 | / | / | IT | ISO |

Vaginitis Test Series

PRODUCT FEATURES



Fungal Vaginitis

* Candida Albicans



Trichomonas

*Trichomonas Vaginitis



Bacterial Vaginitis

*Gardnerella Vaginalis



| | Product name | Catalog No. | Specimen | Format | Specification | Certificate |
|-----------------------|--|-------------|---------------------|----------|---------------|-------------|
| Vaginitis Test Series | Vaginits Multi-test Kit (Dry Chemoenzymatic Method) | 101136 | Urine | Cassette | IT | ISO |
| | Candida Albicans+Trichomonas Vaginalis Antigen Combo Test Cassette | 101133 | Female Vaginal Swab | Cassette | IT/25T | ISO |
| | Candida Albicans+Trichomonas Vaginalis+Gardnerella Vaginalis Antigen Combo Test Cassette | 101132 | Female Vaginal Swab | Cassette | IT/25T | ISO |

Disease Series

Sexually Transmitted Disease Series

Symptom



Chlamydia

Chlamydia



Neisseria

Gonorrhoeae

HIV 1/2

Syphilis


1 drop of specimen
2 drops of buffer



Positive




Negative



Invalid

| | Product name | Catalog No. | Specimen | Format | Specification | Certificate |
|-------------------------------------|--|-------------|-------------------------------------|----------------|---------------|-------------|
| Sexually Transmitted Disease Series | Chlamydia Trachomatis Ag Test | 101038 | Male Urethral Swab/ Female Cervical | Cassette | IT/25T | ISO |
| | Neisseria Gonorrhoeae Ag Test | 101039 | Male Urethral Swab/ Female Cervical | Cassette | IT/25T | CE ISO |
| | Chlamydia+Gonorrhoeae Antigen Combo Test | 101101 | Male Urethral Swab/ Female Cervical | Cassette | IT/25T | ISO |
| | HIV 1/2 Triline Test | 101008 | WB/S/P | Strip/Cassette | 25T/40T/50T | ISO |
| | Syphilis (Anti-treponema Pallidum) Antibody Test | 101015 | WB/S/P | Strip/Cassette | 25T/40T/50T | CE ISO |

Other Disease Test Series



| | Product name | Catalog No. | Specimen | Format | Specification | Certificate |
|---------------------------|---|-------------|-----------------------|----------|---------------|-------------|
| Other Disease Test Series | ToRCH IgG/IgM Test Cassette (Toxo, RV, CMV, HSV I / II) | 109701 | WB/S/P | Cassette | 25T | ISO |
| | Herpes Simplex Virus I Antibody IgG/IgM Test | 101046 | WB/S/P | Cassette | 25T/40T | CE ISO |
| | Herpes Simplex Virus II Antibody IgG/IgM Test | 101047 | WB/S/P | Cassette | 25T/40T | CE ISO |
| | Herpes Simplex Virus I/II Antibody IgG/IgM Test | 101074 | WB/S/P | Cassette | 25T/40T | CE ISO |
| | TSH Thyroid Stimulating Hormone | 109705 | Serum/Plasma | Cassette | 25T/40T | CE ISO |
| | Strep B test | 101134 | Vaginal / Rectal Swab | Cassette | 25T/40T | ISO |
| | IGFBP - 1 (PROM) TEST | 109908 | vaginal secretions | Cassette | 25T/40T | ISO |

PRODUCT CATALOGUE

| | Product name | Catalog No. | Specimen | Format | Specification | Certificate |
|-------------------------------------|--|-------------|------------------------------------|----------------|----------------|-------------|
| Convention Fertility Test Series | HCG Pregnancy Test (Urine) | 109702 | Urine | Strip | 50T | CE ISO |
| | | | | Cassette | 25T/40T | |
| | | | | Midstream | 1T/20T | |
| | HCG Pregnancy Test (Serum) | 109702 | Serum | Cassette | 25T/40T | CE ISO |
| | HCG Pregnancy Test (Serum/Urine) | 109702 | Serum/Urine | Cassette | 25T/40T | CE ISO |
| Digital Fertility Test Series | LH Ovulation Test | 109703 | Urine | Strip | 50T | CE ISO |
| | | | | Cassette | 25T/40T | |
| | | | | Midstream | 1T/20T | |
| | FSH Menopause Test | 109704 | Urine | Strip | 50T | CE ISO |
| | | | | Cassette | 25T/40T | |
| Digital Fertility Test Series | Digital HCG Pregnancy Test | 109702 | Urine | Midstream | 1T/2T/3T/5T | ISO |
| | Digital LH Ovulation Test | 109703 | Urine | Midstream | 7T/10T/15T/20T | ISO |
| | Digital Pregnancy & Ovulation Combination Test Set | 109700 | Urine | Midstream | customized | ISO |
| HPV Test Series | Human Papillomavirus(HPV) Test Midstream | 101105 | Urine | Midstream | 1T/20T | CE ISO |
| | HPV 16/18 E7 Triline Antigen Test Cassette | 101109 | Female Cervical | Cassette | 1T/25T | CE ISO |
| | HPV 16+18 E7 Antigen Test | 101110 | Female Cervical | Cassette | 1T/25T | ISO |
| | HPV L1+16/18 E7 Antigen Combo Test | 101112 | Female Cervical | Cassette | 1T/25T | ISO |
| | HPV 16/18+L1 Combo Antigen Test Cassette | 101116 | Female Cervical | Cassette | 1T/25T | CE ISO |
| | Medical anti-HPV Functional Protein Gynecological Gel | 101135 | / | / | 1T | ISO |
| Vaginitis Test Series | Vaginits Multi-test Kit (Dry Chemoenzymatic Method) | 101136 | Urine | Cassette | 1T | ISO |
| | Candida Albicans+Trichomonas Vaginalis Antigen Combo Test Cassette | 101133 | Female Vaginal Swab | Cassette | 1T/25T | ISO |
| | Candida Albicans+Trichomonas Vaginalis+Gardnerella Vaginalis Antigen Combo Test Cassette | 101132 | Female Vaginal Swab | Cassette | 1T/25T | ISO |
| | Chlamydia Trachomatis Ag Test | 101038 | Male Urethral Swab/Female Cervical | Cassette | 1T/25T | ISO |
| Sexually Transmitted Disease Series | Neisseria Gonorrhoeae Ag Test | 101039 | Male Urethral Swab/Female Cervical | Cassette | 1T/25T | CE ISO |
| | Chlamydia+Gonorrhoeae Antigen Combo Test | 101101 | Male Urethral Swab/Female Cervical | Cassette | 1T/25T | ISO |
| | HIV 1/2 Triline Test | 101008 | WB/S/P | Strip/Cassette | 25T/40T/50T | ISO |
| Other Disease Test Series | Syphilis(Anti-treponemia Pallidum) Antibody Test | 101015 | WB/S/P | Strip/Cassette | 25T/40T/50T | CE ISO |
| | ToRCH IgG/IgM Test Cassette (Toxo, RV, CMV, HSV I / II) | 109701 | WB/S/P | Cassette | 25T | ISO |
| | Herpes Simplex Virus I Antibody IgG/IgM Test | 101046 | WB/S/P | Cassette | 25T/40T | CE ISO |
| | Herpes Simplex Virus II Antibody IgG/IgM Test | 101047 | WB/S/P | Cassette | 25T/40T | CE ISO |
| | Herpes Simplex Virus I/II Antibody IgG/IgM Test | 101074 | WB/S/P | Cassette | 25T/40T | CE ISO |
| | TSH Thyroid Stimulating Hormone | 109705 | Serum/Plasma | Cassette | 25T/40T | CE ISO |
| | Strep B test | 101134 | Vaginal /Rectal Swab | Cassette | 25T/40T | ISO |
| | IGFBP - 1 (PROM) TEST | 109908 | vaginal secretions | Cassette | 25T/40T | ISO |

Hangzhou Testsea Biotechnology Co., Ltd.

Add: No. 13-2 Guanshan Road, Yuhang District, 311115
Hangzhou, Zhejiang, China

Web: www.testsealabs.com

Email: sales@testsealabs.com

Tel: +860571-88600160