



# Favorgen™ High Efficiency (HE) Mini Kit

Create an excellent result for your applications

## Favorgen™ HE Mini Kit – Unleash Efficiency, Elevate Purity!

### High Efficiency – Enhanced Outcomes, Less Effort!

- Fast and streamlined extraction process, economical with time and effort.
- Delivering high-yield nucleic acid even from the demanding samples.
- Available in multiple formats to perfectly match your sample type.

### Ultimate Purity – Reliable Results, You Can Trust!

- Advanced design of HE Columns eliminates buffer retention to avoid the risk of salt contamination and obtain highly concentrated nucleic acid.
- Distinct buffer system efficiently removes PCR inhibitors, producing high-purity nucleic acids for downstream applications.

Column	Favorgen HE Mini Kit	Conventional Mini Kit
Sample Input	1~5×	1×
Column Volume	1000 $\mu$ l	600~700 $\mu$ l
Elution Volume	$\leq$ 30 $\mu$ l	50~100 $\mu$ l
Typical Conc.	3~9×	1×
Typical Yield	1.5~6×	1×

## HE Column characteristic

