

# NEW

# High Efficiency Kit

Efficient manner to eliminate the inhibitors and extract high purity, high concentration and high yields genomic DNA from various tissues and blood samples.

## Feature

### EXTREME



### HE Column

- ✓ Enlarge **30%** sample input
- ✓ Minimum elution volume at **30 µl**
- ✓ Obtain **3~6 times** highly concentrated DNA

## Performance

### DNA Conc.



	Competitor Kit	FABG HE Mini Kit (FABG103)
Input Blood Volume (µl)	200	300 🏆
A260/A230	1.5	2.0 🏆
A260/A280	1.7	1.8 🏆
DNA Conc. (ng/µl)	28.9	166.2 🏆
DNA Yield (µg)	2.9	5.0 🏆

## Ordering Information

Kit	Cat. No.
FavorPrep™ Tissue Genomic DNA Extraction HE Mini Kit	FATG103-050 (50 preps) FATG103-100 (100 preps)
FavorPrep™ Blood Genomic DNA Extraction HE Mini Kit	FABG103-050 (50 preps) FABG103-100 (100 preps)

# FavorPrep™ Endotoxin-Free Plasmid Mini kit




## Features

- ✓ Streamlined endotoxin removal technology to efficiently remove residual endotoxin from plasmid DNA extraction procedure (<0.1 EU/µg DNA).
- ✓ Minimize cytotoxicity and increase transfection efficiency to obtain reproducible and reliable transfection results.

## Ordering Information

Kit	Cat. No.
FavorPrep™ Endotoxin-Free Plasmid Extraction Mini Kit	FSPD302-100 (100 preps)

### Endotoxin Level



FavorPrep™ Endotoxin-Free Plasmid Extraction Mini Kit (FSPD302) effectively removes endotoxin when processing plasmid extraction.



Favorgen Biotech Corporation  
 Tel: +886-8-7621829  
 Customer Service: [order@favorgen.com](mailto:order@favorgen.com)  
 Technical Service: [technical@favorgen.com](mailto:technical@favorgen.com)  
 For more product information, please visit [FAVORGEN.COM](http://FAVORGEN.COM)

