

XFD350 Anti-dog/ chicken/ rabbit/ guinea pig/ horse/ cow/ mouse/ rat/ pig/ non-human primates/ human CD79a Antibody *HM47, XFD350 Same Structure to Alexa Fluor™ 350*

Catalog number: 10790140, 10790141 Unit size: 100 tests, 500 tests

Product Details

Storage Conditions 2-8°C with minimized light exposure. Do not freeze.

Expiration Date 12 months upon receiving

Concentration 0.1 mg/mL

Formulation Phosphate-buffered saline (PBS, pH 7.2), 0.09% sodium azide, 0.2% (w/v) BSA

Antibody Properties

Species Reactivity Dog, chicken, rabbit, guinea pig, horse, cow, mouse, rat, pig, non-human primates, human

Class Primary

Clonality Monoclonal

Host Mouse

Isotype Mouse IgG1 kappa

Immunogen CD79a (Mb-1, IGA)

Clone HM47

Conjugate AF350

Biological Properties

Appearance Yellow liquid

Preparation Antibody purified by affinity chromatography and then conjugated with AF350 under optimal

conditions

Application Flow Cytometry (FACS), Fluorescence Imaging

Spectral Properties

Conjugate AF350

Excitation Wavelength 343 nm

Emission Wavelength 441 nm

Applications

HM47 is an anti-dog/ chicken/ rabbit/ guinea pig/ horse/ cow/ mouse/ rat/ pig/ non-human primates/ human monoclonal antibody that forms an immune complex with the CD79a antigen. CD79a (sometimes referred to as Mb-1 or IGA) is a 47 kD single-pass type I membrane protein that is expressed on the surface of cells such as B cells. CD79a is a component of essential cellular pathways, in particular, the B cell receptor signaling pathway. From a research standpoint, it is of biological interest due to its association with key macromolecules/ligands like CD22, CD79b, CD5 and CD19. CD79a is a fairly uncommon antibody target, with a little more than 4100 publications in the last decade. Even still, CD79a has been widely used in immunology research, commonly serving as a phenotypic marker for differentiating cell types in flow cytometric applications. This antibody was purified through affinity chromatography and conjugated to XFD350 (ex/em = 343/441 nm). XFD350 is manufactured by AAT Bioquest, and it has the same chemical structure of Alexa Fluor® 350 (Alexa Fluor® is the trademark of ThermoFisher). It is compatible with the 355 nm laser and 450/50 nm bandpass filter (for example, as in the BD Special Order LSRFortessa™ Cell Analyzer).