

XFD647 Anti-dog/ chicken/ rabbit/ guinea pig/ horse/ cow/ mouse/ rat/ pig/ non-human primates/ human CD79a Antibody
HM47, XFD647 Same Structure to Alexa Fluor™ 647

Catalog number: 10790180, 10790181
Unit size: 100 tests, 500 tests

Product Details

Storage Conditions	2-8°C with minimized light exposure. Do not freeze.
Expiration Date	12 months upon receiving
Concentration	0.1 mg/mL
Formulation	Phosphate-buffered saline (PBS, pH 7.2), 0.09% sodium azide, 0.2% (w/v) BSA

Antibody Properties

Species Reactivity	Dog, chicken, rabbit, guinea pig, horse, cow, mouse, rat, pig, non-human primates, human
Class	Primary
Clonality	Monoclonal
Host	Mouse
Isotype	Mouse IgG1 kappa
Immunogen	CD79a (Mb-1, IGA)
Clone	HM47
Conjugate	AF647

Biological Properties

Appearance	Dark blue liquid
Preparation	Antibody purified by affinity chromatography and then conjugated with AF647 under optimal conditions
Application	Flow Cytometry (FACS), Fluorescence Imaging

Spectral Properties

Conjugate	AF647
Excitation Wavelength	650 nm
Emission Wavelength	671 nm

Applications

HM47 is an anti-dog/ chicken/ rabbit/ guinea pig/ horse/ cow/ mouse/ rat/ pig/ non-human primates/ human monoclonal antibody that forms an immune complex with the CD79a antigen. CD79a (sometimes referred to as Mb-1 or IGA) is a 47 kD single-pass type I membrane protein that is expressed on the surface of cells such as B cells. CD79a is a component of essential cellular pathways, in particular, the B cell receptor signaling pathway. From a research standpoint, it is of biological interest due to its association with key macromolecules/ligands like CD22, CD79b, CD5 and CD19. CD79a is a fairly uncommon antibody target, with a little more than 4100 publications in the last decade. Even still, CD79a has been widely used in immunology research, commonly serving as a phenotypic marker for differentiating cell types in flow cytometric applications. This antibody was purified through affinity chromatography and conjugated to XFD647 (ex/em = 650/671 nm). XFD647 is manufactured by AAT Bioquest, and it has the same chemical structure of Alexa Fluor® 647 (Alexa Fluor® is the trademark of ThermoFisher). It is compatible with the 642 nm laser and 702/87 nm bandpass filter (for example, as in the Luminex Amnis CellStream).