

This product is for research use only (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

Product no AS14 2822

Anti-V-ATPase, a1 | vacuolar H+-ATPase subunit a isoform 1

Product information

Immunogen KLH-conjugated synthetic peptide derived from *Arabidopsis thaliana* V-ATPase subunit A, UniProt: <u>Q8RWZ7-1</u>, TAIR:

AT2G2852

Host Rabbit

Clonality Polyclonal

Purity Immunogen affinity purified serum in PBS pH 7.4.

Format Lyophilized

Quantity 50 μg

Reconstitution For reconstitution add 50 μl of sterile water

Storage Store lyophilized/reconstituted at -20 °C; once reconstituted make aliquots to avoid repeated freeze-thaw cycles. Please remember to spin the tubes briefly prior to opening them to avoid any losses that might occur from material adhering to

the cap or sides of the tube.

Application information

Recommended dilution 1:1000 (WB)

Expected | apparent 93 kDa (Arabidopsis thaliana)

Confirmed reactivity | Arabidopsis thaliana

MW 95 KDa (Arabidopsis trialiana


Predicted reactivity Capsella rubella, Cucumis sativus, Erythranthe guttata, Glycine soja, Lupinus angustifolius, Morus notabilis, Phaseolus


vulgaris , Vitis vinifera

Species of your interest not listed? Contact us

Not reactive in No confirmed exceptions from predicted reactivity are currently known

Additional information Protein or membrane sample should be treated at 70°C for 10 min before loading on the gel.





10 µg of microsomal membranes were isolated from 6-days-old etiolated (dark-grown) *Arabidopsis thaliana* seedlings with Extraction buffer (450mM sucrose, 50mM HEPES pH7.5, 5mM MgCl2, 1mM DTT, 1x protease inhibitor). Proteins were denaturated in SDS sample buffer for 5 min at 95°C. Proteins were separated on a 10% SDS-PAGE and blotted 1h to a nitrocellulose membrane using tank transfer. Blots were blocked with 4% dry milk in TBS-T for 2h at room temperature (RT) with agitation. Blot was incubated in the primary antibody (#AS14 2822) at a dilution of 1:1000 in blocking buffer for 16h with agitation. The antibody solution was decanted and the blot was rinsed briefly twice, then washed 3 times for 15 min in TBS-T at RT with agitation. Blot was incubated in secondary antibody (anti-rabbit IgG horse radish peroxidase conjugated, from Agrisera #AS09 602) diluted to 1:25 000 in blocking buffer for 2h at RT with agitation. The blot was washed as above and developed for 3 min with ECL according to the manufacturer's instructions peqlab AceGlow Kit & INTAS digital imager. Exposure time was 3 min.

Courtesy of Fabian Fink, University Heidelberg, Germany