

CKR-2 rabbit pAb

Cat No.:ES4394

For research use only

Overview

Product Name	CKR-2 rabbit pAb
Host species	Rabbit
Applications	WB;IF;ELISA
Species Cross-Reactivity	Mouse;Rat
Recommended dilutions	Western Blot: 1/500 - 1/2000. Immunofluorescence:
	1/200 - 1/1000. ELISA: 1/5000. Not yet tested in
	other applications.
Immunogen	Synthesized peptide derived from CKR-2 . at AA
	range: 250-330
Specificity	CKR-2 Polyclonal Antibody detects endogenous
	levels of CKR-2 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	C-C chemokine receptor type 2
Gene Name	CCR2
Cellular localization	
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	46kD
Human Gene ID	
Human Swiss-Prot Number	
Alternative Names	CCR2; CMKBR2; C-C chemokine receptor type 2; C-C
	CKR-2; CC-CKR-2; CCR-2; CCR2; Monocyte
	chemoattractant protein 1 receptor; MCP-1-R; CD
	antigen CD192
Background	CCR2 (C-C Motif Chemokine Receptor 2) is a Protein
	Coding gene. Diseases associated with CCR2 include
	idiopathic anterior uveitis and cd3zeta deficiency.
	Among its related pathways are Signaling by GPCR



+86-27-59760950

ELKbio@ELKbiotech.com


www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



and A-beta Pathways: uptake and Degradation. GO annotations related to this gene include G-protein coupled receptor activity and chemokine receptor activity. An important paralog of this gene is CCR9. eceptor for the CCL2, CCL7 and CCL13 chemokines. Transduces a signal by increasing intracellular calcium ion levels. Alternative coreceptor with CD4 for HIV-1 infection. hemokine CC receptors (CCRs) predominantly recognize CC chemokines. CC chemokines are distinguished by having four conserved cysteines, with the first two cysteines being adjacent to each other. There are 10 chemokine CC receptors.

Western Blot analysis of K562 cells using CKR-2 Polyclonal Antibody





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C