

DiagPoly™ Plain Fluorescent Polystyrene Particles, Red, 250 µm

Cat.No: DCFR-L170

DESCRIPTION

Description	DiagPoly™ Plain Fluorescent Polystyrene Particles, Red, 250 µm (# DCFR-L170) are encapsulated with red fluorophores (excitation: 546 nm, emission: 566 nm) and there is no additional coatings or special functional groups on the surface. The density of liquid suspension is 1.05-1.06 g/mL. Sonicate/Vortex the microspheres before each use.
--------------------	---

APPLICATION

Application Notes	Used as standards or markers, in immunoassays (such as latex agglutination assays and solid phase enzyme immunoassays), or coated with proteins for use in research and diagnostic tests.
--------------------------	---

PRODUCT INFORMATION

Diameter	250 µm
Functional Group	Plain
Concentration	1% w/v
Absorption Max	545 nm
Emission Max	566 nm

STORAGE AND SHIPPING

Storage Buffer	De-ionized water containing a small amount of surfactant and 2mM of sodium azide as an anti-microbial agent.
Storage	Store at 2-8°C. Keep refrigerated. Do not freeze. Keep container closed and fluorescent particles protected from light.