

## DiagPoly™ Hydrazide Polystyrene Particles, 20 µm

Cat.No: DNM-M009

### DESCRIPTION

<b>Description</b>	Hydrazide PS Particles, 20 µm, 25 mg/mL, dispersed in DI water.
--------------------	---

### PRODUCT INFORMATION

<b>Diameter</b>	20 µm
<b>Functional Group</b>	Hydrazide
<b>Concentration</b>	25 mg/mL
<b>Number of Particles</b>	5.8E+6 particles/mL
<b>Size Distribution</b>	< 5 % (C.V.)

### STORAGE AND SHIPPING

<b>Storage Buffer</b>	Suspension in water. Stable in aqueous buffers methanol, ethanol, DMSO. Not stable in halogenated hydrocarbons, toluene, strong acidic solutions, e.g. 10% HCl.
<b>Storage</b>	2-8°C and do not freeze!
<b>Shelf Life</b>	6 months