

## DiagPoly™ Epoxy Polystyrene Particles, 30 μm

Cat.No: DNM-M001

## **DESCRIPTION**

**Description** Epoxy PS Particles, 30 μm, 25 mg/mL, dispersed in DI water.

## PRODUCT INFORMATION

| Diameter            | 30 μm               |
|---------------------|---------------------|
| Functional Group    | Ероху               |
| Concentration       | 25 mg/mL            |
| Number of Particles | 1.7E+6 particles/mL |
| Size Distribution   | < 5 % (C.V.)        |

## STORAGE AND SHIPPING

| Storage Buffer | Suspension in water. Stable in aqueous buffers methanol, ethanol, DMSO. Not stable in halogenated hydrocarbons, toluene, strong acidic solutions, e.g. 10% HCl. |
|----------------|---|
| Storage        | 2-8°C and do not freeze!  |
| Shelf Life     | 2 months  |

45-1 Ramsey Road, Shirley, NY 11967, USA

Email: info@cd-bioparticles.com

Tel: 1-631-633-6938 Fax: 1-631-938-8221 © CD Bioparticles All Rights Reserved