

DiagNano™ Anti-Rabbit IgG conjugated Gold Nanorods, diameter 40 nm, absorption max 750 nm

Cat.No: RRG-40-750

DESCRIPTION

Description 50 OD Anti-Rabbit IgG conjugated gold nanorods with 40 nm diameter, longitudinal SPR absorbance max 750 nm. 250 OD concentrated nanorods are also available.

APPLICATION

Application Notes Gold Nanorods conjugates are suitable for dark-field microscopy, two-photon photoluminescence imaging, photoacoustic imaging, kinetic immunoassay, sensors, drug delivery and photothermal therapy.

PRODUCT INFORMATION

Diameter	40 nm
Length	112 nm
Ligand	Affinity purified anti-Rabbit IgG (H+L)
Number of Particles	1.05E+12 particles/mL
Optical Density	50 OD (250 OD are also available.)
Absorption Max	750 nm

STORAGE AND SHIPPING

Storage Buffer Buffer is customer determined.



Storage This product should be stored at 4°C after opening. Care must be taken to only use sterile glassware when working with this product.

Shelf Life 3 months
