

Absolute Mag™ ROR1 Conjugated Magnetic Particles

Cat.No: WHM-AB14

DESCRIPTION

Description	Absolute Mag™ ROR1 Conjugated Magnetic Particles (# WHM-AB14 #) are obtained by coupling the biotinylated ROR1 protein to streptavidin magnetic particles. The ROR1 coupled particles could capture anything binding to ROR1, and make the following testing easy, such as immunocapture, biopanning and flow cytometry. These magnetic particles will be supplied as lyophilized powder (from 0.22 µm filtered PBS pH7.4 including 0.05% Tween-20 and 10% trehalose).
--------------------	--

PRODUCT INFORMATION

Particle Size	2 µm
Ligand	ROR1 protein
Concentration	N/A
Species Reactivity	Human
Usage Statement	a) Reconstitute the magnetic particles, wash the particles and suspended to 1 mg/mL by ultrapure water. b) Add the prepared particles to your samples. c) Particles can be separated from your samples afterwards using a magnetic plate.
Binding Capacity	>150 pmol (25 µg) antibody/mg particles
Surface Group Density	>439 pmol (20 µg) ROR1/mg particles

STORAGE AND SHIPPING

Storage	Upon receipt, please store the particles at -20°C for 1 year in lyophilized state. Once the particles reconstitution, please use it immediately. Do not to freeze thaw the particles after reconstitution.
----------------	--