

Dry Shipper Warming Rack



- Easily warm cryogenic shippers
- Prevent moisture buildup
- Safe and easy storage
- Available in multiple sizes

The Dry Shipper Warming and Drying Rack is an essential tool for maintaining the efficiency and longevity of your dry shippers. Dry shippers, such as the CX100, are used to transport temperature-sensitive biological samples at cryogenic temperatures.



Dry Shipper Air Dryer

Dry shippers can accumulate ice and moisture, reducing LN2 holding capacity. The Dry Shipper Air Dryer is designed to efficiently remove residual water from your dry shipper.

- Fast shipper drying
- Simple, easy design
- Prevents ice build up
- Extends shipper lifetime



Use Code 25ADWR15 for 15% off these items on your next order!

www.mitegen.com
support@mitegen.com

